IFU Part # 5173000
Rev. 1/1806

For complete information see our
MSDS available at our web-site
www.nordiskadental.se

Box 1082, SE-262 21 Ängelholm, Sweden
Location: Framtidsгатan 1b, SE-262 73 Ängelholm
Int tel: (+46) 431 44 33 60 - Fax: (+46) 431 44 33 99
mail@nordiskadental.se - www.nordiskadental.se

**Nordiska
Dental**

Distributed by:

EC REP

ANA 2000[®]

High Copper Admixed Alloy

Non Gamma 2

Capsules

CE
0344

**Nordiska
Dental**

 **Manufactured by:**
Wykle Research Inc.
E. College Parkway 2222, Carson City, NV 89706, USA

ANA 2000[®]
High Copper Admixed Alloy
Non Gamma 2


Composition alloy:
Ag 43,1%, Sn 30,8%, Cu 26,1%
Lathe cut and spherical particles.

Alloy and Mercury acc to ISO 24234

Net mass;	50-Pack	450-Pack
No. 1 	alloy 20,0 g, mercury 20,0 g.	alloy 180,0 g, mercury 180,0 g.
No. 2 	alloy 30,0 g, mercury 30,0 g.	alloy 270,0 g, mercury 270,0 g.
No. 3 	alloy 40,0 g, mercury 40,0 g.	alloy 360,0 g, mercury 360,0 g.

MIXING TIMES

Mixer r/min	Capsule No. 1	Capsule No. 2&3
3500	12 s.	14 s.
4000	8 s.	10 s.
4500	7 s.	9 s.

The above values are valid for mixers such as M4000/M2000, Silamat S5 and Ultramat 2. The time for mixing may vary depending on quality, age and condition of the mixer. For more details see www.nordiskadental.se

INSTRUCTIONS FOR USE

MIXING

Vibrate the capsule according to the "Mixing times".

MANIPULATION

Insert in small portions and condense well. Overpack the cavity. Remove amalgam with excess of mercury. Carve back for occlusal adaption.

The setting time is somewhat shorter than with conventional amalgams, but the time for occlusal adaption is the same.

Recommended storage:  25°C

For dental use only.

Do not use in patients with known allergies.

Dispose of used capsules in accordance with national regulations.

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Danger

H330 Fatal if inhaled.

H360 May damage fertility. May damage the unborn child.

H372 Causes damage to organs through prolonged or repeated exposure.

H410 Very toxic to aquatic life with long lasting effects.

P201 Obtain special instructions before use.

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P273 Avoid release to the environment.


P284 Wear respiratory protection.

P310 Immediately call a POISON CENTER or doctor/physician.

P501 Dispose of contents/container to an approved waste facility.

Contents Capsules,
Mercury 50%

CAS No. 7439-97-6

Effect of mercury on metals

"Mercury reacts with and embrittles particular metals and their alloys. Avoid unnecessary contact between mercury and those metals (and their alloys)."

Spillages

"Mercury presents a health hazard if handled incorrectly. Mercury is toxic by vapour inhalation and the effect is cumulative. Spillages of mercury should be removed immediately, including from places which are difficult to access. Use a plastic syringe to draw it up. Smaller quantities can be covered by sulfur powder and removed. Individual small droplets can be picked up by tin (i.e Sn) foil and removed. Avoid inhalation of the vapour. During and after a clean-up thoroughly ventilate the area where the spill has occurred. Vacuum cleaners should not be used"

Moisture contamination

"If moisture is introduced into the amalgam before it has set, properties such as strength and corrosion resistance can be affected adversely. If the alloy contains zinc, such contamination can result in an excessive expansion (delayed expansion). Use a dry field whenever it is possible."

Waste disposal

"Waste material and all primary containers that have held mercury shall be disposed of following appropriate management practice."