GB Instructions for Use

Algiloid is a reversible hydrocolloid impression material based on agar. It is used in syringes together with an alginate impression material in the so-called combination impression technique. The special composition of Algiloid ensures a firm bonding with most known brands of alginate impression materials. This method combines the accuracy and authenticity of the impressions obtained with the hydrocolloid with the simple use of the alginate impression materials.

Indications

For impressions of inlays/onlays, crown and bridge preparations.

Directions for use

Algiloid is heated up to 98 °C and boiled for 10 minutes. After that, cool down to the working temperature of 68 °C. This can be done by using a special waterbath or a traditional hydrocolloid conditioner. Algiloid is available as ampoules; the ampoules are used in the special ampoule syringe.

Preparation

1. Place the ampoules filled with Algiloid in the conditioner or dry-heater.
2. Follow the instructions for use of the conditioner or dry-heater. Please take care that the hydrocolloids are heated up optimally and that the boiling temperature is adjusted to 98 °C (boiling temperature = 10 minutes). Algiloid may be stored for up to 48 hours at the storage temperature of approx. 68 °C (see also the directions for use of the conditioner or dry-heater). Just prior to use the consistency may be checked by extruding a small amount from the syringe. The material should flow easily and has to be homogenous.

Taking the Impression

1. Follow the usual procedure for the retraction of the gum tissue from the sulcus using retraction cord. Remove the retraction cord prior to taking the impression. After rinsing the preparation, any liquid should be removed. But the hydrophilic nature of Algiloid does not require a dry surface.
2. Apply stops to a perforated tray or a rimlock tray, using a not too soft wax. Fill the tray with slightly thinner than normal alginate impression material.

- We recommend the use of regularly hardening alginate for working with Algiloid.**
3. Apply ALGILOID generously into the sulcus and over the whole preparation. The needle point must always remain in the material to ensure that no air is trapped in the impression.
 4. While the hydrocolloid is still fluid carefully press the tray with the unset alginate over the preparation. Hold firmly in place for 2-3 minutes until the alginate impression material has set.

N.B. The two impression materials will only bond firmly if they come together while they are in a soft condition!

5. Remove the tray from the mouth with a vertical snap action and rinse it under running water.
6. Within 15 minutes cast the impression using a suitable stone and wrap in a damp cloth to harden.
7. In case this time limit can not be kept, the impression should be stored in a moist towel in a well-closed plastic bag. A few drops of water or a moist cotton wool, attached to the grip of the tray, will maintain an atmosphere of 100% R.H., in which the impression can be kept for a maximum of 3 hours without dimensional change.
8. After 1 to 2 hours the gypsum cast may be removed from the impression. Any Algiloid still adhering to the stone cast can be blown away.

Recommended dental stones: Moldano, Moldastone, Moldasynt

Notes

The ampoule-syringe III (REF 101 405) is supplied with two tips (**Set with 2 needles REF 901 124**). The thin needle is used for impressions of inlays/onlays and single crowns. The needle with the bigger diameter is used for larger impressions. It is advisable to prepare plenty of syringes or cartridges. Left over material which has already cooled off may be reused once more after repeating the warming-up procedure.

Warning notices: The ampoule-syringe must be cleaned, disinfected and sterilised before first use as well as after any further use, according to the enclosed instructions.

Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

Manufacturer: HAGER & WERKEN GmbH & Co. KG
Ackerstraße 1, 47269 Duisburg, Germany
T +49 (203) 99269-0, F +49 (203) 299283
www.hagerwerken.de

Product:
Algiloid Ampoule-Syringe III
REF 101 405
As of: 05/2017 (Rev: 4)

Limitation of reprocessing: Frequent reprocessing has only got a small impact on this instrument. The instrument's shelf-life is mainly determined by wear and damage by use.

Place of handling: Clean surface with a disposable cloth or paper towel.

Storage and transport: No special requirements.

Preparations for cleaning: Disassemble product into its components.

Cleaning: Only use cleansers according to the instructions of the respective manufacturer. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C for mechanical cleaning.

Disinfection: Use disinfecting solution according to the instructions of the respective manufacturer. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C for mechanical disinfection.

Maintenance: No special requirements.

Packaging: Standardised packaging material may and should be used for sterilisation processes.

Sterilisation: Vapour sterilisation at 134 °C for 5 minutes,
Vapour sterilisation at 121 °C for 15 minutes.

Control/Functional check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine controls are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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Hersteller: HÄGER & WERKEN GmbH & Co. KG