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General Processing Instructions For Endodontic Products

**APPLICABLE FOR:
ENDODONTIC FILES, BARBED BROACHES, PROBES, EXCAVATORS,
PLUGGERS, SPREADERS, CONDENSORS, DRILLS, BURS,
ULTRASONIC TIPS, FILLING MATERIALS, PINS & POSTS,
WORKING LENGTH TOOLS, OBTURATOR DIAMETER TOOLS.**

1) GENERAL

The processing steps provided herein:

- apply to the products indicated for treatment of endodontic disease before use and/or reuse;
- are intended to be performed only in a clinical or hospital environment, by dentists or qualified users, following good dental practices.

The processing of products applies to the following situations:

- products supplied in a non-sterile condition (no  symbol on the label). Processing prior to first use is required;
- products that are reusable (no  symbol on the label). Processing is required prior to reuse.

The table below summarizes the various existing situations and when the processing instructions are applicable:

Symbols on the label	Processing
None.	Processing before use and reuse.
	Processing before first use only, since reuse is not permitted.
	Processing before reuse only, since the product is ready to use.
 + 	None.

This instruction for use applies to the products listed in ANNEX 1- List of Products covered by this IFU.

Disclaimer: The instructions for processing products prior to first use/reuse herein have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for processing. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the instructions for use set forth below. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

2) WARNINGS

Single-use products shall not be processed and reused. Reusing these products increases the risk of cross-contamination and/or breakage.

3) PRECAUTIONS

GENERAL

- All products marked as “sterile” have an expiry date. Check the expiry date before use and discard the product if it has expired.
- Sterility is not guaranteed once the package has been opened.
- Inspect the packaging and product before use: if the packaging or the product is damaged, do not use the product and discard it.
- Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.
- Always wear protective clothes for your own safety (gloves, glasses, mask and a waterproof gown).
- Do not use labels or identification markers directly on the product.
- Only use properly maintained processing equipment and materials approved by national laws and regulations and use them according to the manufacturer’s instructions (including calibration, cleaning, loading weight, shelf life, operating time and functional testing).
- Only use a detergent solution with disinfecting action approved for its efficacy (VAH/DGHM listing, CE marking, FDA approval) and in accordance with the IFU of its manufacturer.
- Detergent should be aldehyde-free (to avoid blood impurity fixation) and without di- or tri-ethanolamines as corrosion inhibitors.
- Always follow the instructions and concentrations provided by the manufacturer of the cleaning/ disinfecting agent.

LIMITATION ON RE-USE

- The products shall not be processed more than the number of maximum cycles indicated in the table in section 4) STEP-BY-STEP INSTRUCTIONS. It is the responsibility of the user to monitor the number of processing cycles. Meanwhile re-use is permitted only when the product is defect free and following a visual inspection (see item below). Since certain applications may cause the instruments to prematurely reach the end of their useful life the “Max number of processing cycles” will not always be reached (e.g. when shaping an extremely curved root canal with a file).
- In any case inspect the products before re-use, and discard them in case of defects. These include but not limited to the following:
 - Plastic deformation;
 - Breakage;
 - Loss of colour coding or marking;
 - Bent instrument;
 - Untwisted threads;
 - Damaged cutting surfaces;
 - Dull cutting blades;
 - Missing size marking;
 - Corrosion.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible or/and NiTi files (this will help avoid canal transportation).
 - Visually inspect the active part for all the defects listed in the former paragraph during use (i.e after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

MATERIAL RESISTANCE

The use of cleaning and disinfecting agents can damage the products.

Consult the instructions for use provided by the legal manufacturer of the cleaning/disinfection agent and check compatibility with the material of the products (see list of products and materials in ANNEX 1- List of Products covered by this IFU).

Below, a non-exhaustive list of situations to avoid:

- Do not use Hydrogen Peroxide (H₂O₂) solutions on products with plastic handles, NiTi products, tungsten carbide burs or plastic supports. These products are degraded by hydrogen peroxide.
- NiTi products should not be fully immersed in NaOCl solution concentrate. Only the operative part, which is in contact with the patient, should be immersed. The NaOCl solution concentration shall not exceed 5%.
- Do not use solutions containing phenol or any products which are not compatible with the products.
- For metal devices, use anticorrosion disinfecting and cleaning agents.
- Excessive concentrations or immersion times may cause corrosion or other defects in the products.
- It is recommended to use an alkaline detergent with surfactants, with grease removal, disinfection (against bacteria/fungi) and corrosion inhibition properties.
- Products made of aluminum should not be cleaned with liquids containing soda or mercury-salts.
- Products made of aluminum, tungsten carbide or carbon steel should not be used in a washer-disinfector.

4) STEP-BY-STEP INSTRUCTIONS

Step-by-step instructions for processing prior to use/ reuse vary by product type and are divided into 3 categories (Parts A, B and C).

Type of product	Condition	Processing instructions	Max. number of processing cycles
Endodontic Files, Barbed Broaches, Probes, Excavators, Pluggers, Spreaders, Condensers	Processing before first use.	Part A, Steps 4 to 8.	8*
	Processing before reuse.	Part A, Steps 1 to 8.	
Drills, Burs, Ultrasonic Tips	Processing before first use.	Part A, Steps 4 to 8.	5*
	Processing before reuse.	Part A, Steps 1 to 8.	
Pins & Posts	Processing before first use. Product only for single use.	Part A, Steps 4 to 8.	1
Filling Materials & Calcifiable plastic Posts	Processing before first use. Product only for single use.	Part B.	1
Working Length and Obturator Diameter Tools	Processing before first use.	Part C, Step 1 followed by Steps 5 to 9.	50*
	Processing before reuse.	Part C, Steps 1 to 9.	

(*) number of processing cycles validated without affecting the functionality and safety of the products.

A. Processing instructions for Endodontic Files, Barbed Broaches, Probes, Excavators, Pluggers, Spreaders, Condensers, Drills, Burs, Ultrasonic Tips, Pins & Posts.

Please note: Steps “1.” to “3.” are not applicable prior to first use. They only apply for processing prior to reuse for multiple use products.

For products with a silicon stop and supplied non-sterile, steps “4a.” to “8.” are applicable prior to first use. In this case the silicon stop of the new product should be removed before the cleaning (step “4a.” or “4b.”), stored in a pouch in a clean environment, and placed back on the product during inspection (step “5.”).

	Operation	Operating mode	Precautions in addition to section 3) PRECAUTIONS												
1.	Disassembling	- Disassemble the product: Remove and discard the silicon stop, if fitted.	- None.												
2.	Pre-Cleaning	<p><u>General Instructions:</u></p> <ul style="list-style-type: none"> - Soak all products immediately after use (no later than 2 hours after use) in a pre-cleaning solution according to the manufacturer's instructions (aldehyde-free and intended for pre-cleaning by the manufacturer. Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner at 0.4% for a minimum of 15 minutes was used for validation). Where applicable use a tray made from high density polyethylene or stainless steel. - For visible impurities observed on products, or when needed, mechanical pre-cleaning with a soft brush made from either nylon, polypropylene or acrylic is recommended. Manually brush the product until visible impurities are removed. <p><u>Specific Instructions for Diamond Burs:</u></p> <ul style="list-style-type: none"> - Rinse thoroughly immediately after use (no later than 25 minutes after use) and for at least 1 minute under running water (ambient temperature +15°C/+25°C). - Soak all products immediately after use (no later than 30 minutes after use) in a pre-cleaning solution according to the manufacturer's instructions (aldehyde-free and intended for pre-cleaning by the manufacturer. Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner at 0.4% for a minimum of 15 minutes was used for validation). - Brush the active part of the product (coated with diamonds) for at least 10 seconds (a double Nylon toothbrush was used for the validation). - Visually inspect under appropriate lighting (min. 500 lux), if impurities are visible continue to brush until completely removed. 	<ul style="list-style-type: none"> - The pre-cleaning solution should be changed regularly i.e. when it becomes soiled, or when efficacy is diminished due to exposure to microbial loads. - Only use clean soft brushes designed for this purpose. Do not use metal brushes. <p><u>Additional precautions for Diamond Burs:</u></p> <ul style="list-style-type: none"> - For diamond burs, take special care not to exceed the maximum time before rinsing and soaking. If the maximum time is exceeded, cleaning may be ineffective due to the combination of excessive drying of the soiling and roughness of the active part. 												
3.	Rinsing	- Rinse thoroughly (at least 1 minute) under running tap water (ambient temperature).	- If a pre-cleaning solution contains a corrosion inhibitor, it is recommended to perform the rinsing step just before starting the cleaning step.												
4a.	Automated Cleaning with washer-disinfector	<ul style="list-style-type: none"> - Place the products in a kit, support or container to avoid any contact between products (stainless steel baskets were used for the validation). - Place the products in the washer- disinfector (EN ISO 15883 approved) and perform the defined cycle (Ao value > 3000 or, at least 5 minutes at 90°C (194°F)). The following cycle was used for the validation : <table border="1" data-bbox="375 1288 1098 1523"> <tbody> <tr> <td>Pre-cleaning</td> <td>Tap water <45°C ; 2 minutes</td> </tr> <tr> <td>Cleaning</td> <td>0.4% cleaner ; 55°C ; 5 minutes</td> </tr> <tr> <td>Neutralizing</td> <td>Tap water ; 2 minutes</td> </tr> <tr> <td>Rinsing</td> <td>Tap water ; 2 minutes</td> </tr> <tr> <td>Thermal Disinfection</td> <td>Reverse osmosis water ; 90°C ; 5 minutes</td> </tr> <tr> <td>Drying</td> <td>22 minutes</td> </tr> </tbody> </table> <ul style="list-style-type: none"> - Use a detergent solution with cleaning properties (Neodisher Mediclean Forte at 0.4% was used for validation). 	Pre-cleaning	Tap water <45°C ; 2 minutes	Cleaning	0.4% cleaner ; 55°C ; 5 minutes	Neutralizing	Tap water ; 2 minutes	Rinsing	Tap water ; 2 minutes	Thermal Disinfection	Reverse osmosis water ; 90°C ; 5 minutes	Drying	22 minutes	<ul style="list-style-type: none"> - Pay particular attention to cutting edges, both to avoid injury and damage to the products. - Products made of aluminum, tungsten carbide or carbon steel should not be used in a washer- disinfector. Manual cleaning should be performed instead.
Pre-cleaning	Tap water <45°C ; 2 minutes														
Cleaning	0.4% cleaner ; 55°C ; 5 minutes														
Neutralizing	Tap water ; 2 minutes														
Rinsing	Tap water ; 2 minutes														
Thermal Disinfection	Reverse osmosis water ; 90°C ; 5 minutes														
Drying	22 minutes														
OR															
4b.	Manual Cleaning assisted by an ultrasonic device	<ul style="list-style-type: none"> - Immerse in the detergent solution with cleaning properties according to the manufacturer's instructions (Neodisher Mediclean Forte at 2% was used for validation), assisted by an ultrasonic device for at least 15 minutes. - For visible impurities observed on products: manually brush the product until visible impurities are removed (use a soft brush made from either nylon, polypropylene or acrylic). 	- Always place the products in a kit, support or container to avoid any contact between products.												
	Rinsing	- Rinse thoroughly (at least 1 minute) with deionized water (ambient temperature).	- If the cleaning solution used previously contains a corrosion inhibitor, it is recommended to perform the rinsing step just before starting autoclaving.												
	Drying	- Products should be thoroughly dried with a single-use lint-free cloth until visual traces of moisture are eliminated.	- Take particular care to properly dry joints or cavities within a product.												

5.	Inspection	<ul style="list-style-type: none"> - If applicable, assemble the product using a new, unused Dentsply Sirona silicon stop. - Inspect the product functionality (check the sharpness of cutting parts if applicable). - Visually inspect the product with the naked eye under appropriate lighting (min. 500 lux) and discard any with defects (e.g. cracks, deformations (bent, unwound, twisted), breakage, corrosion, loss of colour coding or marking). - Dirty products should be cleaned again. <p><u>Additional processing information for Carbon steel burs:</u></p> <ul style="list-style-type: none"> - Protect carbon steel burs with corrosion inhibitor according to the manufacturer's instructions before packaging (RS-Spray from Meisinger was used for validation). 	<ul style="list-style-type: none"> - The processing and use of the corrosion inhibitor may cause darkening of the carbon steel burs. This colour change is different from the corrosion observed when no corrosion inhibitor is used and should not be considered a defect.
6.	Packaging	<ul style="list-style-type: none"> - Place the products in a kit, support or container to avoid any contact between products. - Pack the products in "Sterilization pouches" (double-packaged using paper-plastic pouches was used as a worst case for steam sterilization validation). - For sharp products that are not contained within a box, silicone tubes should be placed around the products to prevent piercing of the packaging. - Seal the pouches according to the pouch manufacturer's recommendations. 	<ul style="list-style-type: none"> - Ensure that the pouches are suitable for steam sterilization (141°C, 286°F) and validated and manufactured as per ISO 11607 and EN 868-5. - If a thermo-sealer is used, the process must be validated and the thermosealer must be calibrated and qualified.
7.	Sterilization	<ul style="list-style-type: none"> - Place the pouches in the steam sterilizer according to the manufacturer's recommendations. - Use one of the following sterilization cycles with the Pre-Vacuum air Removal steam sterilizer (saturated steam and compliant with EN 13060 (class B, small sterilizer) and EN 285 (full size sterilizer)): <ul style="list-style-type: none"> • 132°C (269.6°F), 4 minutes; • 134°C (273.2°F), 3 minutes; • 134°C (273.2°F), 18 minutes. - We recommend steam sterilization at 134°C / 273.2°F for 18 minutes for the purpose of deactivating potential prions. - Visually inspect the product with the naked eye under appropriate lighting (min. 500 lux) (packaging integrity, no humidity, no colour change of packaging, positive physical-chemical indicators, conformity of actual cycle parameters with the reference cycle parameters). - If visible signs of moisture are present (damp spots on sterile packaging, pooled water in the load) at the end of the sterilization cycle, repackage and re-sterilize using a longer drying time. - Store traceability records and define shelf life according to the packaging manufacturer's guidelines. 	<ul style="list-style-type: none"> - Use a validated sterilization procedure according to ISO 17665 with a minimum drying time of 20 minutes. - Special attention should be paid to the packaging integrity if the 134°C (273.2°F) 18 minutes sterilization cycle is used. - Check the pouch's validity period indicated by the manufacturer to determine the shelf life. - The owner is responsible for complying with the sterilizer's maintenance procedure, which should be performed in accordance with the requirements on sterilizing medical devices (examples: planning for maintenance, qualification, acceptance criteria condensate and water as per EN 285, annex 2).
8.	Storage	<ul style="list-style-type: none"> - Keep products in sterilization packaging in a clean environment, away from sources of moisture and direct sunlight. Store at ambient temperature (typically 15 - 25°C (59 - 77°F)). - In case of damage to the pouch, a complete processing cycle should be performed. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and use by date). 	<ul style="list-style-type: none"> - After sterilization, the product should be handled with care in order to keep the packaging intact (sterile barrier). - Sterility cannot be guaranteed if the packaging is open, damaged or wet.

B. Processing instructions for Filling material and Calcinable plastic posts

	Operation	Operating mode	Precautions in addition to section 3) PRECAUTIONS
1.	Disinfection	<p>Filling material</p> <ul style="list-style-type: none"> - Immerse the obturation device in 5.25% NaOCl for 1 to 5 minutes at ambient temperature. Remove all the bubbles from the surface of the obturation device. - Gently wipe the gutta-percha with sterile gauze moistened with 70% isopropyl alcohol. - Let the obturation devices air dry. <p>Uniclip and Mooser Calcinable plastic posts</p> <ul style="list-style-type: none"> - Immerse the obturation products in NaOCl (2,5% at least) during 5 minutes at ambient temperature. 	<ul style="list-style-type: none"> - Do not use disinfecting solutions containing phenol or any products which are not compatible with the filling material or plastic post.

C. Processing instructions for Working length Tools and Obturator Diameter Tools

Please note: Steps “2.” to “4.” are not applicable prior to first use. They only apply for processing prior to reuse for multiple use products.

	Operation	Operating mode	Precautions in addition to section 3) PRECAUTIONS								
1.	Disassembling	- Disassemble the product, if required.	- None.								
2.	Pre-Disinfection	- Soak all products immediately after use in a detergent and disinfecting solution combined with proteolytic enzyme (according to the manufacturer's instructions. Perfektan Neu - Dr. Schumacher GmbH at 2% for a minimum of 30 minutes was used for validation).	- Change the pre-disinfection solutions regularly i.e. when it becomes soiled, or when efficacy is diminished due to exposure to microbial loads.								
3.	Rinsing	- Rinse thoroughly (at least 1 minute) under running tap water (ambient temperature).	- If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the products just before the cleaning step.								
4.	Pre-Cleaning	<p><u>Complex products such as measuring tools need to be pre-cleaned:</u></p> <ul style="list-style-type: none"> - Immerse the products in disinfecting solution with cleaning properties according to the manufacturer's instructions (only use a CE-marked and FDA-cleared detergent; Perfektan Neu - Dr. Schumacher GmbH at 2% was used for validation). - Make sure that the products are fully immersed. Then wash and brush them in the cleaning/disinfection solution. - Rinse thoroughly (at least 1 minute) under running tap water (ambient temperature). - Repeat the above 3 steps at least three times in order to remove all visible traces/remnants of soiling. 	<ul style="list-style-type: none"> - The pre-cleaning solution should be changed regularly i.e. when it becomes soiled, or when efficacy is diminished due to exposure to microbial loads. - The solution should be a specific solution targeted by the supplier for pre-cleaning and pre-disinfection. It should be used at the dilution specified by the supplier. It should contain or be combined with a proteolytic enzyme. - Only use clean, soft brushes, or a clean, soft cloth or wipe designed for this purpose. Do not use metal brushes or wire wool. 								
5a.	Automated Cleaning with washer-disinfector	<ul style="list-style-type: none"> - Place the products in the washer- disinfector and execute the defined cycle (Ao value > 3000 or at least 5 minutes at 90°C (194°F)). The following cycle was used for the validation : <table border="1" data-bbox="375 1144 1099 1301"> <tbody> <tr> <td>Pre-cleaning</td> <td>30°C ; 1 minute ; 0.5% cleaner</td> </tr> <tr> <td>Cleaning</td> <td>55°C ; 6 minutes ; 0.5% cleaner</td> </tr> <tr> <td>Rinsing</td> <td>1 minute ; demineralized water</td> </tr> <tr> <td>Thermal Disinfection</td> <td>90°C ; 5 minutes ; demineralized water</td> </tr> </tbody> </table> <ul style="list-style-type: none"> - Use a detergent solution with cleaning properties (Neodisher Mediclean-Forte at 0.5% was used for validation). 	Pre-cleaning	30°C ; 1 minute ; 0.5% cleaner	Cleaning	55°C ; 6 minutes ; 0.5% cleaner	Rinsing	1 minute ; demineralized water	Thermal Disinfection	90°C ; 5 minutes ; demineralized water	- Products made of aluminum (e.g. ENDO-M-BLOC) should not be used in a washer- disinfector. Manual cleaning should be performed instead.
Pre-cleaning	30°C ; 1 minute ; 0.5% cleaner										
Cleaning	55°C ; 6 minutes ; 0.5% cleaner										
Rinsing	1 minute ; demineralized water										
Thermal Disinfection	90°C ; 5 minutes ; demineralized water										
OR											
5b.	Manual Cleaning assisted by an ultrasonic device	- Immerse in the disinfecting solution with cleaning properties according to the manufacturer's instructions, and assisted by an ultrasonic device if suitable for at least 15 minutes (e.g. Perfektan Neu - Dr. Schumacher GmbH at 2% was used for validation).	- None.								
	Rinsing	- Rinse thoroughly (at least 1 minute) with deionized water (ambient temperature).	- If a disinfecting solution contains a corrosion inhibitor, it is recommended to perform the rinsing step just before starting autoclaving.								
	Drying	- Products shall be thoroughly dried with a single-use non-woven cloth until visual traces of moisture are eliminated.	- Take particular care to effectively dry joints or cavities within a product.								
6.	Inspection	<ul style="list-style-type: none"> - Visually inspect the product with naked eye under appropriate lighting (min 500 lux) and discard those with defects (e.g. cracks, deformations, breakage, corrosion, loss of marking). - Dirty products should be cleaned again. - If applicable assemble the products and inspect functionality. 	- None.								
7.	Packaging	See Packaging step “6.” provided for A products above.									
8.	Sterilization	See Sterilization step “7.” provided for A products above.									
9.	Storage	See Storage step “8.” provided for A products above.									

ANNEX 1- List of Products covered by this IFU

ENDODONTIC FILES

Type of product	Brands and/or Ref	Materials	Sterile
Manual endodontic instruments	Readysteel® A012X, A012D, A012C, A012I, A011D, A011C, A011I, A016D. Readysteel® Senseus® A1011, A1012, A1013, A1016.	Stainless steel active part, silicone rubber, and plastic or silicone handle.	STERILE R
Manual endodontic instruments	M•access® A12MA, BNCMAFL, A11MA, A16MA. Micro- Opener A0350. Micro- Debrider A0353.	Stainless steel active part, silicone rubber, and plastic handle.	No.
Manual endodontic instruments	PROTAPER® UNIVERSAL For Hand use A0416, A0417, A0418. K-File NiTiflex®* A012N. VERIFIER A0175.	NiTi active part, silicone rubber, and plastic or silicone handle.	No.
Manual endodontic instruments	PROTAPER® UNIVERSAL For Hand use A0416, A0417, A0418. Protaper Ultimate™* BSTPULH.	NiTi active part, silicone rubber, and silicone handle.	STERILE R
Rotary endodontic instruments	PROFILE* A0344, A0358, A011N. PROTAPER® A0409, A0410, A0411, A0415. PROTAPER® RETREATMENT A1410, A1412, A1415. PATHFILE®* A0015.	NiTi active part, silicone rubber, and metal shank (gold, nickel, or Ruthenium plated).	No.
Rotary endodontic instruments	PROTAPER® A0409, A0410, A0411, A0415. PROTAPER® RETREATMENT A1410, A1412, A1415. ProTaper Gold® A0409, A0410, A0411. ProTaper Next® MULTIPLE USE BSTPTNM. TruNatomy®* BSTTNMM. Protaper Ultimate™* BSTPULR.	NiTi active part, silicone rubber, and metal shank (gold, nickel, or Ruthenium plated).	STERILE R

(*) A sterile single use version is also available for these products.

OTHER INSTRUMENTS

Type of product	Brands and/or Ref	Materials	Sterile
Barbed Broaches	A004B.	Stainless steel active part and plastic handle.	No.
Probes	B093.	Stainless steel active part and silicone handle.	No.
Excavators	B0095, B095S.	Stainless steel active part and metal or silicone handle.	No.
Spreaders	A0198.	Stainless steel active part and metal handle (chromium plated).	No.
Finger Spreaders	A0182*.	Stainless steel active part and plastic handle.	No.
	A182N.	NiTi active part and plastic handle.	No.
	A0206.	Stainless steel active part and plastic handle.	No.
Pluggers	A289S.	NiTi or stainless steel active part and silicone handle.	No.
Gutta-Condensors	A0242, A0244.	Stainless steel active part and brass shank.	No.
Lentulo	A022S*, A0022*.	Stainless steel active part and brass shank.	No.

(*) A sterile single use version is also available for these products.

DRILLS & BURS

Type of product	Brands and/or Ref	Materials	Sterile
Enlarger	ROOT CANAL DRILL A0007. GATE A0008*. LARGO® A0009*.	Stainless steel.	No.
Carbide Burs	CAVITY BURS; METAL TRIMMING; SURGERY / BONE CUTTERS; FINISHING BURS EXXXX.	Tungsten Carbide active part and Stainless steel shank.	No.
CST Burs	CAVITY BURS; FINISHING BURS DXXXX.	Carbon steel.	No.
Diamond Burs	DIAMOND BURS ENDO A0164, A0165. OTHER FXXXX.	Stainless steel with diamond coating for the active part.	No.
Stainless steel Burs	SURGICAL BURS D171I, D172I. POST SPACE BURS A0051. THERMACUT A0050.	Stainless steel.	No.

Burs and Drills for Posts/Pins	MOOSER C0053, C0054. STP RESTORATIVE PIN C0058. RS RADIX-STIFTE C0101, C0102. UNIMETRIC C0212, C0213, C0214. Radix C0220, C0216. Easypost™ C0601. CYTCO®-K C107K, C108K.	Stainless steel, some versions with diamond coating on the active part.	No.
--------------------------------	--	---	-----

(*) A sterile single use version is also available for these products.

ULTRASONIC TIPS

Type of product	Brands and/or Ref	Materials	Sterile
Ultrasonic Tips	START-X™ TIP SATELEC A0660, A0661, A9660.	Stainless steel.	No.

FILLING MATERIALS

Type of product	Brands and/or Ref	Materials	Sterile
Gutta-Percha Points	PROTAPER® A022X. PINK A022C. WHITE A022E. COLOUR CODED A022G, A022M, A022N. AUXILIARY A022F, A022L. WaveOne® Gold Conform Fit B00WGGP. ProTaper Next® Conform Fit B00PNGP. ProTaper Gold® Conform Fit B00P. TruNatomy® Conform Fit GGPB00TNGP.	Synthetic gutta-percha.	No.
Gutta-core®	A1703, A1803.	Synthetic gutta-percha, plastic core and silicone rubber.	No.
Obturators	THERMAFIL® A0166, A0167, A0168, A0348. THERMAFIL® For ProTaper-Gold® A1411.	Synthetic gutta-percha & plastic core.	No.

PINS & POSTS

Type of product	Brands and/or Ref	Materials	Sterile
Impression Pin	MOOSER C055B.	Stainless steel.	No.
Fibre Post	Easypost™ C0600. Radix® Fiberpost C0613, C0614.	Composite.	No.
Metal Post	MOOSER C055C, C055D. RADIX ANKER C0193.	Stainless steel.	No.
Metal Post	CYTCO®-K C106K. RADIX ANKER C193T. RADIX ANKER-LONG C0261. RS RADIX-STIFTE C0100. UNIMETRIC C215T.	Titanium.	No.
Pins	RESTORATIVE PIN C0200, C0210, C0211, C0225.	Titanium.	No.
Calcinable Post	MOOSER C0241. UNICLIP C215U.	Plastic.	No.

WORKING LENGTH AND OBTURATOR DIAMETER TOOLS

Type of product	Brands and/or Ref	Materials	Sterile
Working Length Tool	ENDO-M-BLOC A0184.	Aluminium.	No.
Working Length Tool	MINI ENDO-BLOC A0327.	Plastic.	No.
Obturator Diameter Tools	GAUGE A0186, A186C.	Plastic and bronze insert.	No.

Manufacturer

 (for Working Length Tools and Obturator Diameter Tools only)

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