BeeFill® 2in1

de	Gebrauchsanweisung	02-11
en	Operating Manual	12-20
fr	Mode d'emploi	21-30
it	Manuale d'uso	31-40
es	Manual de instrucciones	41-50
pt	Manual de operação	51-60

Distributor: VDW GmbH, Bayernwaldstr. 15
81737 München, Germany

Manufacturer: Aseptico, Inc., P.O. Box 1548
8333 216th Street SE, Woodinville, WA 98072, USA

EU Authorized Representative:
Advena Ltd.,
Tower Business Centre, 2nd Flr,
Tower Street, Swatar, BKR 4013
Malta

PN 420680 Rev. H, ECO 14654, 02/2019

Printed in the USA

Congratulations!

BeeFill® is engineered to provide many years of reliable service. Please read the instructions provided in this manual to receive the best and longest service from your equipment. We recommend the use of a dental dam during all endodontic procedures.

The operating manual can be requested in other languages.

TABLE OF CONTENTS

PACKAGE CONTENTS	12
SPECIFICATIONS	13
CERTIFICATIONS	13
STANDARD SYMBOLS	13
INTENDED USE	14
CONTRAINDICATIONS	14
WARNINGS	14
PRECAUTIONS	14
ADVERSE REACTIONS	14
STEP-BY-STEP INSTRUCTIONS - SETTING UP THE UNIT	15
OPERATION	15-17
STERILIZATION, DISINFECTION & MAINTENANCE	17-18
CLINICAL TECHNIQUE	18
TROUBLESHOOTING	19
FREQUENTLY ASKED QUESTIONS	19
WARRANTY	19

PACKAGE CONTENTS


- Obturation control device with handpieces for Downpack (Pack) and Backfill (Fill)
- Remote power cord
- Bending tool for the cartridges
- Cleaning tool for the Backfill handpiece
- Spare cartridge nut for the Backfill handpiece
- Heat shield for the Backfill handpiece
- Operating manual

Additional products to be used with the device:

- BeeFill® Pluggers:

Small	Black	ISO 40/.03
Medium	Yellow	ISO 50/.05
Large	Blue	ISO 60/.06
- BeeFill® Thermal Response Tip
- BeeFill® gutta-percha cartridges, assorted sizes available in packs of 10
- Machtou-Plugger size 1-2 and 3-4
- NiTi-Plugger Machtou size 0

Figure 1: BeeFill® 2in1 and available accessories

SPECIFICATIONS

Console Dimensions:	11.2 x 15.0 x 15.2 cm (4.4" x 5.9" x 6.0")
Weight:	1.6 kg (3.6 lbs)
Power Source:	115V/60Hz, 230V/50Hz
Current Rating:	115V/0.6A, 230V/0.3A
Fuses:	115V: 0.6A/250V Slo-Blo® fuse 230V: 0.3A/250V Slo-Blo® fuse
Pluggers:	Stainless Steel
Capsule Contents:	Gutta Percha
Capsule:	Aluminium
Cannula:	Silver

Environmental Conditions:

Operating Temperature:	10° to 28 °C (50° to 82.4 °F)
Storage Temperature:	-20 to 60 °C (-4° to 140 °F)
Relative Humidity:	5 to 95 % non-condensing
Altitude:	0 to 3048 meters (0 to 10.000 feet)
Pack Duty Cycle:	25 %

Temperature range Pack:	100 - 400 °C
Volume range Pack:	0 - 100 %

Temperature range Fill:	160 - 200 °C
Flow rate Fill:	20 % - 100 %

Pack handpiece maximum continuous heating time:

10 seconds for temperatures above 200°C.
15 seconds for temperatures at or below 200°C.
One minute in thermal response mode.

RX: This device is restricted to sale by or on the order of a dentist. For dental use only!

Caution: This device has been tested and found to comply with the emissions requirements of IEC 60601-1-2:2001-09. These requirements provide reasonable protection against harmful electromagnetic interference in a typical medical installation. However, high levels of radio-frequency emissions from electrical devices, such as cellular phones, may disrupt the performance of this device. To mitigate disruptive electromagnetic interference, position this device away from radio frequency transmitters and other sources of electromagnetic energy.

CERTIFICATIONS

CE 2797 This product is classified as a Class IIa product and bears the CE marking (CE 2797). The product complies with standards; Europe: IEC 60601-1 + A1:1991 + A2: 1995, Canada: CAN/CSA-C22.2 No. 601.1 and USA: UL 60601-1, First Edition (2003).

Manufacturer: Aseptico, Inc., Woodinville, WA 98072, USA

EU Authorized Representative:
Advena Ltd., Tower Business Centre, 2nd Flr, Tower Street, Swatar, BKR 4013 Malta

Distributor: VDW GmbH, Bayerwaldstrasse 15, 81737 München, Germany
Phone: +49 (0)89 62734-0 • Fax: +49 (0)89 62734-304
info@vdw-dental.com • www.vdw-dental.com

This medical device was put in circulation after August 18th 2005. The product may not be disposed in the normal house garbage. The directives for the waste of electronic and electrical equipment must be considered.

STANDARD SYMBOLS

Fuse ratings
115V: 0.6A
230V: 0.3A
250V SLO-BLO

Type BF
Applied Part

Attention – Consult
Accompanying Documents

Class II Equipment

Dangerous Voltage

Alternating Current

Hot Surface

Pluggers are autoclavable

1 package content

Do not re-use

Keep away from sunlight and heat

Opened packages are not replaced

Caution, this product contains natural rubber latex which may cause allergic reactions.

Cartridge contains gutta-percha

Cannula: silver

Cartridge: aluminium

Expiry date

EU Authorised Representation

INTENDED USE

BeeFill® 2in1 is an integrated obturation device with handpieces that are used for Downpack (Pack) and Backfill (Fill) of root canals. The Pack handpiece is intended to heat pluggers for warming and softening gutta-percha master cones and searing off gutta-percha cones. This handpiece is also intended to heat accessory tips that subject a tooth to heat in order to determine thermal response of pulp vitality. The Fill handpiece is intended for heating and placement of gutta-percha into the root canal. The single-use cartridges are destined for placement of the heated gutta-percha into the cleaned and shaped canal. The BeeFill® device may be used only with original BeeFill® tips and cartridges! A dental dam should be used with any dental procedure!

CONTRAINDICATIONS

Do not use on patients with a known sensitivity to natural rubber latex, silver or copper.

WARNINGS

Equipment not suitable for use in the presence of a flammable mixture with air, oxygen or nitrous oxide.

The gutta-percha cartridges contain natural rubber latex which may cause an allergic reaction.

Prior to replacing the gutta-percha cartridge, retract the plunger completely and switch off the device. Wait until the cartridge nut and the Fill handpiece have cooled off (cool enough to touch). You can then replace the gutta-percha cartridge.

CAUTION: Cartridge as well as Fill handpiece can be damaged if cartridge is replaced while handpiece is still hot.

When replacing a cartridge during a procedure, the handpiece cartridge nut and expended cartridge may be hot to the touch.

CAUTION: Do not remove a hot cartridge!

Risk of fire: Do not contact the heated tip or handpiece with flammable gases or liquids.

Do not use near standing water. Dropping the unit into water may cause electric shock, will destroy the device and can lead to death.

This equipment provides ordinary protection against harmful ingress of liquids. Do not submerge the handpiece or cartridges in any liquid or spray any fluid directly onto the handpiece.

Chemclaving will not be recommended for sterilization and maintenance of the tips as this would cause corrosion.

Plunger tips become hot during use. Turn off unit and allow to cool before exchanging tips.

PRECAUTIONS

A dental dam should be used with any dental procedure.

Always unplug the unit before changing fuses or adjusting the voltage selector.

Pressing the Fill handpiece too tightly during use, may result in a broken cannula.

Place the cannula gently into the canal. Too much pressure will stop the motor from turning.

Take care when replacing a cartridge during a procedure, the cartridge nut and the expended capsule may be hot to the touch. Let the Fill handpiece cool off before replacing the cartridge.

The forward portion of the Fill handpiece becomes warm during use. The heat shield (included with the system) may be used optionally to reduce the handpiece surface temperature. Avoid contact with the forward portion of the Fill handpiece.


When replacing tips during a procedure, the plunger tips may be hot to the touch.

Do not clean unit with a flammable cleaning solution.

ADVERSE REACTIONS

Use in patients with a known sensitivity to latex, silver or copper may cause an allergic reaction. Such an allergic reaction to latex may result in swollen eyes, lips or face. It may also cause difficulty in breathing. The patient should be advised to notify you immediately if any of these symptoms occur.


Figure 2: Component & Control Identification

STEP-BY-STEP INSTRUCTIONS - SETTING UP THE UNIT

1. Unpack the console and check that the Voltage Selector is set to the proper voltage. Use the 115V position for 110-120V 60Hz voltages, and the 230V position for 220-250V 50Hz voltages. Always unplug the unit before changing fuses or adjusting the voltage selector. To change voltage, use a flat head screwdriver to turn the Voltage Selector on the back of the console to the proper voltage setting. To change voltage, replace the fuse to match voltage (see Specifications). Replace with 0,3 A, 230/250 V slow blow fuse (0,6 A for 115 V). For replacement, press on the black plastic clip of the fuse, pull out and replace with the correct fuse. See Fig. 3

Figure 3: Back of device

2. Remove packaging from handpieces. Clean the exterior of the handpieces with a soft cloth moistened with mild non-chlorinated detergent or disinfecting solution which is approved by the DGHM or FDA or bears a CE mark. Also do not use disinfecting solutions which contain phenol, anticorrosives, very acidic or very alkaline solutions - do not submerge handpieces.
Place the handpieces into their respective holders - Pack on the left side and Fill on the right side.

The forward portion of the Fill handpiece becomes warm during use. The heat shield (included with the system) may be used optionally to reduce the Fill handpiece surface temperature. Avoid contact with the forward portion of the handpiece. Sterilize the heat shield before first use and between each patient use. See Sterilization, Disinfection & Maintenance.

3. Connect handpiece cables to console: The handpiece cables are color-coded and keyed to only fit the proper receptacle on the console. With the arrow on the connector facing up, gently push the connector into the console receptacle.
4. Attach the power cord to the back of the console and plug into a grounded electrical outlet.
5. Sterilize Pluggers and Thermo Tip: Autoclave for 10 minutes at 132°C.
6. Pack: Install plunger into Pack handpiece
Push plunger into handpiece tip and slowly rotate in clockwise direction until it keys in with handpiece tip. Push plunger in fully.
7. Fill: Insertion of cartridge
Retract plunger by pressing Return. Allow the plunger to fully retract and switch device off. Wait until cartridge nut and Fill handpiece have cooled off (cool enough to touch). Unscrew and remove the cartridge nut from the handpiece. Insert a cartridge into the handpiece, cannula facing outward. Slip cartridge nut over cannula and screw cap on lightly and in clockwise direction - do not over-tighten.

CAUTION: Cartridge as well as Fill handpiece can be damaged if cartridge is replaced while Backfill handpiece is still hot.

NOTE: Take care when replacing a cartridge during a procedure, the handpiece cartridge nut and the expended cartridge may be hot to the touch. Please switch device off and let it cool before replacing the cartridge.

NOTE: Do not remove a hot cartridge. If the cartridge doesn't fit completely into the handpiece, press the Power switch on the front of the console to the 'ON' position and press the Return button on the console. **The plunger must be in its retracted position to accept the cartridge.**

8. Gently place the heat shield over the cannula and handpiece as necessary. Then turn the heat shield to make the fill indicator visible.

Figure 4: FILL handpiece with heat shield

9. Use the cannula bending tool to place a smooth radius curve on the cannula so the cannula can extend to within 5 mm of the working length of the canal. Place the cannula between the two raised bending posts. Gently bend the cannula to the desired angle.
10. The device is now ready for operation and can be switched back on.

WARNING: Do not touch hot pluggers or hot cartridges.

OPERATION

After the unit has been set up, you are ready for operation as follows:

1. **Turning the device on:**
Press the Power On/Off switch located on the front panel. Power is applied and removed to the device by alternately pressing and releasing this switch.
Buttons on the left side of the keypad are destined for Down-pack, buttons on the right are for Backfill. To change settings use up and down buttons.
2. **Pack or Fill side activation:**
To activate either side of the device, do one of the following: Press the handpiece cuff or any button for the Pack or Fill side. The blue LED indicator at the top of the device will light indicating which side is active. (At initial power on, the Pack side of the device is turned on by default, regardless of which side of the device was last used.) The last used settings are retained in memory.
3. **PACK Side - Program new temperature and volume settings:**
 - a. Press the Pack Temperature button on the left side of the keypad and then use the up/down buttons to set the desired temperature. Press either button repeatedly to increase/decrease the temperature in 10 °C increments, or hold down the button to rapidly scroll through the temperature range (100 °C minimum to 400 °C maximum). The temperature settings are displayed on the numerical LED Window and the yellow Pack Temperature LED illuminates, indicating that the temperature setting mode has been enabled.

NOTE: The thermal response mode temperature is preset to 90° and is not adjustable.

- b. Press the Volume button and then use the up/down buttons to set the desired audible tone volume. Press either button repeatedly to increase/decrease the volume level in 20 % increments, or hold down a button to rapidly scroll through the volume range (0 % minimum to 100 % maximum). The volume settings are displayed on the numerical LED Window and the yellow Volume LED illuminates, indicating that the volume setting mode has been enabled.

NOTE: The numerical LED Window automatically defaults to displaying the selected temperature setting. If the Volume button is pressed, the numerical LED Window will display the volume for 5 seconds and then revert to the selected temperature setting.

4. FILL Side - Program new temperature and flow rate settings:

- a. Press the Fill Temperature button on the right side of the keypad and then use the up/down buttons to set the desired temperature. The LED diode next to the button will flash, indicating that the handpiece is heating up.

NOTE: At temperatures set higher than 190°C, gutta-percha might weep due to the thermal expansion of the gutta-percha material. To avoid this, it is helpful to heat up the Fill handpiece to the desired temperature just shortly before use and to reduce the temperature to 160°C immediately after use. Fill handpiece is now ready to use when you need it and the weepage of gutta-percha is prevented.

- b. Press the Flow Rate Button and use the up/down buttons to set the desired flow rate.

5. Store your custom settings:

To store your new settings for later recall, press either of the two Preset buttons and hold for approximately two seconds. The specific Preset button's LED will illuminate and you will hear a signal when the settings have been successfully stored. Settings for both the Pack side and Fill side are now stored in memory.

6. Preset Use:

The two Preset buttons store all of the settings for both the Pack and Fill sides of the device. The factory default settings are: 200 °C temperature and 40 % tone volume for the Pack side and 160 °C temperature and 60 % flow rate for the Fill side. The Thermal Response button is set at 90 °C and is fixed at that setting for the Pack side. The presets can be manually programmed to store your personal settings by pressing and holding either Preset button for at least two seconds. The Preset button's LED will illuminate and you will hear a signal indicating that it is active.

7. Operation - PACK:

Press the circumferential cuff on the handpiece to begin heating its tip to the selected temperature setting. If the Temperature LED is illuminated, the selected temperature will display on the numerical LED window. As the tip's temperature rises, one tone frequency (approximately 2500 Hz) is heard. Once the tip is within 20 °C of the selected temperature setting, the tone will change to a lower frequency (approximately 2000 Hz) and remain at this tone until the handpiece cuff-switch is released. When pressing the cuff, the tip will heat for a maximum time of 15 seconds if the selected temperature setting is 200 °C or lower. For temperature settings above 200 °C, the tip will heat for no more than 10 seconds before timing out and turning off. When in the thermal response mode, the tip will heat for up to one minute while the cuff switch is depressed. To restart heating after a timeout has been reached, release the handpiece cuff and then press it again.

NOTE: The temperature setting on the console has been optimized based on clinical review. Therefore, the setting does not directly correlate to the absolute temperature of the heated plunger.

CAUTION: Do not stay in the root canal with a heated tip for more than 4 seconds for safety!

CAUTION: Place gutta-percha on the tip of the thermo response tip before applying the thermo tip to the patient's tooth! Do not place the tip directly on the tooth!

CAUTION: Always press the TR button to set the temperature to 90 °C before starting the pulp sensitivity test!

8. Operation - FILL

Dispense gutta-percha: Press the activation cuff to start the flow of filling material. You will notice a short delay as the plunger engages and pushes gutta-percha to the tip of the needle. Extrude a small amount of gutta-percha from the needle. Wipe the excess gutta-percha from the tip before inserting the needle into the canal. Engage the needle within the root canal in accordance with your preferred technique. Hold the handpiece lightly when placing material to allow the device to readily back out of the canal. As the material is expressed, the indicator will help you estimate how much material remains in the cartridge.

NOTE: At temperatures set higher than 190°C, gutta-percha might weep due to the thermal expansion of the gutta-percha material. To avoid this, it is helpful to heat up the Fill handpiece to the desired temperature just shortly before use and to reduce the temperature to 160°C immediately after use. Fill handpiece is now ready to use when you need it and the weepage of gutta-percha is prevented.

NOTE: Applying excessive downward pressure or not allowing the cannula to back-out of the canal may result in a broken cannula.

9. Standby Mode FILL:

After 20 minutes of inactivity, the FILL heater will turn off and slowly cool to room temperature. The Temperature LED will also go out. Press any button on the FILL side of the console to reactivate the FILL heater.

10. Replacing a cartridge

NOTE: When replacing a cartridge during a procedure, the handpiece cartridge nut and expended cartridge are hot to the touch.

1. Select a matching cartridge.
2. Retract the delivery plunger by pressing the Return button.
3. Turn the unit off.
4. Allow the handpiece to cool.
5. **CAUTION: Do not remove a hot cartridge! Cartridge as well as Fill handpiece can be damaged if cartridge is replaced while handpiece is still hot.**
6. Unscrew and remove the Fill handpiece cartridge nut.
7. Using the hole in the included bending tool, remove the cartridge from the handpiece.
8. Dispose of used cartridge in an appropriate biohazard container.
9. Insert a new cartridge, cannula out, into the handpiece end.
10. Slip cartridge nut over cannula and screw on clockwise lightly - do not tighten.
11. Cool handpiece cartridge nut before removing any excess gutta-percha.

11. Replacing the fuses

NOTE: The BeeFill® 2in1 is manufactured with 250V 300 mA rated fuses installed for a 230VAC power source. If power source is 115 VAC, ensure that 250V 600 mA rated fuses are installed.

WARNING: Turn the power off and unplug the unit before following the steps below.

1. Remove the fuseholder from the power inlet connector.
2. Replace the fuses in the fuseholder.

Replacement Fuses

230V: 300 mA, 250V rated slow-blow type (Fuse size: 5 x 20mm)
115V: 600 mA, 250V rated slow-blow type (Fuse size: 5 x 20mm)

3. Replace the Fuseholder.

AC POWER INLET CONNECTOR

FUSE HOLDER

Figure 5: Fill handpiece components



Figure 6: Pack Handpiece components

STERILIZATION, DISINFECTION & MAINTENANCE

Control Console

Clean the exterior of the console by wiping with a soft cloth moistened with mild non-chlorinated detergent or disinfecting solution which is approved by the DGHM or FDA or bears a CE mark. Also do not use disinfecting solutions which contain phenol, anticorrosives or very acidic or very alkaline solutions.

CAUTION: DO NOT SUBMERGE the device in any fluid.

NOTE: When wiping down the handpiece cable, gently wipe from the middle of the cable out to the handpiece and console. Avoid gripping the cable tightly.

Handpieces

Clean the exterior of the handpiece with a soft cloth moistened with mild non-abrasive, detergent or disinfection solution which is approved by the DGHM or FDA or bears a CE mark. Also do not use disinfecting solutions which contain phenol, anticorrosives, very acidic or very alkaline solutions.

CAUTION: DO NOT SUBMERGE the handpieces in any fluid or spray any fluid directly onto the handpieces.

PACK:

Pluggers and Thermal Response Tip

Pluggers and tips must be cleaned, disinfected and sterilized before every use:

1. Gently use a brush for removal of debris. Clean the disinfected pluggers and tips under running water for a minimum of 1 minute. Clean with water and mild, non-abrasive, detergent or disinfection solution which is approved by the DGHM or FDA or bears a CE mark. Also do not use disinfecting solutions which contain phenol, anticorrosives, very acidic or very alkaline solutions. Clean the disinfected pluggers and tips under running water for a minimum of 1 minute to remove all chemicals. Dry thoroughly.

2. Autoclave plugger tips for 10 minutes at 132 °C before initial use and between each patient use.

3. Observe that the pluggers and thermal response tip you are using are not damaged. A deformed or oxidized tip should be replaced. All tips will gradually lose their heating efficiency over a period of time.

4. The tips can be re-used, if they are treated with care, are not damaged or contaminated. Any further use of damaged or contaminated tips is at the user's own risk. All risk and liability is excluded in this case.

NOTE: Observe the concentration of the cleaning or disinfecting solution stated by the manufacturer!

WARNING: Do not touch hot pluggers!

FILL:

Handpiece Cartridge Heater

If gutta-percha gets into the cartridge heater section of the handpiece, turn on the unit and press the Return button on the console to fully retract the plunger. Allow the heater section to reach 180°C. Turn off unit. Insert the included handpiece cleaning brush into the heating chamber. Rotate the brush several times to remove the gutta-percha from the chamber.

Handpiece Delivery Plunger

Annually: Without a cartridge inserted, press the activation cuff until the gutta-percha indicator moves all the way forward. Then, push the Return button and allow the indicator to completely retract.

Handpiece Cartridge Nut

Allow the Fill handpiece cartridge nut to cool before removing any excess gutta-percha. The cartridge nut may be steam autoclaved for 10 min at 132°C.

Cartridges

Cartridges are intended for single patient use. Prior to using the device on a patient, wipe the cannula with alcohol or a disinfectant. (Make sure that the device has been switched off and cooled down.) Store cartridges at room temperature. Do not immerse the cartridges into any liquid. Dispose of cartridge in a biohazard container. The cartridges must not be used after the expiry date!

Heat Shield

Autoclave for 15 minutes at 132 °C.

NOTE: All liability is disclaimed for failure to follow these instructions or use of non-validated methods for the re-use of instruments!

CLINICAL TECHNIQUE

Only use filling materials compatible with this filling technique such as gutta-percha or root canal sealer. The actual temperature in the root canal largely depends on the amount of filling material used. Although the device is able to precisely control the heat transmission to the tip, it is unable to prevent undesired heating up of the root. We recommend testing on extracted teeth.

ATTENTION: Do not penetrate the root canal with a heated tip for more than 4 seconds for safety!

ATTENTION: Recommended maximum temperature setting is 200°C for all downpack techniques!

NOTE: Do not touch the lips, the gums or the oral mucosa with the plugger the cannula or the handpiece cartridge nut, as it may get hot to the touch after longer exposure.

INSTRUCTIONS FOR USE

Cone Fit & Plugger Selection

1. Prepare the canal optimally recognizing shaping facilitates 3-D cleaning and filling.
2. Select a larger size manual plugger that will work passively and effectively over a range of a few millimeters in the coronal one-third of the canal.
3. Select a medium size manual plugger that will work passively and effectively over a range of a few millimeters in the middle one-third of the canal.
4. Select a smaller size manual plugger that will work passively, effectively and deeper in the straight portion of the canal and to within 4-5 mm of the canal terminus.
5. Select an electric heated plugger that will passively fit through the straight portion of the canal and optimally to within 5 mm from the working length. Set the silicone stop at this depth to promote safety and accuracy.
6. In a fluid-filled canal, fit either a non-standardized, fully tapered gutta-percha master cone or a VDW 4 % or 6 % tapered alpha gutta-percha to working length until it exhibits apical tugback. Please confirm radiographically.
7. Dry the canal with appropriately sized paper points.
8. Trim the master cone back to the canal terminus based on the paper point drying technique.
9. Lubricate the master cone lightly with sealer and gently insert it to length.

Downpack: Continuous Wave Technique

1. Activate the electric heat plugger and seal off the master cone at the level of the orifice.
2. Select the larger manual plugger and compact its working end around the circumference of the canal using short, firm strokes to move gutta-percha apically, clean the canal walls and flatten the material.
3. Use this larger manual plugger and press for five seconds to compact warm gutta-percha vertically and laterally into this region of the root canal system (first wave of condensation).
4. Activate the electric heated plugger and, in one deliberate and continuous motion, press the heated plugger through the thermosoftened gutta-percha until the silicone stop is 2 mm from the reference point. This procedure must be limited to two to four seconds to prevent thermal injury.
5. Deactivate the plugger and continue to maintain firm apical pressure on the cooling instrument until the silicone stop reaches the reference point.
6. Maintain firm apical pressure for ten seconds to compact the mass of warm gutta-percha into the apical one-third of the root canal system and to offset shrinkage during the cooling phase.
7. Activate the electric heated plugger for one second, then deactivate and remove the plugger from the root canal using a back and forth motion. This procedure separates and removes gutta-percha from the coronal two-thirds of the canal without disturbing the gutta-percha in the apical one-third.
8. Select the small manual plugger and compact its working end around the circumference of the canal using short, firm strokes to clean the canal walls and re-condense the most coronal aspect of the gutta-percha within the apical one-third.
9. Use Fill handpiece to optimally fill the rest of the canal.

Backfill

1. Position the tip of the warm cannula against the previously packed filling material for 5 seconds.
2. Press the activation cuff on the Fill handpiece and dispense a small (a few millimeters) segment of warm gutta-percha into this region of the canal. Hold the handpiece lightly so it will back out of the canal during use.
3. Select the smaller manual plugger and compact its working end circumferentially around the canal using short, firm strokes to clean the canal walls and flatten the dispensed warm gutta-percha.
4. Use the same small manual plugger and press for five seconds to three-dimensionally compact warm gutta-percha into this region of the canal and to offset shrinkage during the cooling phase.
5. Position the tip of the warm cannula against the previously packed filling material for five seconds.
6. Press the activation cuff on the Fill handpiece and dispense a few more millimeters of warm gutta-percha into this region of the canal.
7. Select the medium size manual plugger and compact its working end circumferentially around the canal using short, firm strokes to clean the canal walls and flatten the dispensed warm gutta-percha.
8. Use the same medium size manual plugger and press for five seconds to three-dimensionally compact warm gutta-percha into this region of the canal and to offset shrinkage during the cooling phase.
9. Continue the backfilling technique, in the manner described, until the canal is completely filled or stop at any point to accommodate a post to facilitate the restorative needs.

TROUBLESHOOTING

1. **Device does not turn on:**
 - a) Verify that power cord is plugged into the device and into wall outlet.
 - b) Unplug device and check fuse. If fuse is blown, replace it. See Operation, Replacing the fuses
2. **Pack Handpiece does not work:**
 - a) Make sure the tip is sitting in the handpiece correctly.
 - b) Replace tip.
3. **Pack: The tip does not get warm:**
 - a) Check temperature settings.
 - b) Replace tip with a new one.
4. **Fill: Cannot remove handpiece cartridge nut:**
 - a) Press Return button and allow drive mechanism to back out in order to relieve pressure on cartridge nut.
5. **Fill: Gutta-percha does not flow through cannula:**
 - a) Verify that the cannula is not kinked - avoid over bending of cannula.
 - b) Verify that the device has reached operational temperature for gutta-percha to flow.
 - c) Increase temperature if needed.
6. **Fill: Delivery plunger does not retract:**
 - a) If, after pushing the Return button, the gutta-percha indicator does not fully retract to the end of the indicator window, push the Return button again.
7. **Fill: Motor stops turning:**
 - a) If the motor stops turning, the cannula may have been placed too firmly in the canal. Reduce the pressure applied to the cannula and the motor may begin turning again.

FREQUENTLY ASKED QUESTIONS

What size is the smallest tip?

The smallest is the small black plugger with ISO 40 Taper 03. The material of the outer surface is made out of stainless steel.

How hot can a tip get?

At maximum power, a tip can reach 400 °C. You can adjust the heat down to 100 °C by using the lowest power setting. The maximum recommended working temperature for downpack is 200 °C.

Isn't this high temperature dangerous?

Like other dental equipment, the unit should only be used by a trained professional. You can use the high heat setting to sear off gutta-percha in the lateral condensation technique where a cluster of cones are cut through. If you want to remove gutta-percha, with continuous wave technique you should use a lower heat setting than 200 °C.

What is the normal power setting?

The maximum heat setting should be used only when searing off gutta-percha. For the vertical condensation technique and continuous wave technique, we recommend a power setting of max. 200 °C.

How long will the tip last?

It depends on how well you take care of the tip. The life of a tip will be reduced by:

- a. running it too hot for long periods of time
- b. working with high temperatures
- c. bending or putting too much mechanical force on it.

The tips should be used only to transfer heat – not to exert manual force! Our recommendation is that you use a cold manual plugger for that purpose. All tips will gradually lose their heating efficiency over time. This is strongly depending on the user.

WARRANTY

VDW warrants the BeeFill® 2in1 device against defects in material or workmanship for a period of one year from date of original invoice. VDW does not warrant any other products that come along with the device, such as pluggers, thermo tips or cartridges for example, as they are consumables.

VDW's sole obligation under product warranty is (at its sole option and discretion) to repair or replace any defective component or product in part or whole. VDW shall be the sole arbiter of such action.

In the event of alleged defect under warranty, the purchaser is to notify VDW Customer Complaint Department promptly. Customer Complaint will provide instructions, usually directing that the product be returned for service. Shipment to VDW and the cost thereof is always the responsibility of the purchaser.

Accidental misuse, inappropriate installation, or failure to perform directed maintenance voids the warranty.

VDW does not assume, under this warranty, any risks or liabilities arising from the clinical use of its products, whether or not such use involves coincidental utilization of products manufactured by others.

VDW makes no warranty other than that stated above, expressed or implied.

IN CASE YOUR DEVICE IS DEFECTIVE PLEASE RETURN IT TO VDW AND INCLUDE THIS SERVICE FORM

Your address (or stamp): _____

Name: _____

Street: _____

Post code, place, country: _____

Telephone: _____

Fax: _____

E-mail: _____

Contact person: _____

The device was purchased on (date): _____

at (e.g. dealer, or sales rep.): _____

Please clean and sterilize the device before returning it, and include all accessories:

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Description of defect: _____

- Repair within the warranty (invoice necessary)
- If repair costs exceed _____ € please send an estimate.
- Please repair the device, I do not need an estimate.

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
SPECIFICATIONS

Console Dimensions: 11.2 x 15.0 x 15.2 cm (4.4" x 5.9" x 6.0")
Weight: 1.6 kg (3.6 lbs)
Power Source: 115V/60Hz, 230V/50Hz
Current Rating: 115V/0.6A, 230V/0.3A
Fuses: 115V: 0.6A/250V Slo-Blow® fuse
 230V: 0.3A/250V Slo-Blow® fuse
Pluggers: Stainless Steel
Capsule Contents: Gutta Percha
Capsule: Aluminium
Cannula: Silver

Environmental Conditions:
Operating Temperature: 10° to 28 °C (50° to 82.4 °F)
Storage Temperature: -20 to 60 °C (-4° to 140 °F)
Relative Humidity: 5 to 95 % non-condensing
Altitude: 0 to 3048 meters (0 to 10,000 feet)
Pack Duty Cycle: 25 %

Temperature range Pack: 100 - 400 °C
Volume range Pack: 0 - 100 %

Temperature range Fill: 160 - 200 °C
Flow rate Fill: 20 % - 100 %

Pack handpiece maximum continuous heating time:
 10 seconds for temperatures above 200°C.
 15 seconds for temperatures at or below 200°C.
 One minute in thermal response mode.

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Manufacturer: Aseptico, Inc., Woodinville, WA 98072, USA

EU Authorized Representative:
 Advena Ltd., Tower Business Centre, 2nd Flr, Tower Street, Swatar, BKR 4013 Malta

Distributor: VDW GmbH, Bayerwaldstrasse 15, 81737 München, Germany
 Phone: +49 (0)89 62734-0 • Fax: +49 (0)89 62734-304
 info@vdw-dental.com • www.vdw-dental.com

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Fuse ratings
 115V: 0.6A
 230V: 0.3A
 250V SLO-BLO

Type BF
 Applied Part

Attention – Consult
 Accompanying Documents

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1. Turning the device on:

Press the Power On/Off switch located on the front panel. Power is applied and removed to the device by alternately pressing and releasing this switch.

Buttons on the left side of the keypad are destined for Down-pack, buttons on the right are for Backfill. To change settings use up and down buttons.

2. Pack or Fill side activation:

To activate either side of the device, do one of the following: Press the handpiece cuff or any button for the Pack or Fill side. The blue LED indicator at the top of the device will light indicating which side is active. (At initial power on, the Pack side of the device is turned on by default, regardless of which side of the device was last used.) The last used settings are retained in memory.

3. PACK Side - Program new temperature and volume settings:

- a. Press the Pack Temperature button on the left side of the keypad and then use the up/down buttons to set the desired temperature. Press either button repeatedly to increase/decrease the temperature in 10 °C increments, or hold down the button to rapidly scroll through the temperature range (100 °C minimum to 400 °C maximum). The temperature settings are displayed on the numerical LED Window and the yellow Pack Temperature LED illuminates, indicating that the temperature setting mode has been enabled.

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