Celtra® Press investment

Instructions for use

Gebrauchsanweisung
Consignes d'utilisation
Istruzioni per l'uso
Traducción del manual del usuario

Celtra Press investment

EN	Instructions for use.....	05
DE	Gebrauchsanweisung.....	17
FR	Consignes d'utilisation	29
ES	Traducción del manual del usuario	41
IT	Istruzioni per l'uso	53

Celtra Press investment

EN	Instructions for Use	06
1	Technical specifications	07
2	Before you start (hazard and safety warnings)	07
3	Storing the investment	08
4	Investing	09
4.1	Preparations	09
4.2	Measuring the liquid	09
4.3	Mixing	10
4.4	Investing technique	10
5	Setting times	10
6	Pre-heating	11
6.1	Conventional heating	11
6.2	Quick heating	11
7	Pressing	12
8	Divesting	12
	Investing with Celtra Press investment - Quick Reference Guide	13
	Tables	15

Celtra Press investment

Assures good pressing results – and protects your health!

You have chosen the **Celtra Press investment** system. This phosphate-bonded investment was specially developed for pressing lithium disilicate ceramics. Its special composition minimizes the formation of a reaction layer on the ceramic surface. It is easy to air-abrade, protecting your objects during divesting. The object no longer requires acid treatment after air-abrading, saving you valuable time and avoiding the use of potentially dangerous chemicals in the lab. **Celtra Press investment** is also ideal for the conventional ceramic pressing technology.

WHY CHOOSE A DUST REDUCED INVESTMENT?

Respirable crystalline silica (which is found in every investment) is considered carcinogenic and is but one substance on a long list of carcinogenic substances. We recognised the dangers involved in working with such materials and the need to develop a dust reduced investment.

Dental technicians are exposed to hazardous substances on a daily basis. Being under pressure in the laboratory may sometimes lead to the neglect of health and safety precautions. **Celtra Press investment** tackles this problem at its source! Very low dust levels during handling mean one less hazard for the technician at work!

QUALITY TROUGHOUT

Within the framework of a certified quality management system, each batch is checked for correct product properties during the manufacturing process. The final product, too, is subjected to strict quality checks before being released for delivery to the customer. This gives the dental technician a reliable investment material that will yield top results, again and again.

APPLICATIONS

Celtra Press investment may be used for the standard pre-heating process with holding levels or for the quick-heating process (“Quick” or “Speed”) at the final temperature.

1 Technical specifications

Table 1: Technical specifications

Mixing ratio (powder: liquid)	100 g: 22 ml
Mixing time (under vacuum)	60 s
Processing time	4,5 - 9 min
Setting time (Vicat)	5,5 - 12 min
Setting expansion	0,3 - 1,2%
Compressive strength (depending on mixing liquid concentration)	4 - 15 MPa

2 Before you begin...

... we would like to ask you to take two minutes to read our **hazard and safety warnings**.

Even though **Celtra Press investment** has an extremely low level of dust emission, in the interest of health we recommend:

DANGER!

When processing the powder:

- Use a local dust extraction unit.
- In case of insufficient ventilation, wear suitable respiratory equipment.
- Do not inhale dust.
- This material contains quartz, which may cause lung disease or silicosis on prolonged or repeated exposure when inhaled.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Do not mix the powder with other materials.
- Dispose of the contents/container in accordance with local regulations.

When processing the liquid:

- Use local dust removal equipment or
- Respiratory protection
- Caution: Mixing liquids are caustic.
- If mixing liquid has come into contact with the skin, flush with plenty of water immediately.
- When heating above 200 °C (391 °F): Ammonia may form. Ammonia irritates the respiratory system.
- If the liquid comes into contact with the eyes:
Rinse with a copious amount of water while keeping the eyelid open.
If complaints persist, consult an ophthalmologist.
- Do not mix the liquid with other liquids.
- Do not mix the powder with other materials.
- Protect the liquid from temperature below freezing.

When heating above 200 °C /390 °F):

- Formation of ammonia vapors
- Irritation and caustic reactions to the skin and mucosa (eyes, respiratory passages, gastrointestinal tract after swallowing) must be expected on local contact.

Celtra Press investment must never come in contact with gypsum or with substances containing gypsum. Failure to heed this warning may adversely influence the properties of the material. For this reason, separate mixing bowls and spatulas must be used for **Celtra Press investment**.

3 Storing the investment

Celtra Press investment should be stored generally at 5-30°C (41-86 °F). Your experience will have shown you how external factors (such as summer/winter) may influence your pressing results. For your current batch we therefore recommend a cool and dry place that generates temperatures 20 -23°C (41-86 °F) to generate constant results.

Please note the symbols on product labels:

- REF Product code
- LOT Batch or lot number
-  Use before
-  Consult instructions for use

Avoid unpleasant surprises

- Do not use investment that has been stored in a too cold or warm place (warehouse, basement, or similar). Replenish your stock in time - well before the last box/package has been used - in order to ensure that the material used has the recommended temperature.
- Do not use boxes/packages that have just been delivered. Make sure they reach the recommended temperature first.
- Regularly check the temperature in the investment area.
- If you have enough storage space, it is recommended to keep enough material (especially the liquid) in stock to be able to avoid wintertime deliveries (frost hazard).

4 Investing

4.1 Preparations

Our recommendation: *Always use separate mixing bowls and spatulas.*

- Look at your mixing bowl (and the mixing motor). Are they really clean? No investment residue, no other dirt, dust, or stains? Replace old bowls. Use only clean spatulas.
- If you use wax tension-reducing agents, make sure that no residue adheres to the wax pattern.
- The investment ring, bases and covers must be clean.
- Check to ensure that the wax pattern is located in the middle of the investment ring (level, position, sprues, etc.)

4.2 Measuring the liquid

Mix **Celtra Press investment** only with **Celtra Press investment Liquid**. Diluting with distilled water lets you control the expansion of the investment. The higher the concentration of the mixing liquid, the higher the total expansion of the investment.

Table 2: Mixing ratios

Ring size	100g investment ring	200g investment ring
Powder	100 g	200 g
Liquid	22 ml	44 ml

We have compiled a table for the various dilution ratios of the mixing liquid (Table 2). These list different concentrations as percentage. Of course, you might convert these percentages to ml and then create each mix individually for each investing process. This procedure, however, is not without problems:

- *Measuring small quantities in the measuring cylinder is difficult.*
- *Conversions may result in “odd” numbers (e.g. 28.05 ml) that cannot be measured in the measuring cylinder.*
- *If the measuring cylinder is not clean or too old, dried mixing-liquid residue may cause your results to be wrong.*

We therefore recommend pre-mixing your mixing liquid in a calibrated 1-litre cylinder (these cylinders can be ordered from chemical or dental lab suppliers). Your dosage will be more exact, and you will have a constant mixing ratio over the “life” of this one litre of liquid. You might also pre-mix different ratios and store them in plastic bottles.

4.3 Mixing

MIXING RATIO (POWDER : LIQUID) = 100 g : 22 ml

- Pour the exact amount of liquid into the mixing bowl.
Caution: Old mixing bowls with a rough surface can absorb 2–3 ml when dry. Always rinse the bowl with water first and then wipe dry.
- Add the exact amount of powder or the content of the bag. (That is, all of the contents – if anything is spilled, discard the materials and start over.)
- Carefully pre-mix manually for **30 (thirty!) seconds** Close the mixing bowl. If possible, place the bowl in a vacuum for an additional **30 seconds** (*check the vacuum regularly*).
Mix with a mixer under vacuum for another **60 seconds**.

Do not fall for the mixer trap:

The speed of your mixer affects the setting expansion of the investment. All values in these instructions refer to a mixer speed of 360 rpm. If the speed of your machine differs from this value, this may affect your pressing results.

4.4 Investing technique

Pour mixed investment into the investment ring using slight vibration. Avoid trapping air. The material's thin-flowing consistency and its processing time of 6–8 minutes allow easy investing even of difficult patterns without being strapped for time. Once the wax pattern is covered with investment, terminate the vibration action immediately.

5 Setting times

If the investment is to be quick-heated, the investment ring must be inserted at the final temperature, no sooner than 20 minutes after the mixing began. To make sure this time interval is observed exactly, simply set a timer to 20 minutes directly before mixing.

- If the investment is to be entered **conventionally, with holding stages**, we recommend the following setting times, depending on the size of the investment ring:

100g investment ring	=	30 min
200g investment ring	=	30 min

At the end of the setting time, place the investment ring in the cold furnace. Under no circumstances must the investment ring be allowed to dry out.

- If the investment is preheated by **quick heating**, we recommend a setting time of 20 min.

6 Pre-heating

6.1 Conventional heating

At the end of the setting time, place the investment ring to be pre-heated in the **cold furnace** (“cold” meaning less than 100 °C (212 °F)). Depending on the type of furnace used, a program with a heating rate of 7 °C/min (45 °F/min) will be started with the following holding levels and holding times (depending on the size of the investment ring):

- 1st holding level, 45–60 min at 270 °C (518 °F) (cristobalite transition)
- 2nd holding level, 30–60 min at 580 °C (1076 °F) (quartz transition)
- 3rd holding time at final temperature (Celtra Press: 850 °C / 1562 °F):

100g investment ring	=	45 min
200g investment ring	=	60 min

6.2 Quick heating

The following instructions apply to **quick heating**:


- Observe the prescribed setting time of 20 min! The investment ring must not be placed in the hot furnace until the setting time has elapsed.
- While the investment can tolerate a longer setting time (20 to 40 minutes), this may adversely affect the fit of the objects or even lead to cracking.
- Place the ring after the appropriate setting time in the preheated furnace at end temperature acc. the recommendations of the ceramic manufacturer.

We recommend placing the investment ring on a firing support or similar with the opening down during heating to avoid contact with the bottom plate of the chamber. (This ensures that the investment ring is heated uniformly and that the wax can run out freely.)

Caution: During the burnout process, opening the oven door will generate large flames.

The following holding times at the final temperature (Celtra Press: 850 °C / 1562 °F) must be observed:

100g investment ring	=	45 min
200g investment ring	=	60 min

We have summarized this information for you in our Quick Reference Guide (Page 13).

7 Pressing

Meticulously follow the instructions for Celtra Press or, if applicable, the instructions of the pressable ceramic manufacturer and the Instructions for Use for your pressing furnace.

8 Divesting

Caution! Celtra Press investment loses its dust reducing attributes after it has been heated in the furnace!

Allow the ring to cool to room temperature after pressing. *Never quench the investment ring in cold water.* The object can be divested easily after cooling:

- Do not use divesting nippers for divesting ceramic objects.
- With glass beads (50 µm, max. 2 bar).
- Or with aluminium oxide (110 µm, max. 1.5 bar).
- No further processing or treatment with acid is required.

We wish you success with the **Celtra Press investment** compound.

Your **Dentsply Sirona** team

INVESTING WITH CELTRA PRESS INVEST

QUICK REFERENCE GUIDE

1. Mixing ratio: 100 g: 22 ml
2. Prepare and clean tools (mixing bowl, mixer, spatula)
3. Prepare the investment ring.
4. Place sprues as recommended in the Instructions for Use for **Celtra Press**.
5. If using wax tension-relieving agents, blow dry carefully and leave no residue.
6. Quick processing? Get the timer ready.
7. Measure liquid (following Table 1 in chapter 4.2). Temperature?
8. Measure powder.
9. Add the powder to the liquid. CAUTION: For Quick/Speed you should now start the timer (20 minutes!).
10. Premix well for approximately 30 seconds.
11. Close bowl and place in vacuum for 30 seconds. Check vacuum.
12. Start mixing process (60 seconds).
13. Fill investment ring.
14. For Speed/Quick heating: Insert directly when the 20 minutes are over.
15. Check for correct placement in the furnace.

Our tip:

Remove this Quick Reference Guide from your User Manual and display it where you do your investing.

Table 3: Mixing ratios for different concentrations of the mixing liquid

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Table 4: Recommended concentrations of the mixing liquid

		Celtra Press	Your suggestion:
One-surface Inlay	Wax CAD CAM wax	35% - 45%	
MOD inlay	Wax CAD CAM wax	35% - 45%	
Anterior crown	Wax CAD CAM wax Adapta	50% - 60%	
Molar crown	Wax CAD CAM wax Adapta	50% - 60%	
3- to 4-unit bridge	Wax CAD CAM wax Adapta	55% - 65%	

Please note that numerous processing parameters (e.g. temperatures, modelling materials, vacuum mixers, etc.) may influence the fit of the object and that the above recommendations may have to be adapted accordingly.

Table 5: Heating times: Conventionally controlled furnaces

Ring size	100 g	200 g
Setting time	30 min	30 min
Starting temperature	100 °C (212 °F)	100 °C (212 °F)
Temperature increase	7 °C/min (45 °F/min)	7 °C/min (45 °F/min)
1st holding level	270 °C (518 °F)	270 °C (518 °F)
Holding time	30 min	40 min
Temperature increase	7 °C/min (45 °F/min)	7 °C/min (45 °F/min)
2nd holding level	580 °C (1076 °F)	580 °C (1076 °F)
Holding time	20 min	30 min
Temperature increase	7 °C/min (45 °F/min)	7 °C/min (45 °F/min)
3rd holding level*	Final temperature	Final temperature
Holding time	45 min	60 min

* Celtra Press: 850°C/1562°F

Table 5a: Accelerated heating

Ring size	100 g	200 g
Setting time	20 min	20 min
Insert temperature	Up to 600 °C (1112 °F)	Up to 600 °C (1112 °F)
Holding time	N/A	N/A
Heating to final temperature	Immediately	Immediately
Holding time at final temperature	45 min	60 min

Table 6: Heating times: Quick (Speed)

Ring size	100 g	200 g
Setting time	20 min	20 min
Insert temperature	Final temperature	Final temperature
Holding time	45 min	60 min

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DeguDent GmbH