- EN INSTRUCTIONS FOR USE
- DE GEBRAUCHSANWEISUNG
- FR MODE D'EMPLOI
- IT ISTRUZIONI PER L'USO
- ES INSTRUCCIONES DE USO
- AR تعليمات الاستخدام
- CG UPUTSTVA ZA UPOTREBU
- CZ NÁVOD K POUŽITÍ
- DA BRUGSANVISNING
- EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
- ET KASUTUSJUHEND
- FI KÄYTTÖOHJEET
- HU HASZNÁLATI UTASÍTÁS
- LT NAUDOJIMO INSTRUKCIJOS
- LV LIETOŠANAS INSTRUKCIJA
- NL GEBRUIKSAANWIJZING
- NO BRUKSANVISNING
- PL INSTRUKCJA UŻYCIA
- PT INSTRUÇÕES DE UTILIZAÇÃO
- RO INSTRUCȚIUNI DE UTILIZARE
- RU ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
- SK NÁVOD NA POUŽITIE
- SR UPUTSTVO ZA UPOTREBU
- SV BRUKSANVISNING
- TR KULLANIM TALIMATLARI

# CEREC Tessera™

## Advanced Lithium Disilicate

for CEREC® and inLab® machines

### INSTRUCTIONS FOR USE – ENGLISH

**Caution:** This is a medical device. For dental use only.  
USA: Rx Only

#### 1. PRODUCT DESCRIPTION

CEREC Tessera™ CAD/CAM advanced lithium disilicate block for CEREC® and inLab® is a material that can be used for the fabrication of full and partial full contour crowns in the anterior and posterior region using a CAD/CAM procedure. Indirect restorations are fabricated by grinding CEREC Tessera™ CAD/CAM blocks using a Dentsply Sirona CAD/CAM system. Note that CEREC® Version 5.1.1 or higher and inLab® 20.0 or higher are required.

A matrix firing cycle with glaze will increase the flexural strength. The use of stains is optional; however, glaze is required. Matrix firing refers to a post-grinding firing cycle using a CEREC SpeedFire furnace or traditional porcelain furnace with the application of glaze.

##### 1.1 Composition

CEREC Tessera™ CAD/CAM blocks are a high strength advanced lithium disilicate material that is tooth-colored. A matrix firing with glaze will increase the flexural strength. The final strength is achieved after the CEREC Tessera™ CAD/CAM block matrix firing step with glaze (see section 3.3.1).

##### 1.2 Indications

CEREC Tessera™ Advanced Lithium Disilicate CAD/CAM blocks are an all-ceramic system for the creation of:

- Veneers
- Inlays
- Onlays
- Crowns in the anterior and posterior region

Advanced Lithium Disilicate (ALD) dental ceramic material type II, class 4a pursuant to ISO standard 6872 (CTE:  $9.8 \cdot 10^{-6}/K$  (25 – 500°C) and Tg = 582° C).

##### 1.3 Contraindications

- Bridges spanning more than three units
- Temporary restorations
- Parafunction (bruxism)
- Cantilever bridges
- Patients with a substantially reduced residual dentition
- Inlay bridges/Maryland bridges

##### 1.4 Compatible stains and glazes

The use of spray or paint on glaze is required with CEREC Tessera™ CAD/CAM blocks with a matrix firing cycle; however, the use of stains is optional. CEREC Tessera™ CAD/CAM blocks are compatible with Dentsply Sirona's Universal Stains and Glaze System and Dentsply Sirona Universal Spray Glaze or Indenco™ ("Indenco" is not a registered trademark of Dentsply Sirona) Spray Glaze. CEREC Tessera™ can be veneered with the Celtra® Ceram veneering ceramic (cut-back only).

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See Celtra® Ceram IFU for processing instructions.

#### 1.5 Compatible Cements

CEREC Tessera™ CAD/CAM block restorations are compatible with Universal/Self-Adhesive, Adhesive Resin Cement and Conventional Cement systems, including all Dentsply Sirona cement systems (available separately, see complete Instructions for Use). Full coverage crown restorations that have a retentive prep design are compatible with all cements, including all Dentsply Sirona self-adhesive resin cements, adhesive resin cements and conventional cements (available separately, see specific product Instructions for Use). Inlays, onlays, non-retentive crowns should be adhesively bonded using an Adhesive Resin Cement and bonding agent. Use of other cements or cement systems with CEREC Tessera™ CAD/CAM blocks are at the discretion and sole responsibility of the dental practitioner.

#### 2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other sections of these instruction for use.

**Safety alert symbol** This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Follow all safety messages that follow this symbol to avoid possible injury.

#### 2.1 Warnings

- If properly processed and used, adverse effects from this medical device are highly unlikely. However, reactions of the immune system (such as allergies) or localized paresthesia (such as an irritating taste or irritation of the oral mucosa) cannot be completely ruled out. In case of skin sensitization or rash, discontinue use and seek medical attention.
- CEREC Tessera™ CAD/CAM block restorations are not suitable for patients with clinical symptoms of parafunctional habits or bruxism (see Contraindications).
- Do not inhale dust particles during grinding. Wear suitable protective mask.
- **Do not use firing pastes or stains and glazes of other manufacturers as it may impact the performance of the material.**
- In patients with hypersensitivity to any of the ingredients, this medical device may not be used at all or only under the advice of the dentist or physician in charge.

#### 2.2 Precautions

- This product is intended to be used only as specifically outlined in these Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and is the sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
- Contamination of the preparation or margin area with saliva, blood, water, or hemostatic agents during adhesive cementation may lead to an adhesive failure. Ensure adequate isolation and tissue management techniques during adhesive cementation.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- CEREC Tessera™ CAD/CAM block restorations require adequate preparation reduction and restoration thickness. Insufficient wall thickness may lead to premature failure.
- **CEREC Tessera™ CAD/CAM block restorations must be glazed and matrix fired before insertion. Direct insertion without matrix firing and use of glaze may lead to failure.**
- Use only in well ventilated areas.
- Do not cement with provisional cements. Use with provisional/temporary cements could cause fracture of restorations.

#### 2.3 Interactions

CEREC Tessera™ CAD/CAM blocks are designed to be fabricated using a Dentsply Sirona CAD/CAM system. Grinding blocks using non-compatible CAD/CAM systems may lead to inadequate or unacceptable restorations.

#### 2.4 Adverse reactions

No adverse reactions have been reported for CEREC Tessera™ CAD/CAM blocks. Should you hear or receive information about any adverse effects, notify Dentsply Sirona.

#### 2.5. Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Store in a dry place, protect from moisture. Do not use after the expiration date.


### 3. STEP-BY-STEP INSTRUCTIONS

#### 3.1 Preparation


Proper reduction of the tooth during preparation is essential for maximizing the strength, shade and retention of the finished restoration. When preparing anterior or posterior teeth, the anatomical form must be reduced as shown.

**Minimum wall thickness:** The following diagram shows the specified minimum wall thickness for each indication. The minimum wall thickness must still be ensured after all manual adjustments have been made:

#### Minimum Wall Thickness (Conventional Cementation\*)

#### Minimum Wall Thickness (Adhesive Bonding)

#### 3.1.1 Inlays and onlays

- Preparation should be free of undercuts, exhibiting a draw and all internal line angles should be rounded. The preparation should be at least 1.0 mm deep in the central fossa. Ensure margins are placed clear of occlusal contact points.
- Inlays and onlays produced from CEREC Tessera™ CAD/CAM blocks should be adhesively bonded. See complete Instructions for Use of the adhesive and cement products.

#### 3.1.2 Full crowns

- Ensure that there is an axial reduction of 1.0 to 1.5 mm with the walls forming an angle of 4 – 8 degrees with the long axis of the tooth. In centric and dynamic occlusion, reduce incisally/occlusally by 1.5 mm. The lingual shoulders must be extended at least 1.0 mm into the proximal contacts surfaces. It is recommended to use a shoulder preparation without a bevel: All angles must be rounded, and the preparation surfaces must be smooth.
- Full crowns produced from CEREC Tessera™ CAD/CAM blocks can be cemented by adhesive or self-adhesive methods. See complete Instructions for Use supplied with the adhesive and cement products.

#### 3.1.3 Veneers

- The standard reduction is 0.6 mm for the labial surface and 0.4 mm in the gingival area (since the enamel is thinner in this region). Reduce the labiolingual incisal angle by 1.0 to 1.5 mm. The preparation margins should be located in enamel. A chamfer or rounded-shoulder preparation is recommended for all veneer margins. Proximal extensions must be located far enough proximally to conceal preparation margins from site and to avoid proximal gingival undercuts.
- Veneers produced from CEREC Tessera™ CAD/CAM blocks must be adhesively bonded. See complete Instructions for Use of the adhesive and cement products. Self-Adhesive cementation is not recommended for veneer restorations.

#### 3.2 CAD/CAM processing

CEREC Tessera™ CAD/CAM blocks for restorations are produced with CEREC or inLab CAD/CAM systems by Dentsply Sirona. If you have any questions about these systems, please contact Dentsply Sirona.

#### 3.2.1 Software requirements

- CEREC Tessera™ CAD/CAM blocks are supported by CEREC® version 5.1.1 and InLab software 20.1 and above.

#### 3.2.2 Grinding unit and block sizes

- In the CEREC® software, select the CEREC Tessera™ CAD/CAM Block material. The Dentsply Sirona CAD/CAM grinding unit will prompt you to insert a CEREC Tessera™ CAD/CAM Block. For detailed processing, please consult the Instructions for Use and technical manuals of the appropriate CAD/CAM systems. Make sure to follow the manufacturers' recommendations.

#### 3.3 Continue processing the grinded restoration

- CEREC Tessera™ CAD/CAM blocks are grinded, glazed, and matrix fired to increase strength. Instructions for staining and glazing can be found in section 3.3.1. Glaze is required.

#### 3.3.1 Instructions for staining and glazing of CEREC Tessera™ CAD/CAM block restorations

- Custom staining of CEREC Tessera™ CAD/CAM restorations can be performed using Dentsply Sirona Universal Stains and Glaze instructions below.
- For fast matrix firing in the CEREC SpeedFire, use Dentsply Sirona Universal Spray Glaze, Indenco™ Spray Glaze, or Dentsply Sirona Universal Overglaze paste paint-on glaze. If stains will be used with the fast matrix firing program, apply the stains and then apply the glaze over the stains and then matrix fire.**
- For best results, any tool marks on the surface should be removed with a suitable fine rubber wheel.
- Before applying the stains and glaze, the restoration must be clean and free of oil and other surface contaminants. Clean the surface of the restoration with a steam cleaner or in the ultrasonic cleaner with distilled water or a suitable cleaning fluid for 5 minutes. Any contamination after cleaning must be prevented. For disinfection, we recommend wipe disinfection with ethanol followed by immersion (1 min) in ethanol (70%).
- Before dispensing, mix the stains in the container thoroughly with a non-metallic spatula. Withdraw desired amount of Dentsply Sirona Universal stain and place it on a mixing palette. For a thinner consistency dilute the material with Dentsply Sirona Stain and Glaze Liquid.
- Before use, clean the brush with Dentsply Sirona Stain and Glaze Liquid. Do not use water to clean the brush since water can create an "orange peel" appearance.
- Apply a thin layer of stain where desired. Check the shade using the shade tab. Use the brush to increase or decrease the stain amount to achieve the desired shade.

#### Notes:

- A more intensive shade effect can be achieved by repeating cycles of applying and firing the material. However, the use of too many layers of stain may result in an unnatural appearance.
- The cusps and fissures can be individually characterized with stains.
- The basic shade is determined based on the Vita shade groups (A, B, C, D) (see Table 1). Please use the VITA™ classical shade guide.

**Table 1: Assignment of Shade Groups**

VITA™ Shade Group	Bleach	A	B	C	D
Universal Stain Shade	Stain 0	Stain 1	Stain 2	Stain 3	Stain 4

(\* VITA™ is not a registered Trademark of Dentsply Sirona)

**Table 2: Use of Incisal Stains**

Stain	Tooth Shade
Incisal Stain i1	A1, A2, A3, B1, B2, B3, B4, D2, D3, Bleach
Incisal Stain i2	A3.5, A4, C1, C2, C3, C4, D4

#### Use of glazes and stains from third-party manufacturers

Glazes and stains by third-party manufacturers may require other firing temperature and cycles than those recommended for CEREC Tessera™ CAD/CAM blocks. Use only Dentsply Sirona Universal Stain and Glaze and compatible spray glazes.

#### 3.4 Firing the restoration

The final strength is achieved after the matrix firing step with glaze in either a CEREC SpeedFire furnace or a traditional Porcelain furnace. For the traditional porcelain furnace, the CEREC Tessera™ specific programs (see Table 3) will have to be added. **Glaze is required.**

**NOTE: For traditional ceramic furnaces, the CEREC Tessera™ CAD/CAM block program will need to be added. DO NOT USE ANY OTHER FIRING PROGRAM SINCE DOING SO MAY LEAD TO FRACTURE AND/OR FAILURE.**

### 3.4.1 Speed Processing CEREC Tessera™ CAD/CAM Block Restorations in the Dentsply Sirona CEREC SpeedFire Furnace

- Make sure to select the CEREC Tessera™ CAD/CAM Block material under Dentsply Sirona in the CEREC Software prior to scanning as this will ensure the correct firing program will be transferred.
- Be sure to pre-heat the CEREC SpeedFire Furnace by touching the pre-heat button in the lower left section of the screen.

- Once the button has been pressed, the furnace will begin to pre-heat to approximately 400° C. The time to pre-heat to 400° C is about 1 minute.
- Apply Dentsply Universal Sirona Spray Glaze, Indenco™ Spray Glaze, or Dentsply Sirona Universal Overglaze paste paint-on glaze to the restoration. Apply an even amount of glaze to all sides of the restoration.
- Once the furnace is pre-heated to 400° C, push the job to be processed to open the furnace.

#### For Firing with firing pad (Honeycomb + firing pad):

- Place CEREC Tessera™ CAD/CAM restoration directly on the round firing pad; if using a square firing pad, it will need to be trimmed to fit on the round honeycomb.
- If firing an anterior, place the restoration with the lingual facing down and with premolars, place the restoration with the interproximal facing down. If there are any imperfections left on the glaze from the firing pad, this can be easily corrected by polishing the relevant area.

OR

#### For Pin Supported Firing (DS Investment Pins Only):

- Use only Dentsply Sirona Investment Pins (REF# 5365490111).
- Do not use other support devices or refractory putties or processing failure may occur.
- Do not use investment pins with Dentsply Sirona honeycomb or firing pads.
- Place the restoration on the appropriately sized investment pin.
- Place the investment pin in the center of CEREC SpeedFire furnace. Placement in the center is important to ensure the restoration is exposed to the correct temperature to obtain the optimum strength.
- Select the CEREC Tessera™ CAD/CAM job and press the "Start" button on the display screen.
- Please note that due to the hardware, the correct firing temperature of the CEREC SpeedFire is below the temperature listed in this document for conventional furnaces.
- If firing an anterior, place the restoration with the lingual facing down and with premolars, place the restoration with the interproximal facing down. If there are any imperfections left of the glaze from the firing pad, it can be easily corrected by polishing that area.
- **Important to Note that only one crown at a time can be fired in the CEREC SpeedFire with the fast matrix firing cycle.**

### 3.4.2 General firing recommendations in a traditional Porcelain Furnace

- The matrix firing with glaze will increase the flexural strength of the CEREC Tessera™ CAD/CAM block restoration; the application of glaze is required.
- Starting temperature 400°C.
- The holding time is 2.00 minutes.
- Place the restoration on either a firing pad and then on a honeycomb tray OR on a Dentsply Sirona Investment Pin, then place on the firing table of the furnace.
- If firing an anterior, place the restoration with the lingual facing down and with premolars, place the restoration with the interproximal facing down. If there are any imperfections left of the glaze from the firing pad, it can be easily corrected by polishing that area.
- Additional glaze firings may be performed at 760°C in order to accentuate the shade, correct it with glaze, or to increase the gloss. It is necessary to coat the entire surface with Dentsply Sirona Universal Glaze (available separately, see complete Instructions for Use) to obtain a uniform finish.
- Check the shade against a shade tab and adjust if necessary.

Table 3: Firing Charts

General Firing Recommendation								
	Start temperature	Drying	Closing	Pre-heating time	Heating rate	Final temperature	Holding time	VAC
	°C	min	min	min	°C/min	°C	min	min
Paint-On Glaze	400	2:00	2:00	2:00	55	760	2:00	-
2nd & Subsequent Glaze Firing - if needed	400	2:00	2:00	2:00	55	760	2:00	-
Spray Glaze	400	0:00	1:00	1:00	55	760	2:00	-
Stain and Spray Glaze	400	1:00	1:00	1:00	55	760	2:00	-

Programat CS2, EP 5000/5010							
	Standby Temperature B	Closing time S	Heating rate t ↑	Firing temperature T	Holding time H	Vacuum on/off	Long-term Cooling L
	°C	min	°C/min	°C	min	Vac 1 (°C) Vac 2 (°C)	°C
For general recommendation Paint-On Glaze	403	6:00	55	760	2:00	off	0
2nd & Subsequent Glaze Firing - if needed	403	6:00	55	760	2:00	off	0
Spray Glaze	403	2:00	55	760	2:00	off	0
Stain and Spray Glaze	403	3:00	55	760	2:00	off	0

**Note:** Please visit [www.dentsplysirona.com/CERECtessera](http://www.dentsplysirona.com/CERECtessera) for other furnace firing programs.

### 3.5 Try-In and Pre-Cementation Surface Preparation

#### 3.5.1 Preparation of the CEREC Tessera™ CAD/CAM Block restoration

Try-in the restoration for marginal and proximal fit. Take care that the patient does not close down on the restoration at this stage. Make any necessary adjustments.

- Clean the restoration with an ultrasonic or steam cleaner or with alcohol.
- Apply 5 % hydrofluoric acid etching gel (Available separately, see manufacturer's complete Instructions for Use) to the interior of the restoration only. Etching time 30 seconds.
- **CAUTION: Follow manufacturer's precautions. DO NOT ALLOW THE HF GEL TO CONTACT TEETH, TISSUE EYES OR MUCOSAL SURFACES.** Using an air-water spray, rinse the hydrofluoric acid from the restoration surface. Ensure that eye protection is employed during this process.
- Gently dry the restoration with air.
- Apply Calibra® Silane Coupling Agent to the dried intaglio surface of the restoration. Allow the silane to sit on the surface undisturbed for a minimum of 60 seconds. Using compressed air, remove excess silane from the surface. The surface will appear "wet" which is normal.

### 3.5.2 Cementation

CEREC Tessera™ restorations can be cemented using Dentsply Sirona Adhesive Resin Cement (Calibra® Ceram), Universal/Self-Adhesive Cement (Calibra® Universal) and Conventional Cement (Calibra® Bio). Cements from other manufacturers can be used following the Instructions for Use from the respective manufacturer.

#### Adhesive Cementation (Full Coverage Crowns, Bridges, Inlays, Onlays)

Crowns, bridges, inlays and onlays fabricated from CEREC Tessera™ blocks can be cemented adhesively in all clinical situations provided the preparation has at least 1.0 mm of occlusal reduction in the central fossa and all other preparation guidelines are followed. **For successful adhesive cementation, suitable measures (e.g. a cofferdam) must be taken to ensure that no moisture reaches the tooth.** Note that inlays and onlays should only be cemented using a fully adhesive protocol. Cement the CEREC Tessera™ restoration with Dentsply Sirona Calibra® Ceram Adhesive Resin Cement following the product Instructions for Use (IFU). Adhesive resin cement and adhesive from other manufacturers can be used following manufacturer's instructions. Self-adhesive or universal resin cement formulas may also be used to cement CEREC Tessera™ restorations.

#### Conventional Cementation (Full Coverage Crowns and Bridges)

Full coverage crowns and bridges fabricated from CEREC Tessera™ blocks can be conventionally cemented with resin-modified glass ionomer (RMGI) type cements provided that the preparation has at least 1.5 mm reduction in the central fossa and along the axial walls. Additionally, the preparation design should be retentive in design (parallel axial walls, at least 4 mm in height). Inlays and onlays should only be cemented using a **fully adhesive protocol**.

Cement the CEREC Tessera™ restoration with Dentsply Sirona Calibra® Bio Bioceramic Luting Cement following the product Instructions for Use. RMGI type cements from other manufacturers can be used following their respective instructions for use.

#### Cementation with Universal/Self-Adhesive Cements (Full Coverage Crowns and Bridges)

Full coverage crowns and bridges fabricated from CEREC Tessera™ Blocks can be cemented using Universal/Self-Adhesive cements provided the preparation has at least 1.5 mm reduction in the central fossa and all other preparation guidelines are followed.

Cement the CEREC Tessera™ restoration with Dentsply Sirona Calibra® Universal Cement following the product Instructions for Use (IFU). Universal type cements from other manufacturers can be used following their respective instructions for use.

#### Veneer Cementation

Veneers fabricated from CEREC Tessera™ blocks should be cemented using a fully adhesive protocol using Dentsply Sirona Calibra® Veneer Cement (shaded veneer cementation system) following the product Instructions for Use (IFU). Shaded veneer resin cement systems from other manufacturers may be used following their respective instructions for use.

## 4. HYGIENE

CAUTION

### Cross-contamination

Do not reuse single use products.  
Dispose of in accordance with local regulations.

### 4.1. Disposal

Dispose of in accordance with local regulations.

## 5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

### 5.1 Do not use after expiration date.

ISO standard uses: "YYYY-MM" or "YYYY-MM-DD"

### 5.2 The following numbers should be quoted in all correspondences:

- Reorder number
- Lot number
- Expiration date

**5.3** Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations

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DeguDent GmbH