Intraligamentäre Spritze

Syringe for intraligamental anaesthesia

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
Injection system for intraligamental anaesthesia, for use with a 1.8 ml cartridge and metric threaded needle that can be used in intraseptal anaesthesia. (This product is also available with imperial thread under item number 1956-IT.)
Intended for use only by dental professionals.

Intraligamental anaesthesia (ILA) is a special form of infiltration and terminal anaesthesia. ILA differs from terminal anaesthesia in that the injection is placed directly into the gingival sulcus of the affected tooth and not parallel to the tooth axis in the mucobuccal fold of the mucosa adjacent to the tooth to be anaesthetised. This allows selective anaesthesia of only the tooth to be treated, while with terminal anaesthesia at least 2-3 teeth are affected by the anaesthetic.

Instructions for use

1. Construction of the CerviJect syringe

The functional components and assembly of the CerviJect syringe are shown in the following diagram:

- | | |
|--|--|
| 1 barrel | 6 hand piece |
| 2 viewing chamber with cartridge | 7 dosing lever |
| 3 plastic sleeve | 8 resetting key |
| 4 bayonet catch | 9 dosing plunger |
| 5 threaded nozzle for injection needle | 10 disposable needle (not included in the pack contents) |

2. Preparation

To return the plunger (9) to the starting point, depress the resetting key (8) and hold the syringe in an upright position (bayonet catch up).

Place the plastic sleeve (3) in the barrel (1).

Insert a standard cartridge in the barrel. Connect the barrel (1) and hand piece (6) using the bayonet catch (4). Depress the dosing lever (7) until a slight resistance is felt. The dosing plunger is now in contact with the anaesthetic cartridge bung and presses the tip of the cartridge into the inner part of the nozzle.

Insert a disposable needle into the nozzle and screw the hub onto the thread. We recommend the use of the cartridge of disposable needles in sizes 0.3 x 12 mm or 0.3 x 16 mm with the corresponding short needle point.

WARNING: Always insert ampoule into ampoule holder first and then screw on the cannula attachment (needle). In reverse order, there is the risk the needle does not puncture the ampoule and bends out of shape.

Depress the dosing lever until the anaesthetic solution emerges from the needle point. The syringe is ready for use.

After loading and assembling the syringe according to directions, guide the needle along the neck of the tooth into the gingival sulcus with direct access to the periodontal space. The needle is generally inserted from the buccal/labial or palatal/lingual aspect. Make the first press for injecting the anaesthetic as soon as there is perceptible contact with the alveolar bone. There should be a marked resistance, which guarantees that the anaesthetic does not flow into the surrounding soft tissue but that there has been adequate intraosseous infiltration. An external sign is that the mucosa appears pale due to vasoconstriction.

3. Description of applications for ILA

- a. Indications for ILA
- Extraction of single teeth
 - Efficient, additive use in the case of inadequate elimination of pain or failure of anaesthesia
 - Patients with unstable circulation
 - Differential diagnosis of definite pulpal problems
 - Endodontics (especially intracanalicular anaesthesia)
 - Where there is haemorrhagic diathesis
- b. Contraindications for ILA
- Where there are acute infections
 - Where there is intolerance to the vasoconstrictor
 - All types of marginal periodontal diseases
 - Deep periodontal pockets if the tooth does not have to be extracted
- c. Advantages of ILA for the user
- Anaesthesia takes effect after only a short time
 - Only a very small amount needs to be injected
 - Targeted, selective anaesthesia of individual teeth
 - Elimination of the risk of soft tissue haematomas
 - Intercanal effect - e.g. during endodontic treatments
 - Acts in the case of overacidified tissue (inflammation)
- d. Advantages for the patient
- Quick, gentle elimination of pain
 - No feeling of numbness in adjacent anatomical structures (lips, tongue)
 - Only minimal intake of anaesthetic (patients at risk)
 - Anaesthesia only lasts for a short period
 - As the ILA syringe does not have the typical appearance of a cartridge syringe, its slimline, delicate design has a calming effect on the patient. This also makes it generally successful to use with children.

4. Care instructions

The CerviJect syringe is sold nonsterile. Complete a pre-disinfection cycle followed by sterilisation cycle before initial use.

The CerviJect syringe is made from stainless steel and can be easily cleaned, disinfected and sterilised together with other stainless steel instruments in the practice (QSA 313). Adhere strictly to the preparation recommendations of the respective manufacturer.

After each use, loosen the barrel and hand piece of the syringe and remove the cartridge, plastic sleeve and dosing plunger. Thoroughly dry the individual components after cleaning and disinfection. Then place the individual components in a suitable sterilisation medium (sterile packaging, container) and only reopen them shortly before the next use. The individual components should only be assembled when dry.

The CerviJect syringe must not be lubricated.

Only use cartridges with a plunger seal that does not have a hollow.

5. Warranty

This product has been checked for flawless quality prior to dispatch. In case of defects in workmanship and materials, we accept a **2-year warranty** from the date of purchase. This guarantee is void if your CerviJect has been improperly sterilized, poorly stored or incorrectly handled.

All information provided is to the best of our knowledge and is only intended as non-binding information.

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Carl Martin GmbH