VOCO Clip/Clip F

CE 0482

EN Instructions for use
MD EU Medical device

Product description:

Clip is a light-curing material for temporary filling and sealing. The elastic consistency assures an easy removal, even in case of undercuts, and finishing the cavity is not necessary. Thus the fitting of the temporary filling will be improved.

Clip is available in syringes and Caps for direct application into the cavity.

Clip F is available in syringes. **Clip F** contains fluoride and continuously releases fluoride ions.

Indications:

- All types of temporary fillings
- Temporary protection in inlay and onlay technique

Contraindications:

Clip contains acrylates, methacrylates and BHT. **Clip F** contains methacrylates, BHT and fluorides. **Clip/Clip F** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Clip/Clip F is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Clip/Clip F should only be used by a professionally trained dental practitioner.

Use:

Preparation, lining and impression to be carried out as usual. Rinse and dry the cavity. Bring material to room temperature before use.

Clip/Clip F syringe:

Apply **Clip/Clip F** into the cavity with a suitable instrument (round condenser, spatula, etc.).

Clip Caps:

Insert Caps into the opening of the dispenser (follow the instructions for use). Turn Caps into the required direction and remove the protective cap. Apply **Clip** directly into the cavity by slightly and evenly pressing the levers of the dispenser together.

Further handling:

Contour with slight excess. For use in the inlay and onlay technique, keep the approximate contact points free of resin with dental floss or a matrix before curing to maintain tooth mobility.

Cure for 20 - 40 s (depending on layer thickness) with a conventional polymerisation device. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices. As a rule, 5 mm deep layers are sufficiently cured in 40 s. Increased exposure times are required when the light source cannot be brought close enough to the material.

Check occlusion. Excess can be removed with a scalpel. Silicone points are suitable for grinding.

In case of more than two weeks' stay in the mouth, and in occlusal or marginal areas, reinforcement with composite is recommended. Apply a small quantity of a composite material onto the corresponding areas and polymerise after contouring together with **Clip/Clip F**.

The slight expansion of **Clip/Clip F** ensures a tight seal of the temporary restoration. When treating adjacent teeth simultaneously, remove excess material in approximal areas prior to polymerisation to maintain tooth mobility and contact points.

Removal of the temporary filling:

Insert a rigid instrument into the temporary filling and extract it.

Warnings, precautionary measures:

- Adhesion to dental resins or linings (e.g. light-curing materials) may occur. In such cases apply a suitable insulating agent (e.g. glycerine, vaseline).
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Clip/Clip F**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Clip/Clip F**.
- Avoid contact of unpolymerised paste with mucous membranes to avoid sensitivity reactions.
- Each **Clip** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage:

Storage at 4°C-23°C. Refrigeration is advised if the material is not used regularly. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Clip/Clip F** must be reported to VOCO GmbH and the responsible authority.

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VC 60 BB1284 E1 0821 99 © by VOCO