D.T. LIGHT-POST

Translucent quartz fiber

 **VDW**[®]

EN	INSTRUCTION FOR USE
DE	GERBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
NL	GEBRUIKSAANWIJZING
PT	INSTRUÇÕES DE UTILIZAÇÃO
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
NO	BRUKSANVISNING
SV	BRUKSANVISNING
FI	KÄYTTÖOHJE
HE	הוראות שימוש
HI	उपयोग करने के लिए निर्देश
TH	ข้อแนะนำในการใช้งาน
AR	تعليمات الاستخدام

D.T. LIGHT-POST® FOR DENTAL USE ONLY

DESCRIPTION

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D.T. LIGHT-POST® posts:

- D.T. LIGHT-POST® post is single use device. It must not be re-used to avoid the risk of contamination.
- D.T. LIGHT-POST® posts can be cleaned with alcohol. They can neither be sterilized nor disinfected in thermosinfectant.
- The post must be sterilized individually before insertion into the canal with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g EN ISO 11607-1),
 - ✓ Autoclave: type B complying with regulation standard in force (e.g EN 13060),
 - ✓ Sterilization temperature: 134 °C – sterilization time: 18 minutes,
 - ✓ 1 cycle only.
- The post must be immersed in isopropanol for 1 min then dried gently with air syringe for 15 s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- The use of rubber dam is recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.
- For profoundly flared or ovoid canals, the use of accessory post is indicated for best results.
- A D.T. LIGHT-POST® posts re-access kit is available if needed. Contact your RTD distributor.

D.T. LIGHT-POST® drills:

- D.T. LIGHT-POST® drills are supplied non-sterile. Disinfect and sterilize D.T. LIGHT-POST® drills before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermosinfectant.
- Inspect D.T. LIGHT-POST® drills visually before each use: change to new drills when they are damaged, corroded or dull.
- The UNIVERSAL drill may be used only to remove root canal filling partially, and the FINISHING drills may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the D.T. LIGHT-POST® drills stated is 15 uses, following this the dental surgeon is informed that the drill should be replaced.
- D.T. LIGHT-POST® drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,

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

Distributed by: VDW GmbH, Bayerwaldstraße 15, 81737 Munich

- EN** SUGGESTED POST SIZE FOR EACH TOOTH (GUIDELINES ONLY)
- DE** EMPFOHLENE STIFTGRÖSSE FÜR JEDEN ZAHN (NUR ALS RICHTWERT)
- ES** TAMAÑO DE POSTE SUGERIDO PARA CADA DIENTE (SOLO PAUTAS)
- FR** TAILLE DE TENON PROPOSÉE POUR CHAQUE DENT (DONNÉE À TITRE INDICATIF UNIQUEMENT)
- IT** MISURA DEL PERNO PROPOSTA PER OGNI DENTE (VALORI INDICATIVI)
- NL** AANBEVOLEN STIFTMAAT VOOR IEDER GEBITSELEMENT (ALLEEN BEDOELD ALS RICHTLIJN)
- PT** TAMANHO DE ESPIGÃO SUGERIDO PARA CADA DENTE (APENAS PARA ORIENTAÇÃO)

NOT1500VW00EUAREV01-2020-01 -14

D.T. LIGHT-POST® FOR DENTAL USE ONLY

DESCRIPTION

D.T. LIGHT-POST® devices consist of D.T. LIGHT-POST® posts and D.T. LIGHT-POST® drills.

D.T. LIGHT-POST® posts are radiopaque, translucent fiber post. D.T. LIGHT-POST® posts present a smooth surfaced post, apically tapered with a double taper, and a length of 20 mm. The radio-opacity level is around 200% Al equivalent (ISO 4049).

D.T. LIGHT-POST® drills (UNIVERSAL/FINISHING) are intended for mechanical shaping and cleaning of the root canal. They are made of stainless steel and the shape is chosen according to the form of the corresponding post. D.T. LIGHT-POST® drills exist in a similar range of diameters as the associated post, and they are identified by a color code.

D.T. LIGHT-POST drills are associated with:

- D.T. LIGHT-POST posts
- D.T. LIGHT-POST ILLUSION X-RO posts

	■ 0,5	■ 1	■ 2	■ 3
Ø apical tip	0.80	0.90	1.00	1.20
ØPost head	1.25	1.52	1.80	2.18
Double taper	.02/.04	.02/.06	.02/.08	.02/.10
Color coding	■ Black	■ Red	■ Yellow	■ Blue

COMPOSITION

D.T. LIGHT-POST® posts are fiber-reinforced polymer made of radiopaque quartz fiber (80%) and epoxy resin matrix (20%).

D.T. LIGHT-POST® drills are made of stainless steel.

INTENDED USE

D.T. LIGHT-POST® devices have been designed to support and secure the coronal restoration, in case of insufficient residual tooth structure (<4 mm).

D.T. LIGHT-POST® drills are intended to prepare the tooth for insertion of a post that will support and secure the coronal restoration, where insufficient residual tooth structure (<4 mm) exists.

These products should only be used in hospitals, clinics or dental office with qualified dental personnel (dentists).

INDICATIONS FOR USE

D.T. LIGHT-POST® devices are indicated in case of insufficient residual tooth structure (<4 mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

D.T. LIGHT-POST® devices are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of D.T. LIGHT-POST® devices for pregnant women, breastfeeding women and children under 18 years old (for posts only). It is not recommended for such use.

If the device is used in a way different than the one recommended in the DFU, loosening, breakage of the post or extraction of the tooth might occur.

D.T. LIGHT-POST® posts:

- D.T. LIGHT-POST® post is single use device. It must not be re-used to avoid the risk of contamination.
- D.T. LIGHT-POST® posts can be cleaned with alcohol. They can neither be sterilized nor disinfected in thermoisinfectant.
- The post must be sterilized individually before insertion into the canal with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g EN ISO 11607-1),
 - ✓ Autoclave: type B complying with regulation standard in force (e.g EN 13060),
 - ✓ Sterilization temperature: 134 °C – sterilization time: 18 minutes,
 - ✓ 1 cycle only.
- The post must be immersed in isopropanol for 1 min then dried gently with air syringe for 15 s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- The use of rubber dam is recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.
- For profoundly flared or ovoid canals, the use of accessory post is indicated for best results.
- A D.T. LIGHT-POST® posts re-access kit is available if needed. Contact your RTD distributor.

D.T. LIGHT-POST® drills:

- D.T. LIGHT-POST® drills are supplied non-sterile. Disinfect and sterilize D.T. LIGHT-POST® drills before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfectant.
- Inspect D.T. LIGHT-POST® drills visually before each use: change to new drills when they are damaged, corroded or dull.
- The UNIVERSAL drill may be used only to remove root canal filling partially, and the FINISHING drills may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the D.T. LIGHT-POST® drills stated is 15 uses, following this the dental surgeon is informed that the drill should be replaced.
- D.T. LIGHT-POST® drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,

- ✓ Autoclave: type B complying with regulation standard in force,
- ✓ Sterilization temperature: 134 °C – Sterilization: 18 minutes,
- ✓ Store the sterilized components in a dry, dust-free place,
- ✓ If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

ADVERSE REACTIONS

- Debonding of the post
- Root fracture
- Fracture of the post

STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
2. Select the correct **D.T. LIGHT-POST®** post size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the **D.T. LIGHT-POST® UNIVERSAL** drill selected (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the **D.T. LIGHT-POST® FINISHING** drill corresponding to the selected **D.T. LIGHT-POST®** post size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal
6. For **D.T. LIGHT-POST®** post, position the colored O-ring as a cutting mark and then remove post. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post. Remove the O-ring.
7. Immerse or soak the post in isopropanol for 1 minute and then gently air dry for 15 s at least.
8. Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.
9. Apply 2 coats* of the primer with a thin brush (e.g. Composibrush RTD) in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
- 11.a. If using a dual cure resin cement for both cementation and core build up such as **CORECEM™** (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic form, or directly on the tooth.
- 11.b. If using separate cement and a core build up resins: Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
- 11.c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
- 11.c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. RTD recommends that the coronal end of the post be covered by core composite.

STORAGE

Store at ambient temperature. Check the packaging of the product before using it. In case of damage, do not use the product.

DISPOSAL

Dispose in a safe manner in accordance with local/national regulations.

SAFETY AND LIABILITY

WARRANTY

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

LIMITATION OF LIABILITY

Except where prohibited by law, RTD will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only.

** Refer to the manufacturer's instructions for use.*

*** Please refer to your autoclave manufacturer for specific sterilization instructions if needed.*

TERMS AND SYMBOLS IDENTIFICATION

-  Manufacturer
-  Do not re-use
-  Batch code
-  Catalogue number
-  Consult Instructions for use
-  Caution
-  Number of unit / Content

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D.T. LIGHT-POST® posts:

- D.T. LIGHT-POST® post is single use device. It must not be re-used to avoid the risk of contamination.
- D.T. LIGHT-POST® posts can be cleaned with alcohol. They can neither be sterilized nor disinfected in thermosinfectant.
- The post must be sterilized individually before insertion into the canal with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g EN ISO 11607-1),
 - ✓ Autoclave: type B complying with regulation standard in force (e.g EN 13060),
 - ✓ Sterilization temperature: 134 °C – sterilization time: 18 minutes,
 - ✓ 1 cycle only.
- The post must be immersed in isopropanol for 1 min then dried gently with air syringe for 15 s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument,

such as wire cutters, as the pressure can destroy the structure of the post.

- The use of rubber dam is recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.
- For profoundly flared or ovoid canals, the use of accessory post is indicated for best results.
- A **D.T. LIGHT-POST®** posts re-access kit is available if needed. Contact your RTD distributor.

D.T. LIGHT-POST® drills:

- **D.T. LIGHT-POST®** drills are supplied non-sterile. Disinfect and sterilize **D.T. LIGHT-POST®** drills before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfectant.
- Inspect **D.T. LIGHT-POST®** drills visually before each use: change to new drills when they are damaged, corroded or dull.
- The UNIVERSAL drill may be used only to remove root canal filling partially, and the FINISHING drills may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the **D.T. LIGHT-POST®** drills stated is 15 uses, following this the dental surgeon is informed that the drill should be replaced.
- **D.T. LIGHT-POST®** drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,
 - ✓ Autoclave: type B complying with regulation standard in force,
 - ✓ Sterilization temperature: 134 °C – Sterilization: 18 minutes,
 - ✓ Store the sterilized components in a dry, dust-free place,
 - ✓ If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

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STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
2. Select the correct **D.T. LIGHT-POST®** post size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the **D.T. LIGHT-POST®** UNIVERSAL drill selected (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the **D.T. LIGHT-POST® FINISHING** drill corresponding to the selected **D.T. LIGHT-POST®** post size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal
6. For **D.T. LIGHT-POST®** post, position the colored O-ring as a cutting mark and then remove post. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post. Remove the O-ring.
7. Immerse or soak the post in isopropanol for 1 minute and then gently air dry for 15 s at least.
8. Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.
9. Apply 2 coats* of the primer with a thin brush (e.g. Composibrush RTD) in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
- 11.a. If using a dual cure resin cement for both cementation and core build up such as **CORECEM™** (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to

60 seconds. Continue applying the core build up using a plastic form, or directly on the tooth.

- 11.b. If using separate cement and a core build up resins: Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
- 11.c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
- 11.c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. RTD recommends that the coronal end of the post be covered by core composite.

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