



## Duceragold® Kiss

Directions for Use  
 Gebrauchsanweisung  
 Mode d'emploi  
 Istruzioni per l'uso  
 Instrucciones de uso

Instructions for use	1	en
Verarbeitungshinweise	6	de
Indications de mise en oeuvre	11	fr
Note per la lavorazione	16	it
Consejos para su utilización	21	es
تعليمات بخصوص التحضير	26	ar
Pokyny pro zpracová	31	cs
Bearbejdsesindikationer	36	da
Υποδείξεις επεξεργασίας	41	el
Napomene za obradu	46	hr
Feldolgozási útmutató	51	hu
加工について	56	ja
가공처리 주지사항	61	ko
Apdorojimo nuorodos	66	lt
Norādes sakarā ar pārstrādi	71	lv
Instructies voor de bewerking	76	nl
Informacje dotyczące obróbki	81	pl
Instruções de utilização	86	pt
Рекомендации по работе с материалом	91	ru
Beredningsanvisningar	96	sv
İşlem talimatları	101	tr
加工说明	106	zh

Release date: 2017-09

## Duceragold® Kiss

### Indications for use

Duceragold Kiss is indicated for preparation of Crowns and Bridges:

- Veneering metal framework and copings
- Veneering pressable ceramic

### Product Information

- Duceragold Kiss is a ceramic material suitable for veneering metal copings and frameworks for single crowns and bridges with a coefficient of thermal expansion (CTE) of 16.2 to 17.6  $\mu\text{m}/\text{m} \cdot \text{K}$  (25–600 °C) and for veneering the Cergo pressable framework class ceramic.

### Contraindications

- Suitable only for the indications listed above
- Duceragold Kiss is contraindicated in bruxism or other types of parafunction.
- In addition, Duceragold Kiss is contraindicated in situations where the interocclusal dimension is insufficient.

### Warnings

Undesired side effects of these medical products are extremely rarely to be expected with proper processing and application. Immunoreactions (e.g. allergies) and/or local unpleasant perceptions (e.g. taste irritations or oral mucosa irritations), however, cannot completely be excluded on principle. If you know any undesirable side effects – also in doubtful cases –, then please communicate them to us.

In case of hypersensitivity of the patients against Duceragold Kiss veneering porcelain or one of the ingredients, this medical product must not be used or only used under strict supervision by the treating physician/dentist.

Known cross-reactions or interactions of the medical product with other medical products or materials existing in the mouth must be taken into account by the physician/dentist when using the medical product.

Please pass on all the above-mentioned information to the treating physician/dentist, when processing this medical product for an individually manufactured object. Notify the dentist or physician in charge of all factors described above if you use this medical product for a custom construction.

- Do not inhale abrasive dusts
- Pastes and paste fluid are harmful if swallowed.

### Precautions

When working with these materials, make sure to comply with the Instructions for Use and the pertinent Material Safety Data Sheets (MSDS).

- For professional use only

### Adverse reactions

We are not aware of any risks and/or side effects related to Duceragold Kiss veneering porcelain.

### Technical data

- CTE dentine: 15,1  $\mu\text{m}/\text{m}\cdot\text{K}$  (25–500 °C)
- Dental ceramics, type 1, class 1 according to DIN EN ISO 6872
- Metal-ceramic bond, flexural strength and chemical solubility according to DIN EN ISO 9693 / 6872
- Use only with alloys with a solidus temperature of at least 900 °C.

### Alloy selection

- Duceragold Kiss is compatible with high-gold alloys (such as Degunorm and Degunorm supra/pur). Ask the alloy manufacturer about the alloy's composition and coefficient of thermal expansion. If the cooling times are observed, bonding alloys with a CTE of 16,2–17,6  $\mu\text{m}/\text{m}\cdot\text{K}$  (25–600 °C) can be recommended.

Market launch: May 2003

### Transport and storage conditions

- Protect powders and pastes from light and moisture and store them in a vibration-free place.

Keep dry

Keep away from sunlight

### Please pay attention to the following symbols on the product labels:

- REF Catalogue number
- LOT Batch code
-  Use by
-  Pay attention to the instructions for use
-  Do not reuse
-  Date of manufacture

### Combinable Liquids

- Powder Opaques:
  - Ducera® Liquid B
  - Ducera® Liquid OCL universal
- Paste Opaques:
  - Fluid Paste opaque
- Shoulder porcelain:
  - Ducera® Liquid Quick
- Dentines/Incisals, etc.:
  - Ducera® Liquid SD
  - Ducera® Liquid Blend
  - Ducera® Liquid Form
- Stain/Glaze:
  - Ducera® Liquid Stain
- Isolating Fluid:
  - Ducera® Sep Isolating Fluid

### Ceramic furnace

For best results, make sure that the required firing temperatures and times are observed strictly.

If necessary, adjust the firing parameters of the ceramic furnace as required.

## Names/Explanations of important materials

### Power Chroma (PC 1–6)

The Power Chromas are highly chromatic, fluorescent and intensive powders for individual shade adjustments. The power chromas can be used to increase the chroma in the cervical, palatal, and occlusal areas. The powders can be used pure or in a 1:1-mixture. By mixing the Power Chromas with Stand by you can create an excellent mamelon dentine. The shade conversion chart is meant as a guidance.

### Stand by

Strongly opalescent, almost transparent multi functional powder. It is possible to use Stand by pure or to mix it with all other powders of the Kiss concept. Hence the Stand by powder has a key function.

### Opal Effekt Sunrise/Opal Effekt Sunset

Opalescent effect powder for yellow as well orange/reddish incisal areas. Very well suited for increasing the chroma at the 2nd or 3rd dentine firings. Sunset is mainly used for A-shades, whereas Sunrise is used mainly for B-shades. Both powders can be diluted by using Stand by.

### Opal Effekt Sky/Opal Effekt Ocean

Opalescent effect powders for discreet and for strong, deep blue incisal areas – can be diluted by using Stand by.

### Opal Effekt Fog

Opalescent effect powder for grayish incisal areas – can be diluted by using Stand by.

### White Surface

Whitish opalescent effect powder to enhance occlusal cusps in the posterior region as well as for palatal/lingual ridges in the anterior area – can be diluted by using Stand by.

### Final Kiss

Low-fusing, transparent correction powder (Final Kiss) – firing temperature 720 °C.

Table 1: Colour conversion chart, Duceragold® Kiss

Shade	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
<b>Standard layering</b>																
Opaque	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
Dentine	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
Incisal	1	2	3	3	5	1	1	4	6	1	5	5	6	2	4	4
<b>Individual layering</b>																
Opaque	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
Opaque Orange	For characterizing occlusal, cervical and palatal regions.															
Opaque Bleach	For extremely bright/bleached teeth. Normally used in connection with bleach dentines only.															
Opaque Gum	For gingival aspects.															
Shoulder SM/F SM	1	2	2+3	2+4	3+4	1	1+3	3	3+5	1	1+4	2+4	4	1+4	2+4	3+4
Dentine	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
Power Chroma 1																
Power Chroma 2																
Power Chroma 3	1+2	2	2+5	3+5	4+6	1	1+3	2+3	3+6	1+6	2+6	3+6	5+6	1+6	2+6	3+6
Power Chroma 4																
Power Chroma 5																
Power Chroma 6																
Flu Inside 1	x	x	Mix	Mix		x	x	Mix		x	Mix			x	Mix	
Flu Inside 2					x				x			x	x			x
Opal incisal 1	x	x	Mix	Mix		x	x	Mix		Mix	Mix			Mix	Mix	
Opal incisal 2					x				x			x	x			x

The Kiss colour wheel shows you the right powders for the different shades at a glance.

### 1 2 3 4 Framework preparation

The framework should be designed in a small anatomic shape for supporting the porcelain and to ensure an even layering thickness of the veneering material.

Use only cross-cut tungsten-carbide cutters to finish the framework. The use of chip-removing tools is mandatory for this alloy range. Using stones invariably leads to microscopic overlaps, which in turn will lead to bubble formation in the ceramic veneering material.

To avoid tensions in the Duceragold Kiss veneering porcelain, it is necessary to round off any sharp corners and edges of the frame.

Carefully sandblast the framework after finishing and correcting using 110-µm alumina at a recommended pressure of 2 bars.

Then the frameworks are cleaned by steam or in a clean ultrasonic bath (see fig. 1–4).

### Soldering and laser-welding

For soldering and laser-welding metal frameworks, follow the comprehensive instructions and notes of the alloy manufacturer.

### 5 6 Applying the opaque

The Duceragold Kiss paste opaque cover the framework in a layer of uniform thickness in the basic colour selected. It is recommended to perform two opaque firings.

When using alloys containing copper (such as Degunorm), it is necessary to use the Degunorm Classic base paste instead of the standard opaque before the first opaque firing. The Degunorm Classic base paste reliably covers the dark copper oxide, effectively preventing any discolouration of the veneering ceramic. The Duceragold Kiss paste opaque and Classic base paste can be applied easily with a standard paste opaque brush, without diluting. If the paste changes its consistency over time and becomes less easy to apply, the original consistency may be restored by adding a small quantity of paste opaque fluid. Alternatively, the framework may be coated with a very thin paste opaque fluid layer (see fig. 5–6).

Caution: Use only very small quantities of paste opaque fluid. Using too much paste opaque fluid may result in cracks and bubbles during opaque firing.

If dilution becomes necessary, please use the associated paste opaquer fluid. If necessary, please dilute outside the little glass. Moreover, please take note of the longer pre-drying time (now 2 min. instead of 1 min.).

### 7 8 9 Aesthetic Line Basic

In the standard layering technique (opaque/dentine/incisal), highly aesthetic veneers can be created within a minimum of time. Start by building up the dentine core and then use the cut-back technique to obtain room for the incisal layer. Rebuild the cut-back portion of the crown with incisal material. This is followed by the first dentine firing (see table 2). Complete to the original contour using the appropriate dentines and incisals. This is followed by the second dentine firing (see table 2) and the finishing of the restoration. Finally, this is followed by the glaze firing (see table 2), which may be performed with or without glaze or Duceragold Kiss stains for special tooth characteristics (see fig. 7–9).

### 10 11 12 Ceramic shoulder

- When creating a ceramic shoulder, use a carbide bur to reduce the crown margin so it ends 0,5–0,8 mm above the lowest part of the chamfer or ledge. Also, ensure that the crown margin tapers off softly to avoid tension in the ceramic material.
- Sandblast the framework inside (especially the margin) and outside carefully, as described initially, and afterwards clean the framework (steam-cleaner).
- Mark the preparation border with a graphite-free pencil and seal die as usual, e.g. with superglue.
- Apply Ducera-Sep ceramic insulation generously. Expose the insulation to the air and apply a second layer. The insulation must be dry before further processing.
- Mix the shoulder margin material (SM 1–5) with the SD Quick modelling liquid. Afterwards complete the margin area by using the shoulder margin porcelain. Apply the margin porcelain in the cervical quarter. Let the material dry, maybe by using a heat source, and lift the coping carefully from the die. Bake the shoulder margin after the recommendation.
- The gap created by the shrinkage of the ceramic can be complemented by a second firing or after the glaze by using the final shoulder margin materials (F-SM 1–5) (see fig. 10–12).  
With Degunorm logic do not fire over 760 °C.

**13 14 15 16 17 18 Aesthetic Line Individual**

With the individual layering technique, the power chroma and opal effect powders allow you to create natural-looking and highly sophisticated restorations. With the Power Chromas, you can first line the framework with a highly chromatic, fluorescent basic layer for individual tooth shade characteristics. This

is followed by building up dentine layer, as usual. After the cutback, complete the incisal area with the corresponding opal effect enamel and use e.g. Opal Effect Sky or Ocean (in the incisal areas) or OE Sunrise or Sunset (for the cervical and body areas) as supplements (see fig. 13–18).

Table 2: General firing recommendations – Duceragold® Kiss

	Pre-heating temperature °C	Drying time min	Degunorm		Holding time min	Vacuum hPa	Tempering
			Heating rate °C/min	Firing temperature °C			
Oxide firing	575	0:00	55	780	5:00	50	–
Degunorm classic Base	575	7:00	55	780	1:00	50	–
Paste opaque	500	13:00	55	780	2:00	50	–
Margin 1	450	5:00	55	780	1:00	50	–
Margin 2	450	5:00	55	780	1:00	50	–
Dentine 1	450	6:00	55	780	1:00	50	3 min/720°C
Dentine 2	450	5:00	55	780	1:00	50	–
Glaze	450	3:00	55	770	1:00	–	–
Correction (Final Kiss)	450	2:00	55	720	1:00	50	50
Final margin	450	2:00	55	720	1:00	50	–

Generally, a tempering phase is only necessary for the 1st dentine firing. All other firings take place without cooling phase. With Degunorm eco, Econolloy and Econolloy Au, all firings starting with the 1st dentine firing are carried out with a tempering phase. With Degunorm logic, all firings are carried out without a tempering phase.

The values indicated here are recommended values and serve only as a clue. Deviations of the burning results are possible. The burning results depend on the respective output per oven and are due to the manufacturer and to age. Therefore the recommended values must be individually adapted at each burning. We recommend a test burning to control the oven. All indications have been carefully elaborated and tested by us, but they are passed on without any guarantee.

Caution: Note the following procedure for pre-drying try-in objects made of Duceragold Kiss hydrothermal ceramic material:

- Heat objects to about 80–90°C and allow to pre-dry for 30 minutes, or 1 to 1.5 hours for larger objects, until the moisture and organic residue accrued during try-in have dried.
- Depending on the degree of contamination, steam down or sandblast lightly with 50 µm alumina to remove any organic residue.
- Heat the restoration to 450°C at a heating rate of 5–10°C per minute and hold for another hour.
- Conclude by performing any corrections required.

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DeguDent GmbH
Dentsply International Inc.