# Dyract<sup>®</sup> flow

## Flowable Compomer Restorative

Compomère fluide de restauration

Compomero da restauro fluido

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CAUTION: For dental use only.

USA: Rx only.

Content	Page
1 <b>Product description</b> .....	2
2 <b>Safety notes</b> .....	3
3 <b>Step-by-step instructions</b> .....	5
4 <b>Hygiene</b> .....	8
5 <b>Lot number, expiration date and correspondence</b> .....	9

## 1 Product description

**Dyract<sup>®</sup> flow** restorative is a compomer restorative material with flow characteristics which make it ideal for small cavities in anterior and posterior teeth and for cavity lining. **Dyract<sup>®</sup> flow** restorative adapts to the cavity walls without the use of hand instruments.

**Dyract<sup>®</sup> flow** restorative is available in 5 VITA<sup>®1</sup> shades in a regular opacity (A2, A3, A4, B1) and an opaque (O-A3) shade.

### 1.1 Indications

- Restoration of minimally invasive cavity preparations, pits and fissures, including air abrasion and tunnel preparations.
- Small restorations without antagonist contact in anterior and posterior teeth.
- Restoration of minimal, shallow Class V preparations, including incipient Class V abfraction and erosion lesions.
- Base/liner under composite, high-density composite and compomer Class I & Class II restorations.
- Undercut blockout and filling of small defects in indirect restoration preparations.
- Small margin defect repair in otherwise sound restorations.

### 1.2 Contraindications

- Use with patients who have a history of severe allergic reaction to methacrylate resins.

### 1.3 Delivery forms

- Small syringes containing 1.8 g (1 ml) restorative to be used in combination with disposable steel cannulas (applicator needles, dispensing tips) for direct intraoral application.

<sup>1</sup> Not a registered trademark of Dentsply Sirona, Inc.

## 1.4 Composition

- Strontium-alumino-fluoro-phosphor-silicate glass
- Diethylene glycol dimethacrylate (DGDMA)
- Ammonium salt of phosphoric acid modified methacrylate resin
- Carboxylic acid modified methacrylate resin
- Highly dispersed silicon dioxide
- UV stabilizer
- Ethyl-4-dimethylaminobenzoate
- Camphorquinone
- Butylated hydroxy toluene (BHT)
- Iron oxide pigments
- Titanium dioxide

## 1.5 Compatible adhesives

The material is to be used following application of a compatible (meth)acrylate-based dentin/enamel adhesive such as Dentsply Sirona adhesives designed for use with visible light cured composite restoratives (e.g., Prime&Bond® etch and rinse adhesives or Xeno®, self etching adhesives).

## 2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

### Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

## 2.1 Warnings

The material contains methacrylates and polymerizable monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

## 2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- If refrigerated, allow material to reach room temperature prior to use.
- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Only use the cannulas supplied by Dentsply Sirona for placement of the material.
- Cannulas are used for direct intraoral application of restorative. Cannulas are for single use only, discard in order to prevent cross-contamination. Cannula may clog as restorative may set inside.
- Tightly close syringes immediately after use. Do not allow cannula to remain as a syringe cap. Replace original cap tightly after each use.
- Always make sure cannula is properly connected to Luer-Lock before application.
- The material should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may cause material and not fully seated cannulas to eject from the syringe.
- The syringe cannot be reprocessed. To prevent the syringe from exposure to spatula or spray of body fluids or contaminated hands it is mandatory that the syringe is handled with clean/disinfected gloves.
- As additional precautionary measure, the syringe may be protected from gross debris but not from all contamination by applying the single use Disposa-Shield<sup>®</sup> protective barrier<sup>2</sup>.
- Disposa-Shield<sup>®</sup> protective barrier and cannulas are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Interactions:
  - Do not use eugenol- or hydrogen peroxide-containing materials in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
  - The restorative is a light cured material. Therefore, it should be protected from ambient light. Proceed immediately once the material has been placed.
  - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.

## 2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

---

<sup>2</sup> Disposa-Shield<sup>®</sup> Infection Control Barrier for syringe-type delivery systems, Reorder A88013S1.

## 2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store in sealed packs at temperatures between 10 °C and 24 °C (50 °F and 75 °F).
- Keep out of direct sunlight and protect from moisture.
- Do not freeze.
- Do not use after expiration date.

Humidity can adversely affect the properties of unsealed syringes. For unsealed syringes the following is recommended:

- Store in a dry environment (< 70% relative humidity).
- Use unsealed syringes within five months.

## 3 Step-by-step instructions

### 3.1 Shade selection

Shade selection is made whilst the teeth are hydrated and prior to the restorative procedure.

1. Remove any extraneous plaque or surface stain.
2. Use the VITAPAN<sup>®3</sup> classical shade guide. Use the central part of the respective VITA<sup>®</sup> tooth for shade selection.

### 3.2 Cavity preparation

1. Clean the tooth surface with a fluoride free prophylaxis paste (e.g., Nupro<sup>®</sup> prophylaxis paste).
2. Prepare cavity (unless not necessary, e.g., cervical lesion).
3. Use adequate isolation such as rubber dam.
4. Rinse surface with water spray and carefully remove rinsing water. Do not desiccate the tooth structure.

#### 3.2.1 Matrix placement and wedging (Class II restorations)

Place a matrix (e.g., Palodent<sup>®</sup> Matrix System or AutoMatrix<sup>®</sup> Matrix System) and wedge. Burnishing of the matrix band will improve contact and contour.

### 3.3 Pulp protection, tooth conditioning/dentin pre-treatment, adhesive application

Refer to adhesive manufacturer's Instructions for Use for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of the material.

### 3.4 Application

#### **Clogging.**

Separation of cannula from syringe.

1. Do not allow material to dry inside.
2. Do not reuse steel cannulas.
3. Always make sure cannula is properly connected to Luer-Lock before application. The material should extrude easily.
4. Do not use syringe, once removed from pouch after five months.

<sup>3</sup> Not a registered trademark of Dentsply Sirona, Inc.

### **Avoid contamination of the syringe.**

To prevent the syringe from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringe is handled offside the dental unit with clean/disinfected gloves.

Apply material in increments immediately after the application of the adhesive. The material is designed to be cured in increments up to a 2 mm depth/thickness.

1. Remove cap from the end of the syringe.
2. To assure free flow of material from syringe, express a small amount onto a mixing pad.
3. As an additional precautionary measure, the syringe may be protected from gross contamination but not from all contamination by applying the single use Disposa-Shield® protective barrier.
4. Attach disposable steel cannula<sup>4</sup> over the protective barrier to the end of the syringe.
5. Turn cannula clockwise until fully seated, then tug on cannula to assure locking into syringe collar.
6. Dispensing tip may be bent for easy access. Gentle curving of the tip is preferable to an acute bend. Bending may occlude the cannula. Be sure to check that this has not occurred by extruding a small amount through the tip.
7. Material should flow freely with GENTLE PRESSURE. DO NOT USE EXCESSIVE FORCE. If not, remove syringe from patient field and check for obstruction.
8. Apply material through the cannula to surfaces to be restored.

### **Excessive force.**

Injury.

1. Apply slow and steady pressure on the syringe.
2. Always make sure cannula is properly connected to Luer-Lock before application. The material should extrude easily.

#### **3.4.1 Instructions for application as cavity liner**

For most cavity lining applications, a thickness of 1 mm or less is recommended. With the exception of Class II proximal box (gingival) margins, cavity margins/conditioned enamel are usually closed with selected restorative material. Lining is generally confined to internal conditioned dentin surfaces. Avoid lifting the tip while dispensing to minimize air entrapment. Any visible air bubbles should be pierced with a clean, sharp explorer prior to curing. Technique Tip: Use of a clean instrument such as a sharp explorer may be used to evenly spread a thin layer of dispensed Dyract® flow restorative material over prepared cavity (dentin) surfaces prior to curing.

Proceed immediately to complete the restoration per restorative material manufacturer's Instructions for Use. For adhesive restorations, to not disturb or allow contamination of cured liner surface. If contamination occurs, follow adhesive manufacturer's Instructions for Use for re-application of adhesive.

Technique Tip: If surfaces remain undisturbed, additional application of adhesive between increments is neither necessary nor recommended.

<sup>4</sup> Use the disposable steel cannulas included in the Dyract® flow package. Dentsply Caulk Reorder No. of the disposable steel cannulas 685606; DeTrey Reorder No. 606.04.490

### 3.5 Light curing

1. Light cure each area of the restoration surface with a suitable curing light designed to cure materials containing camphorquinone (CQ) initiator, i.e., spectral output containing 470 nm. Minimum light irradiance must be at least 550 mW/cm<sup>2</sup> exposure time according to the shade as listed in the table below. Refer to curing light manufacturer's Instructions for Use for compatibility and curing recommendations.

#### Inadequate curing.

Inadequate polymerization.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check minimum irradiance.
4. Cure each area of each increment for the recommended curing time.
5. Check distance to surface to be cured.

#### Curing light irradiance $\geq 800$ mW/cm<sup>2</sup>

Shades	Increment/layer	Curing time
A2, B1	2 mm	10 sec
A3, A4	2 mm	20 sec
O-A3	2 mm	30 sec

#### Curing light irradiance 550-800 mW/cm<sup>2</sup>

Shades	Increment/layer	Curing time
A2, A3, B1	2 mm	20 sec
A4, O-A3	2 mm	40 sec

### 3.6 Finishing and polishing

1. Contour the restoration using finishing burs or diamonds.
2. Use Enhance<sup>®</sup> finishing and polishing instruments and interproximal finishing strips.

For contouring, finishing, and/or polishing, follow the manufacturer's Instructions for Use.

## 4 Hygiene

### Cross-contamination.

Infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.
- Reprocess reusable products according to the instructions.

### 4.1 Syringes – cross-contamination

#### Cross-contamination.

Infection.

- The syringes cannot be reprocessed.
- To prevent the syringes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringes are handled with clean/disinfected gloves. Do not reuse syringes if contaminated.
- Dispose of contaminated syringes in accordance with local regulations.

To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination.

Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard composite material that has been in contact with any fluid or non-sterile instrument.

Repeated liquid contact may damage label. Dry the syringe with a lint-free single-use cloth.

**NOTE: Vigorous wiping can destroy the label.** Wipe syringe gently.

### 4.2 Disposa-Shield® protective barrier

The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination.

1. Right after use, carefully remove protective barrier together with cannula without contaminating the device.

#### Cross-contamination.

Infection.

1. Protective barriers are intended for single use only. Do not reuse in order to prevent cross-contamination.
2. Discard in accordance with local regulations.

#### 4.3 Steel cannulas – single use

1. Right after use, carefully remove protective barrier together with cannula without contaminating the device.

##### **Cross-contamination.**

Infection.

1. Do not clean, disinfect or reuse steel cannulas.
2. Properly dispose the used steel cannulas in accordance with local regulations.

##### **Clogging.**

Separation of cannula from syringe.

1. Do not allow material to dry inside.
2. Do not reuse steel cannulas.
3. Always make sure cannula is properly connected to Luer-Lock before application. The material should extrude easily.
4. Do not use syringe, once removed from pouch after five months.

#### 4.4 Disposal

Dispose of in accordance with local regulations.

#### 5 Lot number ( ), expiration date ( ) and correspondence

1. Do not use after expiration date.  
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
  - Reorder number
  - Lot number
  - Expiration date
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

© Dentsply Sirona 2019-10-11

[These Instructions for Use are based on Master Version 04]

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Dentsply DeTrey GmbH