EXA'lence™

VINYL POLYETHER SILICONE (VPES™) IMPRESSION MATERIAL

Table with 2 columns: Property and Value. Rows include Viscosity, Extra Light, Light Body, Medium Body, Heavy Body, and Heavy Body Rigid.

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Table with 2 columns: Material and Application. Rows include Extra Light, Light Body, Medium Body, Heavy Body, and Heavy Body Rigid.

EXA'lence™

MATERIALE D'EMPREINTE VINYL POLYETHER SILICONE (VPES™)

Table with 2 columns: Viscosité and Type. Rows include Extra Light, Light Body, Medium Body, Heavy Body, and Heavy Body Rigid.

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
Table with 3 columns: Indication, Impression Technique, and Recommended Viscosities. Rows include Impression for facial, onlay, crown, and Impression for partial denture.

Table with 2 columns: Physical Properties (Average) and Storage Conditions. Rows include Working time, Elasticity, and Compression.

CONTRAINDICATIONS
Avoid use of this product in patients with known silicone and flavoring allergies.

Table with 2 columns: Precautions and Directions for Use. Rows include Cartridge Loading and Dispensing steps.

PRECAUTIONS
When extruding EXA'lence, care should be taken to avoid mixing or contact with the following materials.


DIRECTIONS FOR USE
CARTRIDGE LOADING AND DISPENSING
1. Lift the release lever of GC CARTRIDGE DISPENSER II...

Table with 3 columns: Indication, Application, and Recommended Viscosities. Rows include Abformung für Veneers and Abformung für partielle Prothesen.

Table with 2 columns: Physical Properties (Durchschnittswerte) and Storage Conditions. Rows include Working time, Elasticity, and Compression.

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
Table with 2 columns: Physical Properties (Valeurs moyennes) and Storage Conditions. Rows include Working time, Elasticity, and Compression.

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
Table with 2 columns: Physical Properties (Proprietà fisiche (Promedio)) and Storage Conditions. Rows include Working time, Elasticity, and Compression.

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
Table with 2 columns: Physical Properties (Propiedades físicas (Promedio)) and Storage Conditions. Rows include Working time, Elasticity, and Compression.

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DISINFECTING
The impression can be disinfected with a glutaraldehyde solution or other appropriate disinfectants.

MODEL FABRICATION
After the impression is rinsed with water, wait for 60 minutes before pouring a model.

CAUTION
Avoid contact of base, catalyst, mixed material or adhesive with eyes.

STORAGE
Recommended for optimal performance, store at room temperature (59 - 77°F; 15 - 25°C).

PACKAGES
EXA'lence EXTRA LIGHT BODY, Regular and Fast Set, 59g (48mL) x 2 cartridges.

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