D

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Warning notices: The product must be cleaned, disinfected and sterilised before first use as well as after any further use, according to the enclosed instructions.

Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

Manufacturer: HÄGER & WERKEN GmbH & Co. KG
 Ackerstraße 1
 47269 Duisburg, Germany
 T +49 (203) 99269-0
 F +49 (203) 299283
 www.hagerwerken.de

Product: Fit Rubberdam Frame

REF 355 319; 355 323 – 355 324

Warning notices: No special notices

Limitation of reprocessing:

Frequent reprocessing has only a small impact on this instrument. The instrument's shelf-life is mainly determined by wear and damage by use.

Place of handling: Clean surface with a disposable cloth or paper towel.

Storage and transport: No special requirements.

Preparations for cleaning: No special requirements.

Cleaning: Use cleaning agent suitable for stainless steel according to the manufacturer's instructions. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C in mechanical cleaning.

Disinfection: Use disinfecting solution according to the manufacturer's instructions of. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C in mechanical disinfection.

Maintenance: No special requirements.

Packaging: Standardized packing material for sterilisation can be used (recommendation).

Sterilisation: Vapor sterilisation at 121 °C for 15 minutes
 Vapor sterilisation at 134 °C for 5 minutes.

Control / Functional check: Sight check on damages, wear, deformation.

Storage: No special requirements.

MEDICAL DEVICE

The above-mentioned instructions have been validated as SUITABLE for preparation of a Medical Device and its reuse by the Medical Device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and material applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.