Use only carefully read the instructions for use.

EN

GC Initial™ CAST NP Cobalt-Chrome (CoCr) based casting alloy for metal ceramic veneering (type 4)

GC Initial CAST NP is a cobalt chrome based casting alloy for fabricating crown & bridges with a suitable metal veneering ceramic or acrylic veneering material. Thanks to its thermal expansion of 14,1 x 10⁻⁶ K⁻¹ it is suitable for the most common metal veneering ceramics on the market. It has outstanding properties like an optimal casting performance, an extra fine molecular structure resulting in smooth, compact surfaces offering an easy and straight-forward processing. The alloy offers a reduced oxide formation even after repeated ceramic firings.

GC Initial CAST NP is produced conform to **ISO9693-1:2012** and **ISO 22674:2006** standards. It has passed the cytotoxicity test according to the norm ISO 10993-5. The alloy is highly corrosion resistant and it is free of beryllium, indium and gallium.

COMPOSITION

Cobalt	Co	63%
Chromium	Cr	28%
Tungsten	W	3%
Niobium	Nb	4%
Silicium	Si	1%
Others	Fe / Mo	1%

ALLOY CHARACTERISTICS
Solidus-liquidus temperature: 1253 - 1304°C
Melting point: 1450 °C
Highest firing temperature: 950°C
Thermal expansion coefficient: (25 - 500°C) 14,1 x 10⁻⁶ K⁻¹ / (25 - 600°C) 14,5 x 10⁻⁶ K⁻¹ / (500°C) 8,4 g/cm³
Vickers hardness: 324 HV10
Percentage elongation at fracture: 3,4%
Yield strength: 475 MPa
Modulus of elasticity: 194 GPa
Colour: White
Recommended ceramics: GC Initial MC

INSTRUCTIONS FOR USE

I. WAX-UP
The wax should be 0,5 mm thick to ensure that after the manufacturing process the metal wall thickness is at least 0,3 mm.
Avoid sharp edges and undercuts. Frameworks need to be anatomically reduced ensuring a uniform thickness of the ceramic layer. Connectors should be modelled as strong and as high as possible.

Single crowns: can be directly sprued. Suggested round casting sprues: Ø 3,0 mm to 4,0 mm.
Bridges: should be cast on a crossbar spanning the crown & bridge. A crossbar spanning technique is advised. Connections to crossbar Ø 3 to 3,5 mm.
Distance from modulation to crossbar 2,5mm. Thickness of the crossbars Ø 3 to 5 mm. Thickness of sprues to crossbar Ø 3 to 4 mm.

The higher the volume of the casting the thicker the sprue should be. The casting sprue depends on the object dimensions.
II. INVESTING
Use only phosphate bonded investment material for Crown & Bridges. GC Fijvestium and GC Fijvestil II are optimized for producing high precision fit castings in combination with smooth cast surfaces. Guidelines for optimal use of GC's phosphate bonded investments for Crown & Bridge can be found in the "GC Focus Edition – Phosphate Bonded Investments for Crown & Bridge" manual. Pre-heating temperature of 850-950°C is advised. Follow the manufacturer's Instruction for Use related to the pre-heating times before casting.

III. CASTING

We recommend to only use fresh alloys in own crucibles. Sprues and cones can only be re-used once, given that new metal coming from the same lot is added in equal quantity.
GC Initial CAST NP in pre-heated crucibles used for this alloy only.
Amount of casting alloy to use = wax framework weight with sprues / wax density x metal density.

High frequency melting equipment

When lost ingots have melted in the pre-heated ceramic crucible and the incandescence shadow has disappeared, immediately start the casting process. Do not use any flux.

Open flame melting

Place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux.

Overheating the melted material could cause shrinking hollows, micro-porosity and building of a coarse-grained structure and cause bridges to break or ceramic veneering to crack.

Follow the Instructions of the respective manufacturer of the casting devices for parameters and casting procedures.

After casting bench-cool the casting ring. Remove the investment material with care and sandblast using a non-recycling sandblaster with 100 to 125 µm aluminium oxide at a pressure of 3-4 bar.

IV. PREPARATION FOR VENEERING
Use tungsten burrs to grind and finish the prepared frameworks. Use always the same grinding instruments to avoid contamination.

Use diamond and boron burrs. Grind the surface with slow speed and low pressure. Grind the surface in one direction only to avoid overlying layers (ingot to air bubbles after the first opaque firing). Sandblast the prepared framework using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a maximum pressure of 3-4 bar and then steam clean.

V. OXIDE FIRING

Fire for 5 min. under vacuum at 950-980°C (10°C less than oxide melting firing temperature). After firing, carefully sandblast the oxide layer using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a pressure of 2,5-3 bar, then steam clean. The framework should have a homogeneous grey surface. Sandblast with fine particles using a compatible veneering ceramic following the manufacturer's instruction.

VI. VENEERING

- Veneering with GC Initial MC: please refer to GC Initial CAST NP manual
- Veneering with other ceramics: please refer to the manufacturer's instruction for use
- Veneering with acrylics: please refer to the manufacturer's instruction for use

PACKAGING
GC Initial™ CAST NP 250g
GC Initial™ CAST NP 500g
GC Initial™ CAST NP 1000g

CAUTION

- Cobalt-Chrome based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Before prosthetic application verify if other metal implants in an patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- Metal dusts and smoke are dangerous for health. Use exhaust fans while casting and sandblasting.
- We recommend to keep batch records and link it to patient's file to allowing its complete traceability.
- This product does not need any special preservation precautions. Nonetheless it has to be employed within the expiry date on the packaging
- Dispose of all waste according to local regulations.

GUARANTEE
All recommendations are based on GC findings and GC Initial tests and are therefore only to be used as a general safety claim for compensation GC tests. In any event, shall not exceed the value of the individual order of GC Initial CAST NP.

Undesired effects - Reporting

If you become aware of a kind of undesired effect, reaction or similar events experienced by you of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://eur.europa.eu/growth/sectors/medical-devices/contacts_en or to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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Use the following instructions for use.

(EN)

GC Initial™ CAST NP Cobalt-Chrome (CoCr) based casting alloy for metal ceramic veneering (Type 4)

GC Initial™ CAST NP is a cobalt chrome based casting alloy for fabricating crown & bridges with a suitable metal veneering ceramic or acrylic veneering material. Thanks to its thermal expansion of 14,1 x 10⁻⁶ K⁻¹ it is suitable for the most common metal veneering ceramics on the market. It has outstanding properties like an optimal casting performance, an extra fine molecular structure resulting in smooth, compact surfaces offering an easy and straightforward processing. The alloy offers a reduced oxide formation even after repeated ceramic firings.

GC Initial™ CAST NP is produced conform the **ISO9693-1:2012** and ISO 22674:2006 standard. Es hat den Zytotoxizitätstest according to the norm ISO 10993-5. The alloy is highly corrosion resistant and it is free of beryllium, indium and gallium.

Cobalt	Co	63%
Chromium	Cr	28%
Tungsten	W	3%
Niobium	Nb	4%
Silicon	Si	1%
Others	Fe / Mo	1%

ALLOY CHARACTERISTICS
Solidus-liquidus temperature: 1253 - 1304°C
Casting temperature: 1450 °C
Highest firing temperature: 950°C
Thermal expansion coefficient: (25 - 500°C) 14,1 x 10⁻⁶ K⁻¹/ (25 - 600°C) 14,5 x 10⁻⁶ K⁻¹
Density: 8,4 g/cm³
Vickers Hardness: 324 HV10
Percent elongation at fracture: 3,4%
Yield load strength: 475 MPa
Modulus of elasticity: 194 GPa
Colour: White
Recommended ceramics: GC Initial™ MC

INSTRUCTIONS FOR USE

I. WAX-UP

Wax walls should be 0,5 mm thick to ensure that after the manufacturing process the metal wall thickness is at least 0,3 mm. Avoid sharp edges and undercuts. Frameworks need to be anatomically reduced ensuring a uniform thickness of the ceramic layer. Connectors should be modelled as strong and as high as possible.

Single crowns: can be directly sprued. Suggested round casting sprues: Ø 3,0 mm to 4,0 mm. Bridges: crossbar spruing technique is advised. Connections to crossbar Ø 3 to 3,5mm. Distance from modulation to crossbar 2,5mm. Thickness of the crossbars Ø 5mm. Thickness of sprues to crossbar Ø 3 to 4mm.

The higher the volume of the casting the thicker the sprue should be. The casting sprue depends on the object dimensions.

II. INVESTING

Use only phosphate bonded investment material for Crown & Bridges. GC Fijvest Premium and GC Fijvest II are optimized for producing high precision fit castings in combination with smooth cast surfaces. Guidelines for optimal use of GC's phosphate bonded investments for Crown & Bridge can be found in the "GC Focus Edition – Phosphate Bonded Investments for Crowns and Bridges". A pre-heating and temperature of 850-950°C is advised. Follow the manufacturer's Instruction for Use related to the pre-heating times before casting.

III. CASTING

We recommend to only use fresh alloys in own crucibles. Sprues and cones can only be re-used once, given that new metal coming from the same lot is added in equal quantity.

Melt GC Initial™ CAST NP in pre-heated crucibles used for this alloy only.

Amount of casting alloy to use = wax framework weight with sprues / wax density x metal density

High-frequency melting equipment
When last ingots have melted in the pre-heated crucible and the casting immediately shadow has disappeared, immediately start the casting process. Do not use any flux.

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IV. PREPARATION FOR VENEERING

Use tungsten burs to grind and finish the prepared frameworks. Use always the same grinding instrument to avoid contamination. Do not use diamond burs. Grind the surface with slow speed and low pressure. Avoid the surface in one direction only to avoid overlaying layers (to prevent air bubbles after the first Opaque firing). Sandblast the prepared framework using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a pressure of 2,5-3 bar., then steam clean. The framework should have a homogeneous grey surface. Proceed with firing process using a compatible veneering ceramic following the manufacturer's instructions.

VI. VENEERING

a) Veneering with GC Initial™ MC: please refer to GC Initial™ MC Technical Manual
b) Veneering with other ceramics: please refer to the manufacturer's instruction for use
c) Veneering with acrylics: please refer to the manufacturer's instruction for use.

PACKAGING

GC Initial™ CAST NP, 250g
GC Initial™ CAST NP, 500g
GC Initial™ CAST NP, 1000g

CAUTION

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Before prosthesis application verify if other metal implants are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- Metal dusts and smoke are dangerous for health. Use exhaust fans while casting and sandblasting.
- We recommend to keep batch records and link it to patient's file to allowing its complete traceability.
- This product does not need any special preservation precautions. Nonetheless it has to be employed within the expiry date on the packaging.
- Dispose of all waste according to local regulations.

GUARANTEE

All recommendations are based on GC firings and GC metal tests and are therefore only to be used as guidelines. Any claim for compensation GC receives, in any event, shall not exceed the value of the individual order of GC Initial™ CAST NP.

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