Prior to use, carefully read the instructions **EN**

GC Initial™ INmetalbond

Buffer between GC Initial metal ceramics and dental alloys.

GC Initial INmetalbond INmetalbond is used as a thin layer between the alloy and the first Opaque layer. The INmetalbond blocks the escaping metal oxides and neutralizes differences in the expansion coefficient. The bonder does not increase the bonding strength as such but allows a wider span of CTE compatibility.

Recommended indications Can be used on all PFM alloys, precious as well as non-precious.

Instruction For Use

- Metal framework needs to be prepared according to manufacturer’s instructions.
- Oxidation firing again according to manufacturer’s instructions.
- Stir paste before using.
- Apply the Bonder in a thin layer but mask the metal completely.

IMPORTANT:

- To change the consistency of the INmetalbond paste, use the “GC Initial Paste Opaque Thinner”.
- Only use it in very small quantities.
- To avoid drying out, close the cap after using.
- The fired Bonder should have a yellowish, slightly shiny surface.

NOTE:

- The colour of the fired bonder can vary depending on the composition of the alloys.
- When using NPA (non-precious alloys) with the INmetalbond, the final temperature of the Opaque Washbake does not need to be increased by 20°C anymore. e.g.: 1st Opaque Firing on NPA: GC Initial MC = 940°C GC Initial LF = 830°C

Firing instructions

Preheating Temp.	Drying Time	T* increase	Vacuum	Final Temp.	Holding Time
550°C	6min	80°C/min	Yes	980°C	1 min no vacuum

Storage
Recommended for optimal performance, store in a cool, dry place.

Shelf life
GC Initial INmetalbond: 5 years from date of manufacture

Packages
GC Initial INmetalbond, 2x4gr.

CAUTION
Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: http://www.gceurope.com or for The Americas http://www.gcamerica.com They can also be obtained from your supplier.

Last revised: 11/2017

CE 0086

'GC.'

MANUFACTURED by

GC EUROPE N.V.
 Research-Park Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00

DISTRIBUTED by
GC CORPORATION
 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC EUROPE N.V.
 Research-Park Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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Prior to use, carefully read the instructions

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GC Initial INmetalbond

INmetalbond is used as a thin layer between the alloy and the first Opaque layer. The INmetalbond blocks the escaping metal oxides and neutralizes differences in the expansion coefficient. The bonder does not increase the bonding strength as such but allows a wider span of CTE compatibility.

Recommended indications

Can be used on all PFM alloys, precious as well as non-precious.

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Firing instructions

Preheating Temp.	Drying Time	T° increase	Vacuum	Final Temp.	Holding Time
550°C	6min	80°C/min	Yes	980°C	1 min no vacuum

Storage

Recommended for optimal performance, store in a cool, dry place.

Shelf life

GC Initial INmetalbond: 5 years from date of manufacture

Packages

GC Initial INmetalbond, 2x4gr.

CAUTION

Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

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