Hera

Heraenium® Pw

DE	Gebrauchsanweisung	3
GB	Instructions for use	11
FR	Mode d'emploi.....	19
ES	Instrucciones de uso	28
IT	Istruzioni per l'uso.....	37
PT	Instruções de uso	45
NL	Gebruiksaanwijzing.....	54
SE	Bruksanvisning	62
DK	Brugervejledning.....	70
NO	Bruksanvisning	78
FI	Käyttöohjeet.....	85
GR	Οδηγίες Χρήσης	93
CZ	Návod k použití.....	102
HU	Használati utasítás.....	110
LV	Lietošanas instrukcija	118
LT	Naudojimo instrukcija	126
PL	Instrukcja obsługi.....	134
TR	Kullanma talimatı	143
RU	Инструкция по применению.....	151
BG	Инструкции за употреба.....	160
UA	Інструкція по застосуванню	168

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Instructions for use Heraenium® Pw

According to EN ISO 9693 and EN ISO 22674

Heraenium PW is a cobalt chrome alloy for fabricating crowns and bridges with ceramic or acrylic facings. Heraenium PW is distinctive for bonding well to ceramic as well as for its optimal physical composition and a wide range of applications. The alloy is easy and straightforward to process given its low hardness.

1. Waxing up: Coping patterns must be at least 0.4 mm thick – this ensures that the coping remains at least 0.3 mm after trimming. To ensure that the ceramic forms a uniformly thick layer, the copings and pontics must be waxed up to approximate the reduced dental anatomy. Avoid sharp line angles and undercuts. For acrylic facings copings must also be 0.3 mm thick and retention beads (0.4 mm dia.) must be used.

2. Attaching sprues and investing: We recommend using a reservoir bar sprue system for bridgework. Approx. 5 mm long and 3.5 mm diameter feeder sprues are attached to the palatal or lingual surfaces of the wax patterns at an angle of 45°. Each

bridge unit requires its own sprue. Large molar copings or solid pontics require two. The feeder sprues must always be attached to the thickest section of the pattern. They are connected by a 5 mm diameter horizontal bar sprue. 4 mm diameter sprues are used for single units. The higher the volume of the casting the thicker the sprues should be. Fine-grained, phosphate-bonded, graphite-free and plaster-free investment materials for crowns and bridgeworks such as Heraest Onyx and Moldavest exact are optimised for producing precision fit castings and perfect cast surfaces.

3. Casting: The correct amount of Heraenium PW required for casting is calculated by multiplying the weight of the wax pattern by the density of the alloy divided by the density of the wax (approximately 0.93 g/cm³).

Density of Heraenium PW: 8.9 g/cm³

Formula: Weight of pattern x 8.9 g/cm³ ÷ 0.93 g/cm³.

Casting with induction-heated vacuum/pressure casting machines: Induction-heated vacuum/pressure casting machines by Kulzer are the most suitable for melting and casting. During pre-melting, the ingots are melted until the last ingot has sunk into the melt, and the edges or dull film of the ingot are no longer visible. Melting is then interrupted and the casting chamber is allowed to resume normal pressure and

opened. The casting ring is then inserted and the main melting process commenced. The casting process is started 6 seconds after the dull film from melting has disappeared.

Casting with a motor-driven centrifugal casting machine and melting with a casting torch: Ignite the oxyacetylene gas. Open both valves fully and adjust the flame by reducing the acetylene flow until approximately 3 mm long blue cones appear on the orifices of the casting torch.

Set a neutral flame. Place the Heraenium P ingots in the pre-heated crucible. Hold the torch approximately 40 mm above the top edge of the crucible. Move the torch around the ingots in circles to heat them uniformly until they melt and fuse. Once the ingots have fused, start the centrifugal casting machine immediately.

After casting, bench cool the casting ring. Devest the castings very carefully. Remove the investment material with care. Sand-blast with aluminium oxide abrasive (110–125 μm).

4. Trimming and preparing for ceramic facings: After casting, the copings are trimmed with sharp tungsten carbide cutters in one direction only and without exerting high pressure. Do not use diamond burs. Do not use rubber polishers on the surface to be veneered. Avoid sharp line angles and undercuts.

Sandblast the surfaces with aluminium oxide (100–125 µm, 4 bar) followed by steam cleaning.

 **Metal dust is hazardous to health.**

Attention! When trimming and sandblasting, use a suitable extraction system and facemask (type FFP3-EN 149-2001).

5. Oxide firing: 950°C (1742°F) for 10 minutes in a vacuum. The oxide layer must be stained uniformly. Carefully remove the oxide layer on the facing surfaces with aluminium oxide (110–125 µm, 4 bar) using a non-recycling sandblaster, and then clean (e.g. steam clean).

Note: An oxide firing is not necessary when facing with HeraCeram.

6. Veneering with ceramics:

a) Veneering with HeraCeram: The HeraCeram processing instructions from Kulzer apply.

b) Veneering with other ceramics: When veneering with other ceramics, the other ceramic manufacturers' instruction for use must be observed.

c) Stress-relief cooling: Stress-relief cooling is recommended for large-span and more massive frameworks: setting a cooling period with opening of the ceramic furnace (3–5 minutes)

or leaving the object on the firing tray in the radiated heat from the open firing chamber.

7. Acrylic facings: Process the facing material according to the manufacturer's instructions.

8. Soldering before firing: Special Heraenium P solder 1 is recommended for soldering Heraenium Pw frameworks.

Soldering after firing (furnace soldering recommended):

This soldering can be carried out with cobalt chrome/gold solder 750. Apply a thin coat of flux Hera SLP 99 (Please pay attention to the information on the equivalent Safety Data Sheet!) to the soldering site, and make sure that the flux does not contact the ceramic.

9. Laser welding: The special Kulzer CoCR laser welding wire, dia 0,5 mm is ideal for welding Heraenium Pw frameworks.

10. Recasting: The particular physical and chemical qualities of the alloy can be guaranteed only if new material is used. We recommend using only new material for ceramic facings.

11. Applying Blendgold Neu: Blendgold Neu enhances the shade reproduction qualities and aesthetics of porcelain. Condition the surfaces as described in section 4 before applying and firing a very thin coat of opaquer. Extrude Blendgold Neu

out of the syringe directly onto a clean brush, and apply a thin coat of Blendgold Neu onto the opaquer. Heat the ceramic furnace to 400°C (752°F), and dry the Blendgold Neu in the heat radiating from it. The temperature of the object should reach approximately 180°C (356°F). Raise the temperature to 820°C (1508°F) at a rate of 55°C (131°F)/min., and fire without a vacuum or holding time. Remove the object on the firing tray from the furnace, and allow it to bench cool. Continue processing the ceramic as usual. Simply applying a coat of Blendgold in the cervical region alone may be adequate. Blendgold Neu should not be directly applied on Heraenium Pw. If a warm, yellowish gold colour is needed, we recommend HeraCeram Intensive Opaque gold.

12. Chemical composition in mass percent

Co: 55.2; **Cr:** 24.0; **W:** 15.0; **Fe:** 4.0; **Si:** 1.0; **Mn:** 0.8; **N:** <0.1
(beryllium-free, cadmium-free)

13. Technical data

Type 4 alloy after casting and type 5 after ceramic firing according to EN ISO 9693 and EN ISO 22674

Density: 8.9 g/cm³

Investment material: Phosphate-bonded, e.g. Heravest Onyx, Moldavest exact

Preheating temperature: 950°C (1742°F)

Melting crucible: Ceramic non-precious metal

Melting range: 1320°C– 1400°C (2408°F–2552°F)

Casting temperature: 1550°C (2822°F), 6 seconds after dull film disappears

Oxide firing: 950°C (1742°F)/10 min. in a vacuum (oxide firing is not required for facings with HeraCeram)

CTE 25–500°C (77°F–932°F): 14.3 $\mu\text{m}/\text{m}/\text{K}$

Solder: Heraenium P solder 1,
cobalt chrome gold solder 750

Hardness after casting: 275 HV 10

Hardness after ceramic firing: 290 HV 10

Yield point 0.2%: 530 MPa after ceramic firing
490 MPa after casting

Tensile strength after firing: 770 MPa

Elongation at rupture after firing: 8 %

Elasticity modulus: 208 GPa

Contraindications: If components of this alloys trigger hypersensitivity (allergies) they should not be used.

Side effects: Hypersensitivity reactions (allergies) as well as electrochemically based local sensory disturbances, taste disturbances and irritation of the oral mucosa have been described in isolated cases.

Interactions with other dental alloys: With approximal or antagonist contact with dentures of different alloys galvanic effects may occur. Other materials must be used instead if contact with other alloys causes persistent local paraesthesia on an electrochemical basis.

14. Instructions

Consult the material safety data sheets or national regulations for disposal.

Heraenium® = registered trademark of Kulzer GmbH

We retain the right to make technical changes.

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