en	Instructions for use/Technical description Instruments without tissue contact
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung Instrumente ohne Gewebekontakt
fr	Mode d'emploi/Description technique Instruments sans contact avec les tissus
es	Instrucciones de manejo/Descripción técnica Instrumental sin contacto con tejidos
it	Istruzioni per l'uso/Descrizione tecnica Strumenti senza contatto con i tessuti
pt	Instruções de utilização/Descrição técnica Instrumentos sem contacto com os tecidos
nl	Gebruiksaanwijzing/Technische beschrijving Instrumenten zonder weefselcontact
da	Brugsanvisning/Teknisk beskrivelse Instrumenter uden vævskontakt
sv	Bruksanvisning/Teknisk beskrivning Instrument utan vävnadskontakt
fi	Käyttöohje/Tekninen kuvaus Instrumentit, jotka eivät ole kosketuksissa kudokseen
lv	Lietošanas instrukcijas/tehniskais apraksts Instrumenti bez kontakta ar audiem
lt	Naudojimo instrukcija/techninis aprašas Su audiniais nesiliečiantys instrumentai
ru	Инструкция по применению/Техническое описание Инструменты, не контактирующие с тканями
cs	Návod k použití/Technický popis Nástroje bez kontaktu s tkání
pl	Instrukcja użytkowania/Opis techniczny Instrumenty niestykające się z tkanką
sk	Návod na použitie/Technický opis Nástroje bez dotyku s tkanivom
hu	Használati útmutató/Műszaki leírás Szövettel nem érintkező műszerek
sl	Navodila za uporabo/Tehnični opis Instrumenti brez stika s tkivi
hr	Upute za uporabu/Tehnički opis Instrumenti bez kontakta s tkivom
ro	Manual de utilizare/Descriere tehnică Instrumente fără contact cu țesutul
bg	Упътване за употреба/Техническо описание Инструменти без контакт с тъканите
tr	Kullanım Kılavuzu/Teknik açıklama Doku teması olmayan aletler
el	Οδηγίες χρήσης/Τεχνική περιγραφή Όργανα που δεν έρχονται σε επαφή με ιστό

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply for Instruments without tissue contact.

- ▶ For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eIFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

Instruments without tissue contact are used in universal surgical application in all fields.

2.1.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.1.3 Contraindications

No known contraindications.

2.2 Safety information

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Sterility

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging.

2.3 Application

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the product.

3. Reprocessing procedure

3.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

3.2 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

3.3 Preparation before cleaning

- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Carry out non-fixating/NaCl-free pre-cleaning immediately after use.

3.4 Cleaning/Disinfection

3.4.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
 - that are approved for aluminum, plastics and stainless steel,
 - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water

3.5 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

3.5.1 Manual cleaning and wipe disinfecting

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT	1	-	-	17 % Propan-1-ol, 0.23 % Didecyldimethylammonium chloride*
II	Wipe disinfection	RT	≥1	-	-	17 % Propan-1-ol, 0.23 % Didecyldimethylammonium chloride*

RT: Room temperature

*Recommended: Meliseptol® Wipes sensitive (B. Braun)

Phase I

- ▶ Remove any visible residues with a disposable disinfectant wipe.

Phase II

- ▶ Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- ▶ Observe contact time requirements (at least 1 min).

3.6 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.6.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the product for damage to the spiral element.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.6.2 Functional test

- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.7 Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

4. Technical service

⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

5. Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

⚠ WARNING

Risk of injury due to sharp-edged and/or pointed products!

- ▶ When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

Note

The user institution is obliged to reprocess the product before its disposal, see Reprocessing procedure.

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