
Introcan[®], Introcan[®]-W, Introcan[®] Certo, Introcan[®]-W Certo


B | BRAUN

LLDorder 5362 - grunewald 210261

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PUR, PP, ABS, chrome-nickel steel

Indications

Creation of a secure peripheral venous, central venous or arterial (with FEP only) access. Blood transfusions or infusion of I. V. solutions suitable for administration via peripheral veins. Intermittent intravenous drug administration. Prophylactic creation of a secure venous access in patients who may require urgent intravenous drug administration, in particular prior to diagnostic or therapeutic procedures. Intermittent blood sampling, especially from arterial (only for FEP) sites.

Invasive (intra-arterial) measurement of arterial blood pressure (with FEP only).

14-22 gauge catheters may be used with power injectors for which the maximum pressure setting is 300 psi/21 bar. Facilitate the placement of Vascular Access Devices such as guidewires, indwelling central venous catheters, peripherally inserted central catheters, and midline catheters into the vascular system.

Contraindications

Introcan® I. V. Catheter should not be used in patients with known hypersensitivity to any of the materials employed.

Risks

Depending on how long the cannula is left in situ, on the type and amount of infusions or injections administered, and on individual predisposition, thrombophlebitis may occur in the accessed vein.

For FEP only: In arterial puncture, the artery may occlude in rare cases due to thrombotic or embolic complications, resulting in ischaemia.

Duration of use

Change according to CDC Guidelines and / or Hospital or Institutional protocols. The puncture site should be checked at regular intervals. Introcan® I. V. Catheter should be removed in the event of local or systemic signs of infection.

Warnings

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient. After withdrawal, do not reintroduce the steel needle into the catheter, as the latter may be cut off, leading to catheter embolism.

For FEP only: Always clearly mark arterial lines to avoid inadvertent injection. Verify adequate collateral circulation prior to arterial puncture.

Prior to use with power injectors, ensure that secure connection exists between the catheter and power injector. Patency of catheter must be established prior to use with power injector.

Use only if packaging is intact. Do not re-sterilise.

Application

NOTE:

Do not rotate catheter prior to insertion.

1. After disinfection of the puncture site and removal of the protective cap, puncture a suitable vein. If venipuncture was successful, blood will immediately be visible inside the transparent grip part.
2. Advance the catheter further into the vein, while slightly withdrawing the steel needle.
3. For peripheral venous access:
Use adhesive tape, secure catheter to the skin. The steel needle still in situ minimizes spillage of blood.
4. Before removing the steel needle, compress the vein at the tip of catheter to prevent spillage of blood. Remove needle by pulling needle straight back with a swift, continuous motion parallel to the skin (minimize rotation of needle). Dis-

pose of needle immediately into sharps container.

5. When used for peripheral access: Connect to infusion line and cover puncture site with a sterile dressing.

6. If used to facilitate placement of a Vascular Access:

Device follow manufacturer's instructions for use of the specific device being placed.

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B. Braun Melsungen AG