HOFFMANN'S KUPFER CEMENT

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Instruction for the dentist

Composition:

powder contains: zinc oxide, magnesium oxide, copper rhodanide;

liquid contains: o-phosphoric acid

Areas of Application:

Fixation of crowns, bridges and inlays, orthodontic work, underfilling cement with bacteriostatic effect.

Contra-indications:

Direct application to dentin near the pulp or to exposed pulp.

Side effects:

Irritations of the pulp may occur after cementing.

Interactions with other substances:

No known interactions

Dispensing Instructions:

Underfilling consistency: 2.1 g powder : 1.0 g liquid

Fixation consistency: 1.5 g powder : 1.0 g liquid

Processing Recommendations:

A clean, dry glass plate at 18 - 20°C and a non-discolouring metal spatula are recommended.

Mixing:

Divide the measured quantity of powder into individual portions in order to achieve a homogenous consistency. Begin with the smallest quantity, quickly mix the powder into the measured liquid, a portion at a time. Discard remaining powder or impure liquids. Dentin areas near the pulp must be protected.

Mixing Time:

1 minute 30 seconds

Working Time:

3 minutes (at 23°C)

Net Setting Time:

4.0 - 7.5 minutes (at 37°C). The times indicated apply from the end of mixing. The cement must be applied to the dry tooth surface and must set without contact to moisture.

Notes:

- Close containers immediately after removing material.
- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.

Colour: 3 (yellowish white)

Package Sizes: Powder: 100 g, liquid: 40 ml

Hoffmann's Copper Cement

- Do not use after the expiration date
- Compliant with ISO 9917

HOFFMANN DENTAL MANUFAKTUR GmbH, Wangenheimstraße 37/39, 14193 Berlin, Germany, www.hoffmann-dental.com

04/2007

HOFFMANN'S COPPER CEMENT

Instruction for the dentist

Composition:

powder contains: zinc oxide, magnesium oxide, copper rhodanide;

liquid contains: o-phosphoric acid

Areas of Application:

Fixation of crowns, bridges and inlays, orthodontic work, underfilling cement with bacteriostatic effect.

Contra-indications:

Direct application to dentin near the pulp or to exposed pulp.

Side effects:

Irritations of the pulp may occur after cementing.

Interactions with other substances:

No known interactions

Dispensing Instructions:

Underfilling consistency: 2.1 g powder : 1.0 g liquid

Fixation consistency: 1.5 g powder : 1.0 g liquid

Processing Recommendations:

A clean, dry glass plate at 18 - 20°C and a non-discolouring metal spatula are recommended.

Mixing:

Divide the measured quantity of powder into individual portions in order to achieve a homogenous consistency. Begin with the smallest quantity, quickly mix the powder into the measured liquid, a portion at a time. Discard remaining powder or impure liquids. Dentin areas near the pulp must be protected.

Mixing Time:

1 minute 30 seconds

Working Time:

3 minutes (at 23°C)

Net Setting Time:

4.0 - 7.5 minutes (at 37°C). The times indicated apply from the end of mixing. The cement must be applied to the dry tooth surface and must set without contact to moisture.

Notes:

- Close containers immediately after removing material.
- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.

Colour: 3 (yellowish white)

Package Sizes: Powder: 100 g, liquid: 40 ml

Hoffmann's Copper Cement

- Do not use after the expiration date
- Compliant with ISO 9917

HOFFMANN DENTAL MANUFAKTUR GmbH, Wangenheimstraße 37/39, 14193 Berlin, Germany, www.hoffmann-dental.com

04/2007