Lucitone 199[®]

Denture Base Resin

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Mode d'emploi

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Indications for Use

Lucitone 199[®] Resin is indicated for the fabrication of prosthodontic appliances.

Contraindications

Lucitone 199 Resin is contraindicated for patients and users with a history of allergic reaction to methyl methacrylate monomer.

Warnings

- Material contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wear protective gloves and protective clothing. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice.
- Avoid inhalation or ingestion. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Work in a well ventilated area. Liquid contact with eyes may cause possible corneal damage. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects.
Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.
Ingestion: Contact your regional Poison Control Center immediately.
Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician. Wash skin with soap and water.

Precautions

- When grinding prosthodontic resins, proper ventilation, masks, and vacuum systems should be used.

- Store at 60 °F - 80 °F (16 °C - 27 °C) away from moisture and direct sunlight. The liquid contains unsaturated monomers which may polymerize prematurely if stored at excessively high temperature and/or in sunlight. The powder contains an organic polymerization initiator which may degrade if stored at excessively high temperature.
- Lucitone Liquid and Lucitone 199 Resin uncured resin are hazardous materials. Dispose in accordance with Federal, State, and local regulations.
- Lucitone Liquid contains methyl methacrylate monomer, a flammable liquid with a flash point of 50 °F (10 °C). Keep away from heat, sparks, and open flame.
- Use liquid in a well ventilated area. Replace cap when not in use.
- Do not pack denture resin after work time has elapsed or material has become rubbery.
- Use of wax solvents is not recommended. Residual solvent may cause poor bonding of teeth to acrylic base.
- Corneal damage, headache, nausea, vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)
- Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthesis in warm water for several days.
- Particulates will be generated when grinding acrylic resins. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.
- Dispose of contents/container in accordance with local and national regulations.

Adverse Reactions

- Corneal damage, headache, nausea, vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)
- Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthesis in warm water for several days.
- Particulates will be generated when grinding acrylic resins. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

Symbols on Product Labels

-  Manufacturer
-  Batch Code
-  Reorder Number
-  Manufacture Date
-  Expiration Date
-  Consult Instructions for Use
-  Caution
-  Keep Away from Sunlight
-  CE European Conformity Marking
-  Medical Device
-  U.S. Prescription Only
-  Store Between

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

| | | |
|--|--|---------|
| Residual MMA Content | 2.2% (% mass fraction) | Maximum |
| Type and Class (according to ISO 20795-1:2013)..... | Type 1, Class 1 | |
| Storage temperature for powder and liquid | 60-80 °F (16-27 °C) | |
| Powder/Liquid ratio | 21g (32 cc)/10 ml | |
| Mixing time (time required to wet all particles) | 15-30 seconds | |
| Time to reach packing plasticity @ 73°F ± 2°F(23 ± 1°C)..... | 9 ± 2 minutes | |
| Working time | 10 ± 4 minutes | |
| Material used to prepare mold | gypsum | |
| Temperature of mold when packing | room temperature to 110 °F (43 °C) | |
| Recommended cure time and temperature | | |
| 1st stage | 1 ½ hours @ 163 °F (73 °C) | |
| 2nd stage..... | ½ hour @ 212 °F (100 °C) | |
| Alternate cure time and temperature | | |
| 1st stage | 9 hours @ 163 °F (73 °C) | |
| 2nd stage..... | ½ hour @ 212 °F (100 °C) | |
| Method of cooling flask, time and temperature | | |
| 1st stage | ½ hour in air at 60-80 °F (16-27 °C) | |
| 2nd stage..... | ¼ hour in water at 60-80 °F (16-27 °C) | |

STEP-BY-STEP INSTRUCTIONS

PROPERTIES AND RECOMMENDED PROCEDURES

The preceding table describes typical properties of Lucitone 199 Resin, plus recommended temperatures and other important processing information. Further amplification of this condensed information is found later in the text.

FLASKING

Use conventional detail compression molding methods and gypsum materials for flasking.


WAX ELIMINATION

Soften wax in boiling water for approximately 6 minutes. Separate flask and remove wax by flushing with a solution of boiling water to which a detergent has been added. Use clean boiling water for final flush. **DO NOT USE WAX SOLVENTS.** Be sure case is free of wax and grease.

TINFOIL SUBSTITUTE

Apply separator liberally to areas of the warm mold (approx. 120 °F (49 °C) that will contact the resin. Remove excess film with a dry brush. The film of the separator should be thoroughly dry (approx. 3 minutes following application) before the resin is packed.

IMPORTANT - Measure powder and liquid accurately. Do not whip air into mix by excessive spatulation.


MIXING

To assure color uniformity, shake powder jar (with a rocking motion) to prevent concentration of small particles at the bottom of jar. Add 1 unit powder 21 g (32cc) to 10 ml liquid.

Stir sufficiently (15 seconds) to assure wetting of all powder particles. Do not over-mix. Cover mixing jar and allow material to reach packing consistency (approx. 9 minutes at room temperature of 73 °F ± 2 °F (23 °C ± 1 °C).


Powder measure and liquid graduate are supplied in packages of Lucitone 199 Resin. The following table provides the correct ratios for various cases.

| CASE | POWDER | LIQUID |
|-------------|--------|--------|
| Extra Large | 42 cc | 13 ml |
| Large | 32 cc | 10 ml |
| Average | 26 cc | 8 ml |
| Partial | 19 cc | 6 ml |


PACKING

Lucitone 199 Resin should be packed in warm flask (room temperature to 110 °F (43 °C)). The correct packing consistency is identified when the resin dough is free from stickiness but not rubbery. Remove resin dough from jar and condense with finger pressure into mold.

TRIAL PACKING

Lucitone 199 Resin does not require trial packing except in extra-thick cases. However, trial packing may be employed if the operator desires. The flask should be closed **slowly** in all cases.

Close flask slowly

If the no trial packing technique is employed, it is essential that a sufficient amount of material be used to assure a full mould. Excess resin squeezing out of the complete periphery of the two halves of the flask during closure is a good indication of sufficient volume of material.

WORK TIME

The time of Lucitone 199 Resin is approximately 10 minutes. Never pack after work time has elapsed or material has become rubbery.

CURING

Submerge closed flask (locked by compress or spring clamp) in water at 163 °F \pm 2 °F (73 °C \pm 1 °C) for 1 ½ hours. Follow by ½ hour in boiling water. A periodic check of water bath temperature with an accurate thermometer is recommended. Alternate cure: 9 hours in water bath at 163 °F \pm 2 °F (73 °C \pm 1 °C). Follow by ½ hour in boiling water.

CURING OF EXTREMELY THICK CASES

Extremely thick cases require special handling to prevent porosity and excessive surface shrinkage. The following procedure should be employed:

1. Pack material when it reaches a heavy, putty-like consistency (approx. 12 minutes at $73\text{ }^{\circ}\text{F} \pm 2\text{ }^{\circ}\text{F}$ ($23\text{ }^{\circ}\text{C} \pm 1\text{ }^{\circ}\text{C}$)).
2. Condense well during packing by closing flask slowly and trial packing at least 3 times.
3. Closed flask should remain under pressure for 30 minutes prior to curing.
4. Use alternate cure method as specified.

COOLING

Bench cool the clamped flask at room temperature for 30 minutes. Then immerse in cool water $60\text{-}80\text{ }^{\circ}\text{F}$ ($16\text{-}27\text{ }^{\circ}\text{C}$) for 15 minutes before deflasking.

FINISHING AND POLISHING

The denture should be finished and polished in the usual manner.

REPAIRING AND RELINING

Relines or repairs to the completed denture may be made with LUCITONE HIPA (or other repair materials) in accordance with instructions provided. If desired, repairing or relining can also be done by the usual technique using Lucitone 199 Resin itself. Cure in water at $163\text{ }^{\circ}\text{F}$ ($73\text{ }^{\circ}\text{C}$) for 9 hours. Follow by 1/2 hour in boiling water.

INFECTION CONTROL PROCEDURE

INSTRUCTIONS FOR DENTURE DISINFECTION

Note: *Instructions are not intended for routine denture care.*

| | Method | Warning |
|--|---|---|
| Cleaning & Disinfection: Manual | <ul style="list-style-type: none">• Saturate a lint-free wipe with an alcohol quaternary-based solution.• Wipe the entire appliance ensuring to wipe all grooves and edges well.• Spray the entire surface with the alcohol quaternary-based solution to thoroughly wet the appliance - <i>maintain wet for the contact time as recommended per the manufacturer's instructions.</i>• After the manufacturer recommended contact time has elapsed, rinse the appliance under potable water for a minimum of 30 seconds.• Allow the appliance to air dry completely. | *Follow the manufacturer's instructions for solution contact time. |

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Registered Trademark - Chile, Guatemala, Mexico, The Trademark under which this product is sold is registered in various countries.
York, PA 17404 U.S.A.
**DeguDent GmbH**