Deutsch (de)  
English (en)  
Español (es)  
Français (fr)  
Italiano (it)

Stirnleuchte MedLED® Chrome

Headlight MedLED® Chrome

Lámpara frontal MedLED® Chrome

Lampe frontale MedLED® Chrome

Lampada frontale MedLED® Chrome

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Instructions for Use

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Mode d'emploi

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**REF** 90-085-58-10

Revision 01

Date of Release: 2021-04

Deutsch .....	9
English .....	32
Español .....	54
Français .....	77
Italiano .....	100

**Symbolerklärung**

**Symbol Explanation**

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**Explications des symboles**

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## Table of Contents

1	General Information .....	32
1.1	Manufacturer .....	32
1.2	Hotline.....	32
1.3	Adverse Incident Reporting Requirement .....	33
1.4	Notes on this Document .....	33
1.5	Terms & Acronyms.....	33
1.6	Validity of this Document.....	33
1.7	Accompanying Documents .....	33
2	Scope of Delivery.....	34
2.1	Inspecting the Delivery for Completeness and Correctness .....	34
3	Intended Use .....	35
3.1	Intended Use .....	35
3.2	Indications.....	35
3.3	Contraindications .....	35
3.4	Possible Side Effects.....	35
3.5	Patient Target Groups .....	35
3.6	Users .....	35
3.7	Ambient Conditions During Application .....	35
3.8	Limitations of Use .....	36
3.9	Warnings .....	36
4	Application .....	37
4.1	Description of the Components .....	37
4.1.1	Structure, Operating Principle and Performance Specifications .....	37
4.1.2	Materials Used.....	37
4.1.3	Variants / Combinations .....	37
4.1.4	Combination Products and Accessories.....	38
4.2	Before Initial Use.....	38
4.3	Application .....	38
4.3.1	Headlight.....	38
4.3.2	Batteries.....	41
4.3.3	Charging Unit .....	42
4.3.4	Charging the Batteries on the Headlight .....	42
4.3.5	Charging the Batteries in the 4-Bay Desktop Charger .....	43
4.4	After Completion of the Application.....	43
5	Cleaning, Disinfection and Sterilization.....	44
5.1	Limitations and Restrictions to Cleaning, Disinfecting and Sterilizing .....	44
5.2	Pre-Treatment at the Application Location prior to Cleaning .....	44
5.3	Preparation before Cleaning.....	44

5.4	Cleaning and Disinfection .....	45
5.4.1	Manual Cleaning and Disinfection .....	45
5.4.2	Machine Cleaning and Disinfection .....	46
5.5	Inspection, Functional Check, Maintenance .....	46
5.5.1	Inspection and Functional Check .....	46
5.5.2	Maintenance .....	46
5.6	Packaging .....	46
5.7	Sterilization .....	46
5.8	Storage and Transport .....	46
6	Maintenance .....	47
6.1	General Notices .....	47
6.2	Regular Checks to Be Carried Out by the Operator .....	47
6.2.1	Safety Check (SC) .....	47
6.3	Maintenance by the Manufacturer / Periodic Maintenance .....	47
6.4	Format of the Serial Number .....	47
6.5	Headlights, Accessories & Spare Parts .....	48
7	Ecological Information / Disposal .....	50
7.1	Packaging .....	50
7.2	Batteries .....	50
7.3	Disposal of the Product .....	50
7.4	National Regulations .....	50
8	Technical Data .....	51

## 1 General Information

### 1.1 Manufacturer

We thank you for buying one of our products.

This product bears the CE mark. This means that it meets the basic requirements on the safety and the performance of medical devices set forth by the respective applicable European regulations.

We are the manufacturer of this product:

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info@klsmartin.com · www.klsmartin.com

### 1.2 Hotline

Should you have any questions on how to handle the device or product or use it for clinical applications, please contact the Product Management:

Phone: +49 7461 706-352

Fax: +49 7461 706-312

---

### **NOTICE**

Each packaging and most products are marked with a serial number or a production lot number (LOT) and a catalog number (REF). Please always provide serial number or LOT and REF in case of a complaint.

---

### 1.3 Adverse Incident Reporting Requirement

All serious adverse events occurring in connection with the product must be reported to Gebrüder Martin and the respective authorities without delay.

### 1.4 Notes on this Document

**Possible danger to the life of patients, users and others when these Instructions for Use are disregarded**

Read Instructions for Use completely and follow them carefully. In particular, be sure to note all cautionary and warning notices.

This document refers to persons of all genders alike. Reference to various genders is avoided purely for reasons of improved readability.

The electronic version of these Instructions for Use can be requested at [www.klsmartin.com](http://www.klsmartin.com).

### 1.5 Terms & Acronyms

Acronym	Description
dpt.	Diopter
LED	Light-emitting diode(s)
OR	Operating room

### 1.6 Validity of this Document

This document applies to MedLED® Chrome Headlights.

### 1.7 Accompanying Documents

n/a

## 2 Scope of Delivery

Designation	REF
<b>MedLED® Chrome Set</b> consisting of:	See brochure
Headband: MedLED Chrome hardtop or softtop	
Bezel: MC3 or MC6 or MC7 or all three models together	
2 rechargeable lithium-polymer batteries	
USB charging cable	
power supply unit incl. adapter for EU, US, UK, AU	
carrying case	
these Instructions for Use	

### 2.1 Inspecting the Delivery for Completeness and Correctness

Immediately upon receipt, the goods must be checked for completeness and damages.

Notice of any damage must be given immediately.

After delivery of the product, check the original packaging and the packaging seal for damages.

### **3 Intended Use**

#### **3.1 Intended Use**

The headlight is designed to illuminate a surgical area or examination field.

It must not come into contact with the patient. Any other usage counts as misuse.

#### **3.2 Indications**

The MedLED® headlight is indicated to increase visibility only and is not intended to come into contact with the patient.

#### **3.3 Contraindications**

The MedLED® headlight system is contraindicated for use in neonatal transillumination, ophthalmic procedures, photosensitive patients or patients who have received photosensitizing agents.

#### **3.4 Possible Side Effects**

n/a

#### **3.5 Patient Target Groups**

No direct patient contact.

#### **3.6 Users**

Professional users:

The products are to be used exclusively by medical specialists.

Cleaning and disinfection must be performed by trained specialists.

#### **3.7 Ambient Conditions During Application**

The application is carried out in the intended areas.

Where applicable, the specific parameters for ambient conditions during the application (temperature range and air humidity) as stated on the product label must be observed.

For detailed information, see chapter 4.3 "Application", page 38.

---

### 3.8 Limitations of Use

---

**⚠ CAUTION****Possible risk of injury due to glare!**

Do not shine the headlight directly into anybody's eyes.

Never hold microscopes or other optical lenses in front of the headlight. The LED light source is intense and has a high secondary exposure potential due to its blinding effect. Temporary reduction in visual acuity and afterimages can occur, leading to irritation, annoyance, visual impairment, and even accidents, depending on the situation.

Use the rheostat for dimming to the appropriate setting.

---

### 3.9 Warnings

---

**⚠ WARNING****Possible danger to the life of others by shipping contaminated products!**

In case of return shipments, send only clean and disinfected products in sterile packaging.

---

**⚠ CAUTION****Possible danger of injury to the patient or user due to unsafe use!**

Switch off the device using the rheostat dimmer when not in use.

When in use, the aperture of the headlight must not be covered with a towel or any similar material, as it could otherwise become very warm. The bezel should not come into contact with any skin during use.

---


## 4 Application

### 4.1 Description of the Components

#### 4.1.1 Structure, Operating Principle and Performance Specifications

Detailed information can be found in the special brochure by the manufacturer.

#### MedLED® Chrome

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>1. Headstrap</li> <li>2. Headband</li> <li>3. Battery compartment</li> <li>4. Rheostat: On/off switch and brightness controller</li> </ul> | <ul style="list-style-type: none"> <li>5. Micro-B USB port for plugging in the bezel (🔌)</li> <li>6. Micro-B USB port for charging the batteries (🔋)</li> <li>7. Pivot joints</li> <li>8. Bezel with aluminum ring for aperture adjustment (beam and spot size)</li> </ul> |
|---|--|

#### 4.1.2 Materials Used

Designation	Material(s)
Plastics	Polypropylen, TPE
Metal	Aluminium
Others	Nylon bristles, Neoprene, Velcro fastener, LED

#### 4.1.3 Variants / Combinations

The products are available in various variants. Please use the special KLS Martin brochure for this purpose.

#### 4.1.4 Combination Products and Accessories

Combination products and/or accessories can be found in the special KLS Martin brochure or in chapter 6.5 "Headlights, Accessories & Spare Parts", page 48.

## 4.2 Before Initial Use

Remove protective caps or protective covers if present.

The products can be incorporated in the usual processing routine. In order to do this, follow the instructions from chapter 5 "Cleaning, Disinfection and Sterilization", page 44.

Before first use, switch off the headlight -Lever down rheostat- (④), insert rechargeable batteries (③) and charge overnight (⑥).

## 4.3 Application

### 4.3.1 Headlight

---

**⚠ CAUTION****Possible risk of injury due to glare!**

The LED light source is extremely bright. Its glare effect can lead to temporary loss of visual acuity and afterimages. Depending on the situation, this can be irritating or annoying, hamper the view or even lead to accidents.


- Do not point the headlight directly into the eyes of others.
  - Under no circumstances hold microscopes or other optical lenses in front of the headlight.
- 



**NOTICE**

The aluminum bezel / light engine may get warm during use.

This is normal and not a malfunction. Avoid skin contact to these areas during these times.

---

Designation	Picture
<p>Switch off the headlight (④).            Lever down (Off).</p>	
<p>Insert batteries (③) and charge overnight (⑥), see section 4.3.4 “Charging the Batteries on the Headlight”, page 42 or section 4.3.5 “Charging the Batteries in the 4-Bay Desktop Charger”, page 43.</p>	
<p>Remove Micro-B USB plug again (⑥).</p>	
<p>Turn on the headlight (④).            Lever up (On).</p>	
<p>Put on headlight and adjust to head size (①).</p>	
<p>Adjust brightness with the rheostat (④) as required. Stepless adjustment from “Off” (lever down) to “maximum brightness” (lever up).</p>	

Designation	Picture
<p>Adjust the light beam using the pivot joints (7).</p>	
<p>Adjust the beam and spot size on the lamp body (8). Do not cover the lamp body (8) while in use (with a cloth or other object), as it can get very hot.</p>	
<p>Check the LEDs of the battery status indicator, see section 4.3.2 "Batteries", page 41.</p>	

- We recommend connecting the headlight to the charging unit while it is not in use.
- To disconnect the mains cable, always pull on the plug and not the cable.
- If the device is to be stored for more than two weeks, remove the batteries and place them in the original packaging.
- When using the lamp in dental medicine: A yellow filter can be used to avoid premature polymerization of composite material.

### 4.3.2 Batteries

#### **Possible danger to life due to electric shock, fire, ignition, explosion, leaking or bursting of the batteries!**

To avoid injuries or property damage, e.g. due to electric shock, fire, ignition, explosion, leaking or bursting:

- Use only the batteries and charging units provided in the delivery for the headlight system purchased.
- Charge and use batteries only in ambient temperatures from 0°C to +30°C (+32°F to +86°F).
- Use the batteries only at a relative humidity of 45 - 85%.
- Store batteries only at ambient temperatures from 0°C to +25°C (+32°F to +77°F).
- Do not open the batteries, throw them into fire, insert them inverted, use them together with used batteries or batteries of another type, short-circuit them, disassemble them, heat them to over 60°C (140°F) or dispose of them into waste incineration plants.
- If the batteries swell, crack or get hot, do not charge or use them. They must be disposed of properly and be replaced by new ones. The battery used in this device may present a fire or chemical burn hazard if mistreated.
- Do not touch damaged/leaking batteries. Dispose of the batteries immediately and properly, observing the handling and safety measures for lithium-polymer batteries.

Pay attention to the 5 LED indicators on the front of the left battery compartment:

- Each of the 5 LEDs indicates a battery charging level of approx. 20%. Note that it can take a few seconds until the battery charging level is indicated.
- Replace the batteries as soon as only one LED is lit. As soon as all 5 LEDs are dark, only a short period remains until the batteries are fully depleted and the headlight is shut off.
- To ensure optimum performance of the headlight, the batteries should be replaced every 12 months.
- If a battery can be inserted or removed only with difficulty, it has expanded and must be put out of operation immediately.
- Do not charge the batteries in the operation field.

### 4.3.3 Charging Unit

- Use the charging unit only in dry environments in a closed area. Protect from rain, snow, inflammable materials, excessive heat and/or moisture.
- Do not connect the charging unit to an extension cable.
- Do not excessively bend the mains cable or place any heavy objects on it; do not pull, drag or strain the cable.
- Do not use the charging unit if the cable has been severed or the batteries and/or the 4-chamber charging unit have been dropped or damaged.
- Do not use damaged charging units, do not disassemble the charging unit or modify any parts of the charging unit.
- Do not use unapproved power supply units, and do not charge damaged batteries.
- Do not bring metallic objects into contact with the metal parts of the DC connectors on the charging unit or the battery contacts.

### 4.3.4 Charging the Batteries on the Headlight


#### **NOTICE**

The headlight battery packs or 4-bay charger console may warm up during charge mode. This is normal and not a malfunction. Avoid skin contact to these areas during these times.

On either sides of the headlight, one single battery can be used or charged. If two batteries are inserted, they are discharged or charged, respectively, at the same time. If only one battery is connected (active), the LED indicators flash.

If the batteries are being charged in the headlight, the USB charging unit requires approx. 6.5 h.

#### **Application**

- Insert the Micro-B USB connector (2) into the socket (with the marking ⚡) of the headlight (1).
- Connect the normal USB connector (3) either to the 5 V plug-in power supply (4), or to a USB port on a computer.

### 4.3.5 Charging the Batteries in the 4-Bay Desktop Charger

#### **NOTICE**

The headlight battery packs or 4-bay charger console may warm up during charge mode. This is normal and not a malfunction. Avoid skin contact to these areas during these times.

- The 4-bay desktop charger is intended for keeping fully charged replacement batteries ready for use.
- The batteries are fully charged after approx. 10 h.
- Position the console in a safe position on a table top that is protected from dangers of spillage and fire near a wall box.

#### **Application**

- Insert the charging cable.
- The red/orange LEDs light up.
- Insert the batteries into the chambers.  
Ensure that the batteries are properly aligned and fully inserted into the chamber.
- Indicator LED reference:

Indicator	Meaning
Green flashing light	Battery is being charged.
Green continuous light	Battery is fully charged.
Red/orange continuous light	No battery inserted or no contact. Fully re-insert the battery.
Red flashing light	Battery defective. Try another battery.

### 4.4 After Completion of the Application

#### **NOTICE**

Even improper depositing can cause damage to the products. Care should therefore be taken when handling the product.

Deposit the products properly and transfer to processing in accordance with chapter 5 "Cleaning, Disinfection and Sterilization", page 44.

---

## 5 Cleaning, Disinfection and Sterilization

---

** WARNING****Possible danger to life of the patient as a result of non-disinfected handling!**

All products, including non-sterile products, must be cleaned and disinfected before each use.

---

- The operator/preparer is responsible for the cleaning, disinfecting and sterilizing of the used products. It is essential that national regulations, including restrictions, be observed.

### 5.1 Limitations and Restrictions to Cleaning, Disinfecting and Sterilizing

The product lifetime is primarily determined by wear and possible tear during use.

A general statement with respect to a maximum of cycles is therefore not possible.

### 5.2 Pre-Treatment at the Application Location prior to Cleaning

Conduct the cleaning and disinfecting of contaminated products as soon as possible after use.

### 5.3 Preparation before Cleaning

---

** WARNING****Possible danger to life due to electric shock!**

Before cleaning, disconnect the charging unit from the headlight and power supply and power off completely.

---

When choosing a cleaning detergent, the compatibility with the material as well as the suitability for and the efficiency in the cleaning of medical devices must be considered.

The manufacturer's instructions for the cleaning or combined cleaning and disinfecting agent regarding concentrations, temperature, and application time must be observed.

## 5.4 Cleaning and Disinfection

### 5.4.1 Manual Cleaning and Disinfection

#### **NOTICE**

The optical areas of the headlights must be treated with the utmost caution to avoid risks.

Ensure that you wipe the equipment with a damp non-abrasive cloth during cleaning and that no parts of the product are sprayed.

If too much liquid is sprayed onto a cloth or used, it may permeate areas which could lead to damage to batteries, circuits and optics.

Always remove the mains plug for the headlight before cleaning the product.

#### **Lens cleaning:**

To clean the lens, clean as needed with a soft microfiber lens cloth and/or a lens cloth moistened with optical lens cleaner.

#### **Headlight cleaning:**

The device casing can be cleaned using a soft cloth that has been moistened with alcohol or a non-aggressive disinfectant solution.

- Use cleaning cloths or a lint-free cloth combined with a disinfectant based on isopropyl alcohol or simply 70 % isopropyl alcohol.
- Use a dry, lint-free cloth to remove residue and fingerprints. Never use paper fiber materials such as a paper tissue or a paper towel, as these can leave scratches.
- Wring out the cloth before use.
- Do not immerse in or spray with disinfection agents.
- Avoid excessive moisture, as this can lead to damage on the electrics, safety risks and voidance of the warranty.
- When choosing the cleaning agents and disinfectants to be used, it must be ensured that they are generally suitable for the cleaning and disinfecting of medical devices made from metals and plastics.
- Only approved cleaning solutions may be used, as there is otherwise a risk of damaging the equipment and voiding the warranty.

The black padding strips are cleaned using the same method. The padding strips can be removed for cleaning. However, we recommend leaving them on the headlight for cleaning. New padding strips are available as accessories.

Additional occupational safety measures (e.g. protective clothing, protective goggles, gloves; air filtering) in line with national specifications must be properly performed during manual cleaning/disinfection whenever there is a possible risk of injury or infection.

#### 5.4.2 Machine Cleaning and Disinfection

Automated cleaning and disinfection is not possible.

### 5.5 Inspection, Functional Check, Maintenance

#### 5.5.1 Inspection and Functional Check

---

**NOTICE**

Further use is confirmed by the successful inspection of the product.

---

- The products must be macroscopically clean, i.e. free from visible contamination, after each cleaning process.
- Inspect the products for fractures, cracks, deformation, damage and proper function.
- Worn, deformed, porous or otherwise damaged products must be replaced.

#### 5.5.2 Maintenance

n/a

### 5.6 Packaging

n/a

### 5.7 Sterilization

Cannot be sterilized.

### 5.8 Storage and Transport

- Store and transport cleanly.
- Avoid severe temperature fluctuations.
- Protect the products against mechanical damage.
- If the headlight is to be stored for more than two weeks, remove the batteries and place them in the original packaging.
- Store batteries only at ambient temperatures from 0°C to +25°C (+32°F to +77°F).
- Use the charging unit only in dry environments in a closed area. Protect from rain, snow, inflammable materials, excessive heat and/or moisture.

Where applicable, the specific parameters for ambient conditions during the application (temperature range and air humidity) as stated on the product label must be observed.

## 6 Maintenance

### 6.1 General Notices

The product may be repaired only by Gebrüder Martin or a qualified person or company expressly authorized by Gebrüder Martin to perform such work.

Modification of the product may lead to unforeseen risks and is therefore not permitted.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product are never allowed and shall void any claims of liability against Gebrüder Martin.

### 6.2 Regular Checks to Be Carried Out by the Operator

The product must be checked for proper functioning and operational safety prior to each use, including a visual inspection for mechanical damage that might impair the function:

- Function according to section 4.3 "Application", page 38.
- Intactness of the batteries

#### 6.2.1 Safety Check (SC)

n/a

### 6.3 Maintenance by the Manufacturer / Periodic Maintenance

The headlight itself is maintenance-free.

### 6.4 Format of the Serial Number

As a matter of principle, indicate the complete serial number in all cases of queries and complaints.

The serial number is indicated on the adjustment button of the headband (MedLED® Chrome):

	<p>The serial number is structured as follows:</p> <table border="1"><tr><td><b>MC</b></td><td>Gebrüder Martin MedLED® headlight</td></tr><tr><td><b>150520</b></td><td>Date (YYMMDD)</td></tr><tr><td><b>700411</b></td><td>Consecutive number per year of manufacture *</td></tr></table> <p>* Each year of manufacture starts with "0001" again.</p>	<b>MC</b>	Gebrüder Martin MedLED® headlight	<b>150520</b>	Date (YYMMDD)	<b>700411</b>	Consecutive number per year of manufacture *
<b>MC</b>	Gebrüder Martin MedLED® headlight						
<b>150520</b>	Date (YYMMDD)						
<b>700411</b>	Consecutive number per year of manufacture *						

## 6.5 Headlights, Accessories & Spare Parts

Designation	REF
<b>Headlights</b>	
MedLED Chrome MC3, Hardtop, cpl.	17-500-40-04
MedLED Chrome MC3, Softtop, cpl.	17-500-42-04
MedLED Chrome MC6, Hardtop, cpl.	17-500-50-04
MedLED Chrome MC6, Softtop, cpl.	17-500-52-04
MedLED Chrome MC7 PRO, Hardtop, cpl.	17-500-60-04
MedLED Chrome MC7 PRO, Softtop, cpl.	17-500-62-04
MedLED Chrome MC3/6/7 PRO, Hardtop, cpl.	17-500-70-04
MedLED Chrome MC3/6/7 PRO, Softtop, cpl.	17-500-72-04
<b>Accessories</b>	
Power Plug EU f. MedLED	17-500-01-04
Power Plug US f. MedLED	17-500-02-04
Power Plug UK f. MedLED	17-500-03-04
Power Plug AU f. MedLED	17-500-04-04
Magnifier, 2D, Sapphire/Onyx, MC3	17-500-12-04
Magnifier, 3D, Sapphire/Onyx, MC3	17-500-13-04
Magnifier, 4D, Sapphire/Onyx, MC3	17-500-14-04
Rechar. Station f. 4 Batt., MedLED	17-500-17-04
MedLED Chrome, USB Cable (only)	17-500-19-04
MedLED Chrome, Spare batteries, 2pc.	17-500-21-04
MedLED Chrome, Magnifier, 2dpt.	17-500-22-04
MedLED Chrome, Magnifier, 3dpt.	17-500-23-04
MedLED Chrome, Magnifier, 4dpt.	17-500-24-04
MedLED Chrome, Case, alone	17-500-29-04
MedLED Face Shield, for Headlight	17-500-30-04
MedLED Chrome, GoPro adapter	17-500-31-04
MedLED Onyx/MC3, amber filter	17-500-34-04
MedLED MC6, amber filter	17-500-35-04
<b>Spare Parts</b>	
MedLED Chrome, Headstrap, HT	17-500-32-04
MedLED Chrome, Headstrap, ST	17-500-33-04
MedLED Chrome, spare head band, HT, NCV	17-500-36-04

Designation	REF
MedLED Chrome, spare head band, ST, NCV	17-500-37-04
Bezel MC3	17-500-41-04
Bezel MC6	17-500-51-04
Bezel MC7 PRO	17-500-61-04

## 7 Ecological Information / Disposal

### 7.1 Packaging

Gebrüder Martin will take back the full packaging upon request. Whenever possible, parts of the packaging will be reused.

If this does not apply, the packaging can be disposed of in the regular paper- and ordinary municipal waste.


### 7.2 Batteries

Immediately and properly dispose of batteries after the end of their service life. Depending on your region, lithium batteries can be disposed of in normal household waste or industrial waste. Please dispose of batteries responsibly after the end of their service life.

### 7.3 Disposal of the Product

This product may not be disposed of as normal household waste. This product must be disposed of separately as an electronic device in accordance with the local regulations.

We therefore offer you to take the product back for proper disposal and recycling.

**Marking of electric and electronic equipment in accordance with Directive 2012/19/EC (WEEE Directive) and the German Electrical and Electronic Equipment Act (ElektroG)**

This symbol on the product or its packaging indicates that the product may not be disposed of as normal household garbage.

### 7.4 National Regulations

The national regulations and disposal provisions must be observed for all disposal measures.

## 8 Technical Data

Designation	Headlight MedLED® Chrome
<b>Headlight</b>	
Brightness	MC3: 100.000 lux MC6: 200.000 lux MC7 PRO: 275.000 lux
Weight	MC3: 335 g MC6: 335 g MC7 PRO: 350 g
Spot size	MC3: 10-100 mm MC6: 70-90 mm MC7 PRO: 100-200 mm
Light source	High-performance diodes
<b>LED</b>	
Service life	Up to 100,000 h. Replacement of the LED light source is not required
Color temperature	5.000 Kelvin
<b>Batteries</b>	
Ambient temperature (charge/use)	0 °C to +30 °C (+32°F to +86°F)
Humidity (use)	45–85 % relative humidity
Storage temperature (storage)	0 °C to +25 °C (+32°F to +77°F)
Charging time while in the headlight	6,5 h
Charging time in the 4-bay desktop charger	10 h
Total service life	More than 500 charging/discharging cycles; replace every 12 months
<b>Conformity to standards</b>	
USB plug-in power supply	IEC 61000
Battery	UL 2054, EN 55014-2:1997+A1:2001, meets the requirements of the battery directive
LED	IEC 62471:2006 – group “without conditions”
<b>Other information</b>	
Application part	Type A
Protection class	Class II
Maximum current intensity	1800 mA

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### **Martin Nederland/Marned B.V.**
### **Gebrüder Martin GmbH & Co. KG**
### **KLS Martin Taiwan Ltd.**
### **KLS Martin UK Ltd.**