E

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GB

Manufacturer Information on Reprocessing of Medical Devices according to DIN ISO 17664:2004

Manufacturer:	HAGER & WERKEN GmbH & Co. KG Ackerstraße 1, 47269 Duisburg, Germany T +49 (203) 99269-0 F +49 (203) 299283 www.hagerwerken.de	valid as of 07/2014 (Rev: 2)	Product: Miraclean Implant Black REF 355 220 – REF 355 224
---------------	---	------------------------------	--

Warning notices: No special notices

Limitation of reprocessing: max. 20 x

Place of handling: Clean surface with a disposable cloth or paper towel.

Storage and transport: No special requirements.

Preparations for cleaning: None

Cleaning: Use cleaning agent suitable for stainless steel according to the manufacturer's instructions. After cleaning, remove cleanser thoroughly under running water. When disinfecting mechanically do not exceed 93 °C

Disinfection: Use disinfecting solution suitable for stainless steel according to the manufacturer's instructions. Remove disinfecting solution thoroughly under running water. When disinfecting mechanically do not exceed 93 °C.

Maintenance: No special requirements

Packaging: Standardized packing material for sterilisation can be used (recommendation).

Sterilisation: Vapor sterilisation at 134 °C for 5 minutes or vapor sterilisation at 121 °C for 15 minutes.

Control/functional check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The above-mentioned instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine controls are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.