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Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

Manufacturer: HÄGER & WERKEN GmbH & Co. KG Ackerstraße 1 47269 Duisburg, Germany T +49 (203) 99269-0 F +49 (203) 299283 www.hagerwerken.de	valid as of: 02/2011 (Rev:2) Product: Mirahold-Block	REF 605 210 REF 605 238 REF 605 220
---	---	--

Warning notices: no special notices
--

Limitation of reprocessing: 10 x

Place of Handling: Clean surface with a disposable cloth or paper towel.

Storage and Transport: No special requirements.

Preparations for Cleaning: Separate lip retractor from biting block.

Cleaning:

Only use cleansers suitable for synthetics/plastic according to the instructions of the respective manufacturer. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93° C for mechanical cleaning.

Disinfection:

Use disinfecting solution suitable for synthetics/plastic according to the instructions of the respective manufacturer. Remove disinfecting solution thoroughly under running water. Do not exceed 93° C for mechanical cleaning.

Maintenance: No special requirements.

Packaging: Standardized packing material for sterilisation can be used (recommendation).

Sterilization:

Vapor sterilization at 134° C for 5 minutes or vapor sterilization at 121° C for 15 minutes.

Control/Functional Check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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GB Instructions for Use

- Spring bite blocks for mouth opening
- Silicone block is positioned between upper and lower jaw and keeps the mouth reliably open by spring force
- Fixation possible by cheek retractor
- Autoclavable

Bite-Block MIRAHOLD-BLOCK for mouth opening with bite block and retractor. The bite block is used as shown in pictures 1 – 4.

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Mirahold-Block Intro-Set **REF 605 238**

Mirahold-Block Set **REF 605 210**

For other languages of the Manufacturer Information on Reprocessing of Medical Devices see www.hagerwerken.de.

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Fabricante: HÄGER & WERKEN GmbH & Co. KG Versión: 02/2011 (Rev: 2)