Mirawhite® Laser Bleaching

Gebruiksaanleitung

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Bestimmen Sie den tooth colour by comparing it with the Vita® Shade Guide. To evaluate the colours, please re-sort the Vita® Shade Guide as follows:

Values sorted by colour hue: A1 A2 A3 A3,5 A4 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4
Sorted by brightness value: B1 A1 B2 D2 A2 C1 C2 D3 A3 D4 B3 A3,5 B4 C3 A4 C4

B2 – B1 = Values are increased by twofold / A3 – A2 = Values are increased by fourfold

Use air-flow or a cleaning paste to clean the patient's teeth. Use a cheek/lip retractor or similar instrument. Make sure to protect the gingiva (e.g. with Gingiva Protector or by applying a rubber dam). Do not anaesthetise the patient. **The mixed bleaching gel irritates the skin and mucous membranes, and causes severe eye damage. Always wear protective gloves and safety glasses when using the material! Strictly avoid contact of the material with eyes, skin, and mucous membranes.**

Application with LaserHF or another diode laser

Prepare the laser system following the settings given in the manufacturer's handbook and turn on the laser in order to make sure it is functioning properly.

While facing away from the patient, start by removing the clamp. Then remove the cap and replace it with a mixing tip. Now squeeze the gel from the syringe until a homogeneous colour mass with gel-like consistency is coming out of the tip. Discard the initial inhomogeneous mass. Storage for long periods at high temperatures may sometimes cause bubbles to form in the gel. This can also lead to inhomogeneous and more liquid mixtures, particularly in the last quarter of the dual syringe. In such a case, do not use the affected parts of the material. The gel should be applied in a thickness of about 1 – 2 mm to the labial surfaces of the teeth to be bleached. The point of the mixing tip is well-suited to spread the gel evenly. **Make sure the patient wears protective glasses during the whole length of treatment. Please check the Laser Unit's manual for the length of treatment.** Following the laser treatment, aspirate the gel, wipe off the residue, and rinse thoroughly with water. Now dry the tooth surface and polish with polishing paste. Subsequently, a remineralisation with After-Bleaching-Care or with a fluoride solution or fluoride gel is recommended. This is also suitable for treating pain if sensitisation occurs. If there is still enough material for another bleaching treatment, remove the mixing tip and replace the blue cap. Store the material in a refrigerator. Use a new mixing tip for the next treatment. The patient should not eat, drink, or smoke for 1 hour after treatment. Substances with strong colours (e.g. red wine, coffee) and food containing acid should also be avoided in the days following the treatment.

Conventional application

Following a reaction time of about 10 minutes, aspirate the gel and wipe off any residue. Should the bleaching effect be insufficient, the procedure may be repeated twice until the desired effect has been achieved. **Prior to each treatment, ensure that the gingiva protection is still sufficient.** Then determine the tooth colour again. After the final bleaching procedure, aspirate the gel, wipe it off, and rinse thoroughly with water. Now dry the tooth surface and polish with polishing paste. Subsequently, a remineralisation with After-Bleaching-Care with a fluoride solution or fluoride gel is recommended. This is also suitable for treating pain if sensitisation occurs. If there is still enough material for another bleaching treatment, remove the mixing tip and replace the blue cap. Store the material in a refrigerator. Use a new mixing tip for the next treatment. The patient should not eat, drink, or smoke for 1 hour after treatment. Substances with strong colours (e.g. red wine, coffee) and food containing acid should also be avoided in the days following the treatment.

Note: Any composite fillings and restorative materials made of ceramic or metal are not bleached. Fillings should, therefore, be replaced 2 weeks after treatment.

Devital/endodontically treated teeth

The gel is additionally applied through a lingual embrasure. First, check the height of the alveolar ridge. Then remove the complete filling material from the pulp chamber. Apply a layer of glass ionomer filling with a thickness of 1-2 mm to the bottom of the pulp chamber on either the alveolar ridge level or coronally. Apply the gingiva protection to both, the labial and the lingual sides of the tooth. Then apply the bleaching gel to the labial surface as described above under "Conventional application and additionally through the access opening. Let the gel act for 10 minutes. Afterwards, carefully aspirate the gel and lightly rinse with water. If the bleaching effect is insufficient, this procedure can also be repeated twice until the desired effect has been achieved. Following the final bleaching process, aspirate the gel and rinse thoroughly with water. Carefully remove the gingiva protection. Then, close the access opening with a temporary filling material free from eugenol. Wait 2 weeks before preparing the permanent filling.

Shelf life / storage

As an adequate bleaching effect can no longer be guaranteed, **Mirawhite Laser Bleaching** must not be used past its expiration date. **Refrigerate the bleaching gel (3 – 9 °C), but do not freeze! The product contains hydrogen peroxide, potassium nitrate, sodium fluoride, sodium hydroxide, glycol derivatives, thickening agent and coloring.**

Medical device

Troubleshooting

Insufficient bleaching	Repeat the bleaching process (several times if necessary)
Patient is experiencing pain	The gel has come into contact with exposed dentine or a loose filling (please cover before treatment). The patient is particularly sensitive (depending on the individual tooth structure).
White spots on the enamel following the bleaching process	Local demineralisation or dehydration. Reversible. Disappears after fluoridation.
Gingiva is coloured white after bleaching (procedure can be reversed)	The gingiva has not been protected adequately, or the protection is loose.

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Mirawhite® Laser Bleaching

Instructions for use

Mirawhite Laser Bleaching is an in-office bleaching gel intended for medical treatment in dental practices only. The material is used by the dentist to brighten individual teeth or devital teeth including endodontically treated ones. The gel is suitable for use with diode laser systems with a wavelength between 810 and 980 µm. It contains bleaching gel and an activator and comes in 4:1 self mixing syringes (cartridge). One syringe contains 5 g of material. The hydroperoxide content of the ready-mixed gel is 35 %.

The set includes:

1 syringe 5 g / 3 mixing tips / 1 syringe Gingiva-Protector 1.5 g / 1 syringe After Bleaching Care 2 ml / 3 luer tips for individual syringe.

Indications

- Brightening of discoloured, devital teeth
- Brightening of teeth discoloured due to illness, injury, as well as iatrogenic teeth as an alternative to a crown or veneer
- Treatment of tooth discolouration caused by prescription drugs (i.e. tetracycline)

Contraindications

- Exposed dentine
- Defective/ faulty filling margins
- Hypersensitive patients
- Allergy to or intolerance of the ingredients
- Pregnancy and lactation (due to lack of experience)

May cause slight tooth sensitisation in some patients during or shortly after treatment. It will usually disappear after a short while. **Stop treatment immediately if pain out of proportion is being experienced!**

The gel from the cartridge has an inhomogeneous colour (colourless / red)	There are air bubbles in the gel, resulting from being stored too long or at temperatures too high.
Reduced composite-enamel adhesion following the bleaching process	Peroxide residue. Wait 1 – 2 weeks before restoration.

 GDF GmbH, Dieselstr. 5 – 6, D-61191 Rosbach, www.gdfmbh.com
Distribution: Hager & Werken GmbH & Co.KG – Ackerstr. 1 – D 47269 Duisburg – www.hagerwerken.de

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