PERMA SHARP SUTURES

PERMA SHARP SUTURES
Distributed by:
Hu-Friedig Mfg. Co., LLC
3322 N. Rockwell St.
Chicago, IL 60618 USA
www.hu-friedig.com

CE 0120
EN 03-012102-1018

DESCRIPTION
PGA™ FA (Polyglycolide) is a synthetic, resorbable, sterile surgical suture composed of a homopolymer of glycolic acid.

INDICATIONS
This suture is indicated only for use in soft tissue approximation of the skin and mucosa, where only short term wound support (7-10 days) is required.

CONTRAINDICATIONS
Do not use this suture where extended periods of tissues under stress is required or where wound support beyond 7 days is required.

WARNINGS
Do not resterilize. Discard open, unused sutures and associated surgical needles.

PRECAUTIONS
Care should be taken to avoid damage when handling. Avoid crushing or crimping the material with surgical instruments such as needle holders and forceps.

ADVERSE REACTIONS
Adverse effects associated with the use of this device may include wound dehiscence, failure to provide adequate wound support in closure of sites where expansion, stretching or distension occur, etc.

STERILITY
PGA™ FA (Polyglycolide) synthetic absorbable sutures are sterilized using gamma radiation. Do not resterilize. Do not use if package is opened or damaged.

HOW SUPPLIED
PGA™ FA (Polyglycolide) synthetic absorbable surgical sutures are available sterile as braided multifilament, uncoated (clear), coated strands in USP sizes 5-0 to 2-0 (metric sizes 1.0 through 3.0) in a variety of lengths adapted to various needs.

SYMBOLS USED FOR LABELING
Product Code
Caution - Refer to Accompanying Documents
Do Not Reuse
Use by
Sterilized using irradiation
Do Not Resterilize
Lot Code
Do not use if package is damaged

CAUTION
Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.
CE mark and identification number of notified body. Product conforms to the essential requirements of The Medical Device Directive 93/42/EEC.

Authorized Representative in the European Community
Pearsalls Limited
Tancred Street, Taunton
Somerset, TA1 1RY, England

DESCRIPTION
PGA™ FA (Polyglycolic Acid) suture is a synthetic absorbable, sterile, surgical suture composed of a copolymer of glycolic acid and L-lactide.

INDICATIONS
PGA™ FA (Polyglycolic Acid) synthetic absorbable suture is indicated only for use in soft tissue approximation of the skin and mucosa, where only short term wound support (7-10 days) is required.

CONTRAINDICATIONS
Do not use this suture where extended periods of tissues under stress is required or where wound support beyond 7 days is required.

WARNINGS
Do not resterilize. Discard open, unused sutures and associated surgical needles.

PRECAUTIONS
Care should be taken to avoid damage when handling. Avoid crushing or crimping the material with surgical instruments such as needle holders and forceps.

ADVERSE REACTIONS
Adverse effects associated with the use of this device may include wound dehiscence, failure to provide adequate wound support in closure of sites where expansion, stretching or distension occur, etc.

STERILITY
PGA™ FA (Polyglycolic Acid) synthetic absorbable sutures are sterilized using gamma radiation. Do not resterilize. Do not use if package is opened or damaged.

HOW SUPPLIED
PGA™ FA (Polyglycolic Acid) synthetic absorbable surgical sutures are available sterile as braided multifilament, uncoated (clear), coated strands in USP sizes 5-0 to 2-0 (metric sizes 1.0 through 3.0) in a variety of lengths adapted to various needs.

SYMBOLS USED FOR LABELING
Product Code
Caution - Refer to Accompanying Documents
Do Not Reuse
Use by
Sterilized using irradiation
Do Not Resterilize
Lot Code
Do not use if package is damaged

CAUTION
Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.
CE mark and identification number of notified body. Product conforms to the essential requirements of The Medical Device Directive 93/42/EEC.

Authorized Representative in the European Community
Pearsalls Limited
Tancred Street, Taunton
Somerset, TA1 1RY, England

Authorized Representative in the European Community
Pearsalls Limited
Tancred Street, Taunton
Somerset, TA1 1RY, England

Authorized Representative in the European Community
Pearsalls Limited
Tancred Street, Taunton
Somerset, TA1 1RY, England

Authorized Representative in the European Community
Pearsalls Limited
Tancred Street, Taunton
Somerset, TA1 1RY, England

Authorized Representative in the European Community
Pearsalls Limited
Tancred Street, Taunton
Somerset, TA1 1RY, England

Authorized Representative in the European Community
Pearsalls Limited
Tancred Street, Taunton
Somerset, TA1 1RY, England

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