Deutsch
English
Español
Français
Italiano
Português

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Table of Contents

1	General.....	32
1.1	Manufacturer	32
1.2	Hotline	32
1.3	Adverse incident reporting requirement	32
1.4	Summary of safety and clinical performance	32
1.5	About this document	33
1.6	Abbreviations and Terms	33
1.7	Validity of this Document.....	33
2	Scope of Delivery.....	34
2.1	Checking the delivery for completeness and correctness	34
3	Intended use	34
3.1	Intended purpose.....	34
3.2	Indications.....	34
3.3	Contraindications	34
3.4	Possible adverse effects	34
3.5	Residual risks.....	34
3.6	Patient target group.....	35
3.7	Users	35
3.8	Environmental conditions during use.....	35
3.9	Restrictions on use.....	35
3.10	Warnings.....	36
4	Application	43
4.1	Description of the components.....	43
4.1.1	Structure, Functionality and Performance Characteristics.....	43
4.1.2	Variants/Combinations	44
4.1.3	Combination products and accessories.....	44
4.2	Before First Use and Before Each Further Use	44
4.3	Application / Intraoperative Procedure	45
4.4	After Completion of the Application	45
5	Cleaning, disinfection, and sterilization.....	46
5.1	Limitations and restrictions on cleaning, disinfection, and sterilization.....	46
5.2	Pre-Treatment at the Application Location prior to Cleaning	47
5.3	Preparation before Cleaning	47
5.4	Cleaning and Disinfection.....	48
5.4.1	Detergents and Disinfectants.....	48
5.4.2	Machine Processing	48
5.4.3	Manual Processing.....	49
5.4.4	Final Rinsing	49

5.5	Inspection, Functional Check, Maintenance	49
5.5.1	NE with Non-Split Conductive Surface	49
5.5.2	NE with Split Conductive Surface	49
5.6	Packaging	49
5.7	Sterilization	50
5.7.1	Before Sterilization	50
5.8	Storage and transport	51
6	Maintenance	51
6.1	General Notes	51
7	Environmental information/disposal	52
7.1	Packaging	52
7.2	Disposal.....	52
7.3	National regulations.....	52

1 General

1.1 Manufacturer

Thank you for choosing a KLS Martin product.

This product has the CE mark, which means that it is compliant with the essential safety and performance requirements for medical devices pursuant to applicable European regulations.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group
KLS Martin Platz 1 · 78532 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

1.2 Hotline

If you have questions about how to use this device or product or questions about clinical applications, please contact the product management team:

Tel: +49 7461 706-243

Fax: +49 7461 706-190

If you have technical questions or questions about maintenance contracts and training, please contact the Martin Service Center:

Tel: +49 7461 706-343

Fax: +49 7461 706-484

E-Mail: service@klsmartin.com

NOTICE

All packaging, and in some cases also the product itself, is marked with a batch number (LOT) and a catalogue number (REF). When making a product complaint, please include the LOT and REF.

1.3 Adverse incident reporting requirement

All serious incidents occurring in connection with the product shall be reported immediately to Gebrüder Martin and to the competent authority.

1.4 Summary of safety and clinical performance

The Summary of Safety and Clinical Performance can be downloaded from the European Database for Medical Devices (EUDAMED).

1.5 About this document

Could result in the death of or serious injury to the patient, user or a third party if these instructions are not followed

Read and follow all of the instructions for use. In particular, take all precautionary measures and heed all warnings.

This text refers to male, female and non-binary individuals equally. but not all pronouns have been used, solely for reasons of readability.

The electronic version of these instructions for use can be downloaded from www.klsmartin.com.

1.6 Abbreviations and Terms

Abbreviation	Description
MRI	Magnetic resonance imaging
OR	operating room
W/D	Washer-disinfector

1.7 Validity of this Document

These Instructions for Use apply to the following products:

Item no.	Name	Cable	Comments
80-310-04-04	Hand cylinder neutral electrode	4 m	Hand cylinder NE as grounding aid for monoterminial applications with the MD62
80-332-03-04	Rubber neutral electrode 8 cm x 16 cm	4 m	
80-342-05-04	Rubber neutral electrode 8 cm x 16 cm	4 m	
80-342-06-04	Twin-Pad neutral electrode 15 cm x 25 cm	4 m	2 x item no. 80-371-02-04 required!
80-342-07-04	Twin-Pad II neutral electrode		

2 Scope of Delivery

Name	REF
1x of following depending on order:	
80-310-04-04	Hand cylinder neutral electrode
80-332-03-04	Rubber neutral electrode 8 cm x 16 cm
80-342-05-04	Rubber neutral electrode 8 cm x 16 cm
80-342-06-04	Twin-Pad neutral electrode 15 cm x 25 cm
80-342-07-04	Twin-Pad II neutral electrode
these Instructions for Use	

2.1 Checking the delivery for completeness and correctness

Immediately after receiving the delivery, check that it is complete and in good condition.

Report any damage that occurs during transit promptly.

After delivery of the product, check the original packaging and the packaging seal for intactness.

3 Intended use

3.1 Intended purpose

The neutral electrodes by Gebrüder Martin are designed for connection to electrosurgical units by Gebrüder Martin and are intended for closing an HF circuit to be used for monopolar coagulation and monopolar cutting in biological tissue.

3.2 Indications

Not applicable.

3.3 Contraindications

Not applicable.

3.4 Possible adverse effects

Not applicable.

3.5 Residual risks

No product-specific risks are known when used properly. The user is responsible for informing the patient about possible risks and complications.

The user is responsible for informing the patient about existing residual risks in connection with the use of this product.

3.6 Patient target group

There are no restrictions regarding the patient target group.

3.7 Users

The products are to be used exclusively by qualified surgical staff.

Cleaning, disinfection and sterilization are carried out by trained qualified staff in the processing unit for medical devices.

3.8 Environmental conditions during use

Use only in ORs under conditions appropriate for surgery or in medical facilities designed for this.

3.9 Restrictions on use

Danger of injury in case of insufficient qualification of the user!

Instruments for electrosurgery may be used only by persons who have been specially trained or instructed in their use!

3.10 Warnings

 **WARNING**

Peril to life or danger of serious physical injury due to prion infection!

Neutral electrodes that may have come in contact with prions must be destroyed at once (danger of CJD contamination)! Never reuse them in such case!

 **WARNING**

Risk of infection in case of non-sterile handling!

Improper sterilization and non-sterile use of the neutral electrodes can pose severe health hazards to patients.

Sterilization must be carried out according to a validated steam sterilization process, for example in a sterilizer complying with EN 285:2009 and DIN EN 13060, validated in accordance with ISO 17665-1:2006.

NOTICE

Risk of neutral electrode damage!

- Do not use hot air for neutral electrode sterilization!
 - Do not use the “STERRAD®” system for neutral electrode sterilization! (STERRAD is a registered trademark of Johnson & Johnson Inc.).
 - Do not use “flash” autoclaving!
 - Do not use irradiation sterilization, gas sterilization, or formaldehyde sterilization!
-
- Note that high-temperature sterilization and long cycle times reduce the life of the neutral electrodes.
 - Ensure sufficient drying.
 - Follow the sterilization recommendations of the manufacturer of your sterilizer regarding handling, loading and parameter setting.
 - The operator must ensure that sterility is maintained after the sterilization process.
 - Only use the sterilization methods listed below. Do not use other sterilization procedures.

Method	Exposure temperature	Exposure time	Drying time
Fractionated vacuum (dynamic air removal)	min. 132 °C	min. 3 min	min. 30 min

 **CAUTION**

Danger of burns!

The products are delivered ex factory slightly twisted due to the manufacturing process. On level ground, the distance between the contact surface and the NE should not exceed 1.5 cm.

Significantly twisted neutral electrodes (>1.5 cm) or those with hardened silicone pads are faulty and can cause burns.

- Make sure that there is a distance of <1.5 cm between the contact surface and the highest point of the neutral electrode.

 **WARNING**

Danger of burns in case on non-observance of these instructions!

Non-observance of these instructions can lead to user or patient injury!

-
- A **non-split** cannot be monitored. In case of insufficient contact, **no warning signal** will be triggered! We therefore generally recommend use of a split disposable neutral electrode.
 - The maximum output power specified above applies to normal application times of a few seconds. Intermittent cooling phases may be required.
 - Avoid output power settings on the electrosurgical unit where the maximum output voltage might lead to exceedance of the rated voltage of the accessories used.
 - The temperatures under the neutral electrode can rise above permissible increase of 6 K if high currents (>700 mA) are transmitted via the NE for ≥ 60 seconds. The danger of burns increases. In such a case, current and application duration should be minimized immediately.
 - Contact quality monitoring is possible only when using a split neutral electrode. However, intraoperative checks are required in these cases as well!
 - The hand cylinder neutral electrode (item no. 80-310-04-04) is an exception. The neutral electrode may be used only for monoterminal applications with the MD62. The patient is conscious and holds the NE in his/her hand. The palms of the hands must be cleaned.
 - When selecting a monitorable neutral electrode, make sure that it is compatible with the contact quality monitor of the electrosurgical unit used.
 - In order to create safety reserves, please select a neutral electrode with maximum surface and with contact quality monitoring. The selection depends on the possible application site and the electrosurgical unit used.
 - Be sure to observe the safety instructions provided in the user manual of the electrosurgical unit, as well as the instructions received during training!

- Neutral electrodes by KLS Martin are delivered in an unprocessed condition and are not packed in sterile packaging. The neutral electrode must be processed before use. At least disinfection with silicone-compatible agents should be carried out in order to avoid cross-contamination between patients. Hot steam sterilization at 134 °C with fractionated vacuum is possible for all the items mentioned here. Deviating reprocessing procedures must be avoided.
- Inspect all accessories regularly, especially the electrode connection cables and endoscopically used active accessories, for possible damage (e.g. with a magnifying glass). Damaged accessories must not be used!
- Unlike disposable neutral electrodes, reusable neutral electrodes do not have their own adhesive force. Therefore, mechanical attachment must be provided, preferably on the thigh (or the upper arm) by inserting the neutral electrode under a support stocking according to compression class RAL-GZ 387 2000-09, class 1=KKL 1, or by attachment with a Flexoversal bandage. Be aware that the pressure applied to the neutral electrode – and thus to the well-vascularized tissue underneath – should be uniformly distributed across the entire NE or skin surface.
- There are risks that can be caused by neuromuscular stimulation. This can particularly occur in operating modes that generate arcs between the active electrode and the tissue.
- **Always hold the connector when disconnecting the HF cable. Never pull on the HF cable, as this will cause damage to the cable which is difficult to detect and will result in premature failure of the product.**
- Insufficient contact between the neutral electrode and the patient's skin will be indicated by an audible alarm only if you use a compatible, monitorable neutral electrode such as item no. 80-342-06-04 together with a contact quality monitor (e.g. called PCS system or NE monitor).
- To prevent burns, regular intraoperative checks are required, especially when neutral electrodes with a non-split conductive surface are being used.
- Care must be taken when handling reusable neutral electrodes, as bending of the conductive surface, cable connections and tight wrapping with the cable can significantly reduce their function and service life. Furthermore, there is always a danger of tripping when handling cables. For avoidance of functional defects or burns, the items must not be subjected to mechanical stress during storage (e.g. pressure from above) either.

⚠ WARNING

Risk of injury!

Non-observance of the following notices can lead to death or serious injury!

- **Pacemakers:**
Can be disturbed or damaged by RF current application. Consult a cardiologist before performing the surgery.
- **Metal shunts in RF current path:**
e.g., implants can increase current density to the extent of causing unnoticed endogenous burns.
- **Cables run in parallel:**
Interference can be induced into cables laid out parallel to the RF cable, thus causing e.g. monitor image troubles.
- **Grounded objects:**
e.g., if patient insulation is insufficient, unwanted leakage currents can be increased by contact with an operating table in a way that can compromise the proper functioning of the unit or lead to burns.
- **Disinfectants and detergents:**
Cleaning or disinfecting agents that have not been allowed to dry off sufficiently can cause burns.

NOTICE

Improper processing can lead to damage!

Non-observance of the processing instructions can destroy the neutral electrode or lead to damage such as

- **loss of all electrical and mechanical design properties**
- **deformation of plastic parts**
- **increasing porosity of the soldered connection inside the plug**

NOTICE

Danger of burns!

Applies only to use of the neutral electrode TWIN-PAD II, REF 80-342-07-04: Incorrect application may result in burns.

- **When using the neutral electrode TWIN-PAD II, pre-treat the application surface with a damp cloth soaked in NaCl solution.**
- **Observe the application rules for neutral electrodes.**

⚠ WARNING**Danger of burns!**

- There is a danger of ignition of endogenous gases. If saturated with oxygen, some materials such as cotton and gauze can be ignited by the sparks generated during the proper use of the electrosurgical unit.
 - Cables with defective/broken insulation or broken conductors can cause user or patient burns or even fire. If a neutral electrode is damaged in any way, withdraw it from service.
 - Never use defective NEs! Non-compliance can lead to burns!
-

⚠ CAUTION**Danger of injury!**

Non-compliance with the following instructions can cause injury!

- The use of ignitable anesthetics and oxidizing gases, such as nitrous oxide (N₂O) and oxygen, should be avoided when performing an operation in the thoracic or head area, unless these substances are removed with an evacuator.
 - Wherever possible, non-flammable substances should be used for cleaning and disinfection. Flammable substances used as detergents or disinfectants or as solvents for glues should be allowed to evaporate prior to using electrosurgery.
 - Ground the operating table.
-

 **CAUTION**

Danger of injury!

Non-compliance with the following instructions can cause injury!

- To prevent monitor image trouble, avoid running RF cables closely parallel to camera cables. In any case, we recommend using monitoring systems incorporating devices for limitation of the high-frequency current.
- RF cables must never be allowed to get in touch with the patient's skin.
- Never run RF cables in loops.
- Place the patient on a dry and electrically insulated support.
- Skin-to-skin contact (e.g. between the arms and the body of the patient) should be avoided. For example, use dry gauze as insulating layers.
- To prevent unwanted tissue damage, the bipolar technique can be preferable for surgical interventions where the RF current may be expected to flow through body parts with a relatively small cross section.
- Use an output power as low as possible for the intended purpose.
- Select the place of application for the neutral electrode with care:
 - Use a muscular or well-vascularized, convex skin surface showing no anomalies (such as hair growth, punctures, inflammation, injuries, scars, warts, spots, bony prominences) in the vicinity of the surgical site, but at least 20 cm (preferably approx. 30 cm) away from it, for adults preferably on the upper arm or thigh, for children preferably on the trunk.
 - Make sure that the expected current paths inside the patient's body will be as short as possible. Never allow such paths to run through the entire body.
 - Never allow current paths to run across the heart.
 - Avoid places where liquids could accumulate. There is a danger of flammable liquids accumulating underneath the patient, in pits such as the navel or in body cavities such as the vagina. Any liquid that has accumulated in such places should be wiped off before using the electrosurgical unit.
 - The neutral electrode must be attached closer to the surgical site than any ECG electrode.
- Prepare the application site carefully:
 - shave evenly
 - clean (remove any grease, e.g. cosmetic substances)
 - drying.
- Attach the neutral electrode using a suitable means of fixation that ensures highly uniform pressure distribution across the entire skin surface.
- Full-surface application of the NE to the patient's skin is important! Preferably on the thigh by inserting it into a compression class I elastic stocking, or on the upper arm, or by using silicone bandages for fixation.
- Be aware that the pressure applied to the neutral electrode – and thus to the well-vascularized tissue underneath – should be uniformly distributed across the entire NE or skin surface.
- Run the NE connecting cable so that it touches neither the patient nor other cables. Avoid cable loops.
- Now connect the neutral electrode, the foot switch and the RF cable to the electrosurgical unit.
- Verify that the neutral electrode has been reliably attached to the patient's skin across its entire surface.

⚠ CAUTION**Risk of injury!**

- Repositioning a securely attached neutral electrode should be avoided, but is permitted. If necessary, proceed as follows:
 - Remove the NE
 - Select a suitable skin area and clean it
 - Reattach the NE as described above!
 - If you need to reposition the patient during the surgery, make sure afterwards that the neutral electrode is still securely attached across its entire surface!
-

NOTICE

- When using electrosurgical units in conjunction with handles featuring two separate finger switches for optional activation of the CUT and COAGULATION high-frequency operating modes, different settings can be used for each mode regarding type of current and output power. In this case, operating the CUT and COAG activation buttons during the normal use of the electrosurgical unit can lead to variations in the effective output power. This does not pose a hazard, however, and therefore is no reason to stop the required surgical procedure.
 - If the cutting or coagulation effect diminishes during an operation, be sure to check the neutral electrode for correct, full-surface skin contact before increasing the output power on your electrosurgical unit!
Besides, be sure to observe the safety instructions provided in the user manual of your electrosurgical unit!
-

⚠ WARNING**Danger of infection from contamination!**


Avoid skin contact during disposal! Use only clearly marked disposal containers.

4 Application

4.1 Description of the components

4.1.1 Structure, Functionality and Performance Characteristics

The reusable neutral electrodes by Gebrüder Martin are intended for connection to electrosurgical units by Gebrüder Martin. In addition, they are attached to the outside of the patient's skin. They are used to return the current that flows through biological tissue during monopolar coagulation and cutting to the electrosurgical unit. The patient is part of the circuit and is connected between the active electrode and the neutral electrode:

No.	Description
1	Patient
2	Active electrode
3	Flow field
4	Neutral electrode

The reusable neutral electrodes are attached to the patient by straps or inserted into a support stocking to ensure secure attachment and uniform pressure on the tissue.

Despite the transfer of electromagnetic energy within a direct spatial field around the contact surfaces of the neutral electrodes, the reusable neutral electrodes do not cause any physiological effects in the patient tissue. The field strength around the neutral electrode decreases strongly in its spatial extent with increasing distance from the neutral electrode.

The split neutral electrodes have an insulating clearance between the two contact surfaces. The insulating clearance enables monitoring of the quality of the skin contact using an NE monitor in the electrosurgical unit, as the two contact surfaces are separated from each other only by the skin and the lower tissue layers. The NE monitor operates in a much lower frequency range and at much lower voltages than the HF energy. The insulating clearance crosses the entire contact surface of the split neutral electrode and thus enables efficient monitoring of the quality of the skin contact. A change of the current flowing through the insulating clearance indicates a change in the quality of the skin contact. A change towards high impedance detected by the NE monitor indicates poor quality of the skin contact, and the NE monitor deactivates the HF energy. Split, reusable neutral electrodes and the NE monitor must be matched to each other. Inductors (e.g. L=180 µH) which were arranged in parallel above the conductive surfaces and thus, among other things, reduced the so-called run-in times, were also used for matching split neutral electrodes in the past. Twin-Pad II dispenses with this solution in favor of an application specification which specifies how the perspiration process can be supported.

4.1.2 Variants/Combinations

Several different variants of reusable neutral electrodes are available:

Item no.	Name	Cable	Comments
80-310-04-04	Hand cylinder neutral electrode	4 m	Hand cylinder NE as grounding aid for monoterminal applications with the MD62
80-332-03-04	Rubber neutral electrode 8 cm x 16 cm	4 m	
80-342-05-04	Rubber neutral electrode 8 cm x 16 cm	4 m	
80-342-06-04	Twin-Pad neutral electrode 15 cm x 25 cm	4 m	2 x item no. 80-371-02-04 required!
80-342-07-04	Twin-Pad II neutral electrode		

4.1.3 Combination products and accessories

A high-frequency electrosurgery unit with compatible plug connectors is required for all neutral electrodes.

An additional connection cable is needed for the "Twin-Pad II" neutral electrode, REF 80-342-07-04:

Item no.	Name
80-294-40-04	NE connection cable with KLS Martin connector
80-294-43-04	NE connection cable with Erbe connector
80-294-44-04	NE connection cable with international connector
80-372-00-04	Flexible fabric band with hook-and-loop fastening 90 cm for Twin-Pad II

4.2 Before First Use and Before Each Further Use

Carry out a functional check of the NE before the start of the operation:

- Non-split NE: Insert the device connector into the HF device.
NE monitor switches from red to green. When the cable is moved, the NE monitor remains green.

- Split NE: Insert the device connector into the HF device.
NE monitor does not switch from red to green. If the two conductive surfaces are short-circuited over a large area, e.g. after correct application to the patient and run-in time, the NE monitor switches from red to green. When the cable is moved, the NE monitor remains green.

4.3 Application / Intraoperative Procedure

- Apply the neutral electrode with fixing bandage or by insertion into a support stocking according to compression class RAL-GZ 387 2000-09, class 1=KKL 1.

As experience shows, unexpected output power changes pose an extremely high danger of burns!

If you use your standard setting on your electrosurgical unit, but you do not get the usual coagulation or cutting result, or the power output decreases during operation, or the electrosurgical unit even fails (with or without NE alarm), be sure to heed the following:

- Never increase the output power on the unit unless you have checked the surgical situation!
- In particular, check
 - all RF plugs and cables for proper contact and the foot switch or handle pushbutton for proper functioning
 - the insulation of all RF components for visible damage
 - the distal end of the active electrode(s) for cleanness and wear
 - in particular, the neutral electrode for correct application!
- Do not use conductive gel to improve skin contact. Instead, attach a new NE if necessary!
- Note that using conductive gel could lead to so-called “hot spots” due to uneven attachment of the NE. This is because the thickness of the gel layer applied would be uncontrollable and uniform drying of the gel might not be guaranteed.
- Never reduce the conductive surface by cutting it to size! If necessary, use a new, smaller NE! Cutting to size is not permitted because this could lead to exposure of the metal inserts, thus creating a risk of the bare metal (= shunt) coming in direct contact with the skin.
- Avoid repositioning the NE. If necessary, be sure to reattach it correctly!

4.4 After Completion of the Application

- Following use, simply unfasten and remove the neutral electrode, but avoid pulling at the connecting cable during the process.
- Do not roll up the NE or wrap the connecting cable around it. Non-compliance can lead to damage and premature failure of the product.

5 Cleaning, disinfection, and sterilization

Possible danger to life of the patient in case of non-sterile handling!

Products supplied in non-sterile condition must be cleaned, disinfected and sterilized before the first and before any further use.

- In case of patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible CJD variants, processing of the products must be carried out according to the applicable national regulations.
- The operator/processor is responsible for cleaning, disinfecting and sterilizing the products used. It is essential that national regulations, including restrictions, be observed.

The processing procedure must be validated prior to use of the medical device. The operator/processor is responsible for this.

Because of potentially inadequate labeling on the product, products that have been removed from their sterile packaging but have not been used and are not contaminated cannot be considered equivalent to products that have been obtained as non-sterile products. In this case, the operator/processor is responsible for use including labeling, cleaning, disinfection and sterilization.

The specific information for each reference number for the following section is available at:

www.klsmartin.com/processing

5.1 Limitations and restrictions on cleaning, disinfection, and sterilization

The product life of a surgical instrument is essentially determined by wear and possible damage during use.

This is recognized in the functional check and visual inspection to be performed in each individual processing cycle. If required, appropriate corrective measure shall be taken. This functional check and visual inspection are the basic prerequisite for the release of the medical device for the next processing step, for the storage of the medical device in the validated sterile barrier system and for the application of the medical device after sterilization.

For this reason, a general specification of maximum cycles is not possible.

5.2 Pre-Treatment at the Application Location prior to Cleaning

Conduct the cleaning and disinfecting of contaminated instruments as soon as possible after use.

The following generally applies:

The sooner the inner lumina are cleaned, the better the result will be. If timely cleaning is not possible due to the duration of the application or as a result of organizational aspects, the users must determine and validate measures on their own responsibility to prevent the contamination from drying completely or to ensure successful cleaning despite drying.

We generally recommend dry return transport in a closed transport container without the addition of liquids, cleaning agents or disinfectants.

5.3 Preparation before Cleaning

Possible danger of injury to the patient from residues on the medical device!

Particles from cloths and brushes may get stuck on products with rough surfaces, threads, sharp edges, narrow gaps or the like.

- Use only soft, clean and lint-free cloths for drying.
- The brushes for ducts must be a little bit bigger than the inner diameter of the duct, and the shaft of the brush must be at least as long as the duct.

A machine-based process (W/D) must be used for cleaning and disinfecting. The information whether the product has been validated for a manual procedure can be found at www.klsmartin.com/processing. Due to the significantly lower effectiveness and reproducibility, a manual procedure – also using an ultrasonic bath – should only be used if a mechanical procedure is not available.

When selecting the cleaning agent used, the material compatibility, suitability and effectiveness for cleaning medical devices must be taken into account.

The concentrations, temperatures and exposure times specified by the manufacturer of the cleaning agent and/or disinfectant as well as the specifications for rinsing must be observed. Use only freshly prepared solutions, only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water) or only a soft, clean and lint-free cloth and/or filtered air for drying.

Store the medical devices on carriers suitable for cleaning, e.g. trays or baskets. Observe the operating and loading instructions of the W/D manufacturer.

5.4 Cleaning and Disinfection

5.4.1 Detergents and Disinfectants

It is mandatory to follow the directions provided by the manufacturers of the detergents and disinfectants used. When using machine cleaning, the instructions provided by the manufacturer of your washer-disinfector must be observed as well.

The warranty excludes any malfunction related to the disinfectants and processes used, including the effectiveness of the disinfectants.

Using disinfectants containing Glucoprotamin® can adversely affect silicone materials. Be sure to observe the manufacturer's instructions prior to use!

Silicone oils and greases can impair the function and must not come into contact with the silicone product.

Amines can have a negative effect on the performance of the product. It is recommended to avoid these during cleaning and sterilization.

During cleaning, it should be ensured that there are no residues of the cleaning or disinfection agent on the product. Otherwise the material may harden, and the product may be damaged.

Use of disinfection agents containing chlorhexidine is not recommended.

Care should be taken to avoid soap solutions or alkaline cleaning agents if a citric acid-based neutralizing agent is to be used subsequently.

Disinfecting agents comprising phenol derivatives must not be used.

We recommend using detergents and disinfectants by the following company:

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We recommend SinAlc sensitive as an additional alcohol-free disinfecting agent.

5.4.2 Machine Processing

Machine cleaning (program with thermal disinfection during the final rinse) should be preferred.

For example, "neodisher® MediClean forte" can be used for machine cleaning.

The following procedure is recommended for machine cleaning:

1. Prewash with cold mains water for 2 minutes.
2. Dry product.
3. Wash the product for 5 minutes using warm mains water with a temperature of 55 °C and 0.5 % cleaning agent (MediClean forte).
4. Dry product.
5. Rinse product for 3 minutes using cold, unfiltered water.

6. Dry product.
7. Rinse product for 2 minutes using cold, unfiltered water.
8. Dry product.

5.4.3 Manual Processing

If only manual cleaning is possible, the neutral electrodes should be immersed in a disinfecting solution containing cleaning components. Be sure to observe the specified immersion time (do not exceed it either!).

The neutral electrode can be “wipe-disinfected”. Be sure, however, to use silicone-compatible disinfectants only.

5.4.4 Final Rinsing

We recommend using demineralized water for the final rinse both for machine cleaning and for manual cleaning.

5.5 Inspection, Functional Check, Maintenance

5.5.1 NE with Non-Split Conductive Surface

NE with non-split conductive surface offer no possibility of monitoring the skin contact resistance between neutral electrode and patient. This means that the contact quality cannot be checked. Nevertheless, the NE monitor of the electrosurgical unit shows a color change from red to green (red = no NE connected / green = NE connected) already upon insertion of a NE with non-split conductive surface. However, this means only that the inserted NE with a non-split conductive surface is in a technically perfect condition and has no line break – it does not mean that the application on the patient is of good quality!

The NE must be carefully applied to the patient, which is not supported by the NE monitor.

The neutral electrode is ready for operation only when applied to the patient.

5.5.2 NE with Split Conductive Surface

In case of NE with a split conductive surface, the quality of the contact on the patient’s skin is checked. Only when the application is perfect does the color change from red to green (red = insufficient contact quality / green = good contact quality).

Subsequently, this neutral electrode is ready for operation.

Only the split NE together with a compatible NE monitor allows detection of unintentional removal of the NE from the patient. If the alarm surface is reached during detachment, the HF energy is deactivated only in case of a compatible system.

5.6 Packaging

Approved sterilization packaging (e.g. conforming to EN 868, ISO 11607) must be used for sterilization, subsequent transportation and storage.

5.7 Sterilization

5.7.1 Before Sterilization

- Inspect the neutral electrodes visually and check them for insulation damage, cleanness and integrity.
- Clean and disinfect the neutral electrodes, rinse them with distilled water and dry them carefully.

! WARNING**Peril to life or danger of serious physical injury due to prion infection!**

- Neutral electrodes that may have come in contact with prions must be destroyed at once (danger of CJD contamination)! Never reuse them in such case!

! WARNING**Risk of infection in case of non-sterile handling!**

Improper sterilization and non-sterile use of the neutral electrodes can pose severe health hazards to patients.

Sterilization must be carried out according to a validated steam sterilization process, for example in a sterilizer complying with EN 285:2009 and DIN EN 13060, validated in accordance with ISO 17665-1:2006.

NOTICE**Risk of neutral electrode damage!**

- Do not use hot air for neutral electrode sterilization!
 - Do not use the “STERRAD®” system for neutral electrode sterilization! (STERRAD is a registered trademark of Johnson & Johnson Inc.).
 - Do not use “flash” autoclaving!
 - Do not use irradiation sterilization, gas sterilization, or formaldehyde sterilization!
-
- Note that high-temperature sterilization and long cycle times reduce the life of the neutral electrodes.
 - Ensure sufficient drying.
 - Follow the sterilization recommendations of the manufacturer of your sterilizer regarding handling, loading and parameter setting.
 - The operator must ensure that sterility is maintained after the sterilization process.

- Only use the sterilization methods listed below. Do not use other sterilization procedures.

Method	Exposure temperature	Exposure time	Drying time
Fractionated vacuum (dynamic air removal)	min. 132 °C	min. 3 min	min. 30 min

- Steam sterilizer complying with EN 13060 or EN 285
- validated according to DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance evaluation (PQ))

5.8 Storage and transport

The items must not be subjected to mechanical stress during transport and storage (e.g. no pressure from above, no storage together with objects that have sharp edges) in order to avoid functional defects or burns. Silicone products must not be stored together with PVC or rubber products.

Parameter	Range
Air humidity	no specifications
Storage temperature	no specifications
Direct sunlight	Protect from sunlight and artificial light during storage.
Shelf life	unlimited

6 Maintenance

6.1 General Notes

The product may only be repaired by Gebrüder Martin or a qualified person or firm expressly authorized by Gebrüder Martin to perform such work.

Modification of the product can lead to unforeseeable risks and is therefore not permitted.

If the repair is carried out by a person or firm specially authorized by Gebrüder Martin, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must show the date of the repair and the details of the person or firm carrying out the work and must be signed.

In all cases where a party other than Gebrüder Martin performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Performing unauthorized actions on the product is strictly prohibited at any time. Non-compliance will void any liability claim against Gebrüder Martin.

7 Environmental information/disposal

7.1 Packaging

Gebrüder Martin will take back the full packaging upon request. Whenever possible, parts of the packaging will be reused.

If this option is not used, the packaging can be disposed in the general paper and household waste.

7.2 Disposal

The product was designed to use as few composite materials as possible. This design approach ensures that it can be recycled to a large degree. This is why we offer to take back the product and to dispose of it properly.

7.3 National regulations

The national regulations and disposal provisions must be observed for all disposal measures.

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Martin Nederland/Marned B.V.
KLS Martin Taiwan Ltd.
KLS Martin UK Ltd.