NOBETEC[®]

ZOE-cement

106913-2020-08-12

0413

DIRECTA AB
P.O. Box 723, Finvids väg 8
SE-194 27 Upplands Väsby,
Sweden, Tel: +46 8 506 505 75

DIRECTA

Device description

Nobetec is a ZOE cement, temporary cement for restorative dental treatment. Powder and liquid are mixed into a cement.

Composition

Powder: Zinc oxide 65%, colophony 35%

Liquid: Eugenol 60%, balsam of Canada 35%, Balsam of Peru 5%

REF

1250100 – Powder, 100 gram

1250200 – Liquid, 20 ml

EN ISO 3107:2011

Classification: Type 1.

Indications for use

For temporary cementation of temporary crowns and bridges made of polycarbonate, composite, tin and aluminum. To be used up to 30 days.

Contraindications

- Temporary filling
- Temporary cementation longer than 30 days
- Temporary cementation of permanent bridges
- Patients or dental personnel with known allergy to eugenol
- Patients or dental personnel with known allergy to other included materials
- Direct contact with pulp

Warning

- Direct exposure/contact to the pulp tissue may harm/irritate the pulp because of the eugenol. If Nobetec unintentionally gets in contact with the pulp, and the pulp get irritated, keep the pulp under control
- Do not use as a permanent solution, since the high precision of the construction/power of adhesion cases difficulties in removal of cement.
- May be dangerous by ingestion and may cause sensitization by skin contact. If on skin: wash with plenty of soap and water. If skin irritation or rash occurs, get medical advice/attention.

Pre-caution

- Wear protective gloves and mouthguard when working with the product
- Before temporary cementation close to the pulp, place a layer of normally used isolation material to protect the pulp

Recommended storage conditions

- Keep the containers tightly closed
- Store at <25°

Shelf life

- The shelf-life of powder and liquid is 3 years from the date of manufacture under the recommended conditions of storage
- Expiry date: see product packaging

Disposal

- Used material shall be disposed as risk waste
- Unused powder and liquid shall be disposed according to information in safety data sheets.

Instructions

Working time: minimum 1 minute

Setting time: 5-10 minutes

- When preparing a temporary post crown, that is an polycarbonate crown combined with a screw post, the screw must go up into the crown about 5mm, the Nobetec cement there can be of rather firm consistency.

Mixing

- Use the supplied measurement spoon. A half measurement spoon corresponds approx. 0,5 ml.
- Start with ½ measurements of powder + 2 drops of liquid.
- If required, increase the amount of liquid to adjust for desired consistency.
- Mix minimum 1 minute on a glass slab or a mixing pad with a spatula
- The paste stays mouldable for 5-10 minutes after mixing.
- To shorten the setting time; add a few drops of water to the powder before mixing with liquid.
- A high degree of humidity in the air or contact with saliva may also shorten the setting time.
- By adding of a few drops of sesame oil or Vaseline the setting time is lengthened.

Application

- Dry the prepared teeth and the surface of the restorations.
- Apply a thin layer of the mixed cement to the internal surfaces of the temporary restoration.
- Firmly seat the restoration in the mouth.
- After the material has set remove excess material.

Removal and cleaning of site

- To prevent weak bond strength of the permanent cementation, when resin based cement is used, the surface needs to be cleaned from any possible residue of eugenol with hydrogen peroxide/ethanol and etched with phosphoric acid, to prepare the permanent cementation.

Cleaning of measurement spoon

Use pure alcohol or alcohol containing detergents.

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