KULZER
MITSUI CHEMICALS GROUP

Palajet®

DE

GB

FR

ES

IT

PT

RU

- DE** Betriebsanleitung
- GB** Operating Instruction
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- RU** Инструкция по эксплуатации

PALA®

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

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Operating instruction – method and instructions for use

Injection process for manufacturing dentures

Contents

1	Application	2
1.1	In general	2
1.2	EC – Declaration of Conformity	2
1.3	Explanation of symbols and accompanying word(s)	2
2	Safety Notes	3
3	Designated Use	3
4	Contents of delivery and accessories	3
5	Equipment Characteristics	3
5.1	Operating Elements	4
6	Starting the apparatus	4
6.1	Installation	4
7	Description of Operation	4
7.1	Preparation	4
7.2	Investment	4
7.3	Isolation	5
7.4	Preparing acrylic	5
7.5	Injecting	5
7.6	Polymerising	6
7.7	Deflasking	6
7.8	Preparation for following processes	6
8	Cleaning – Maintenance – Repairs	7
8.1	Disposal	7
9	Technical data	7
10	Sources of error	7
11	Name plate	7
12	Service	8
13	Document history	8

1 Application

1.1 In general

Palajet® is a registered trademark of Kulzer GmbH, 63450 Hanau (Germany).

These operating instructions apply to:

Order no.	Type Features	Date	Valid from Serial no.
66020450	Palajet – Injection unit	2017-07 / 66035914/11	2017-07-1000

GB

1.2 EC – Declaration of Conformity

We, Kulzer GmbH, Leipziger Straße 2, 63450 Hanau (Germany), hereby declare that the design, structure and version introduced to the market of the unit described below correspond to the relevant safety and health requirements of the EU Directive.

This declaration shall become invalid in the event of an alteration made to the unit without our prior authorisation.

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
1.3 Explanation of symbols and accompanying word(s)

Symbols	Accompanying word(s)	Explanation
	WARNING	Identifies conditions or practises that present risk of injury to the patient and or user.
	NOTE	Instructions on using the unit optimally.
		Read the Operating instruction – method and instructions for use.
		Device manufactured in conformity with European directives and guidelines and is labeled with CE mark.

Operating instruction – method and instructions for use

Injection process for manufacturing dentures

2 Safety Notes

WARNING

The safety of persons, operation and processing quality with this laboratory apparatus depends primarily on the actions of the persons who work with the apparatus.

Read the operating instructions carefully before using the apparatus and observe the directions for avoiding errors and resulting damage, particularly damage to health.

In addition to the directions in the operating instructions, the applicable national laws, regulations and guidelines must also be observed for installation and operation.

Maximum lift speed 50 mm / sec. Flip the switch upwards only when the fill cylinder is inserted.

Caution: danger of injury! Before working on the pressure module (pressure plunger) disconnect the apparatus from the compressed-air supply!

NOTE

Keep the operating instructions in an accessible position near the apparatus!

3 Designated Use

The Palajet is a pneumatic injection apparatus for the production of dentures with PalaXpress® ultra and PalaXpress® cold-curing polymer or heat-curing polymers such as Paladon® 65.

Area of application: dental technology, such as laboratories, schools, universities.

The manufacturer does not accept liability resulting from use that is not in accordance with the instructions for use or unauthorised changes.

The apparatus is designed for continuous operation. Personal protective equipment, such as safety gloves, safety glasses etc., must be used at all times for protection from contamination with the processed materials.

4 Contents of delivery and accessories

Palajet contents of delivery

- Injection cylinder including cover
- 1 Palajet Duoflask (incl. allen-wrench (13 mm), opening key for flasks, distance holder, placeholder)
- 12 injection cylinder cartridges
- 12 cover cartridges
- 50 sealing disks
- Measurement cup
- Mixing cup
- Mixing spatula
- Channel wax (33 strips Ø 7 mm and 30 strips Ø 3 mm)
- Timer
- Flask filling basket

Palajet accessories

- Palajet singleflask (incl. allen-wrench (13 mm), opening key for flasks, distance holder)
- Palajet Duoflask (incl. allen-wrench (13 mm), opening key for flasks, distance holder, placeholder)
- Injection cylinder including cover
- Timer
- Flask filling basket


Timer

5 Equipment Characteristics

- The Palajet is a pneumatic injection apparatus for injecting denture material into denture flasks.
- The apparatus is factory-set to a working pressure of 4 bar, which can be controlled with the pressure meter (see paragraph 10 Sources of error).
- The pneumatic feed pressure of max. 10 bar must not be exceeded.
- The plunger speed during injection is 8 – 10 mm / sec.
- The Palajet is designed and manufactured according to the most common safety regulations and guidelines.

5.1 Operating Elements

- 1) Flask holder knob
- 2) Injection cylinder
- 3) Pressure plunger
- 4) Vent seal
- 5) Pressure meter
- 6) Switch for operating the pressure plunger
- 7) Maintenance set
- 8) Compressed air connection
- 9) Adjusting valve for piston speed
- 10) Cylinder bolt for flask wrench

6 Starting the apparatus

On receipt of the apparatus, please check it for transport damage and report any damage to the transport company no later than 24 hours after delivery.

WARNING

Do not operate damaged apparatus under any circumstances!

6.1 Installation

WARNING

The apparatus must be installed in dry areas only and on a solid foundation so that it cannot be tipped.

Weight capacity of foundation: min. 20 kg.

Area: approx. 300 x 280 mm.

Before connecting the compressed air, check that the primary pressure is not more than 10 bar and not less than 5 bar. Also ensure that the compressed air does not contain oil or water. The apparatus is connected to the air supply with a hose (internal \varnothing 8 mm).

7 Description of Operation

During the injection process the mixed acrylic dough is pressed through the injection cylinder into the denture flask cavity by the pressure plunger.

After the filling process is completed, the remaining denture acrylic will escape through the ventilation channel.

The ventilation channel is then sealed by the vent seal.

NOTE

Please watch for deviations when using heat-curing polymers.

The injection process runs automatically after setting the switch. We recommend the Palamat® elite polymerisation apparatus for polymerizing the dentures.

7.1 Preparation

After connecting the compressed air, flip the switch down so that the pressure plunger releases the fill cylinder (pressure plunger moves down). Move the pressure module to the rear as far as it will go. Remove the fill cylinder from the apparatus by screwing anticlockwise and remove it. Unscrew the fill cylinder cover from the fill cylinder.

7.2 Investment

NOTE

Before embedding, make sure that both flask halves fit exactly (check for the same identification number).

Operating instruction – method and instructions for use

Injection process for manufacturing dentures

GB

- 5 -

Lubricate the inside with Vaseline so that the plaster can later be easily removed. The Palajet singleflask (1 model) or the Palajet Duoflask (1 – 2 models) can be used depending on the number of models. The placeholder can be used to save plaster if the Palajet Duoflask is used with only one model.

The fabricated wax denture is embedded with the injection channel positioned at the maxillary tuberculum on the upper jaw or the retromolar triangle on the lower jaw. Remove excess plaster and smooth the surface. Clean excess plaster from the edges of the flask. Use plaster type 3 for embedding.

Applying channel wax:

After the plaster has hardened, attach the wax for the injection channel (Ø 7 mm PalaXpress® channel wax) and the ventilation channel (not required for Paladon® 65).

Position the ventilation channel (with Ø 3 mm PalaXpress® channel wax) in the area of the maxillary tuberculum or retromolar triangle (when using the single flask opposite the injection channel), choosing the shortest distance to the ventilation channel. If two models are used in the Palajet Duoflask, the dentures must be connected with an injection channel on the opposite side of the injection opening in the tuber area.

Insulate the plaster surface with Aislar® (for more detailed information please see Aislar® instruction for use).

Place the flask halves on top of each other and firmly screw them together with the included spanner (do not use a hammer on brass flasks). Fill in plaster through the flask opening so the teeth are covered and leave to harden. Remove excess material so the flask fits between the centring disks of the Palajet.

7.3 Isolation

After the plaster has hardened, preheat the flask in hot water (approximately 90°C (194°F)) without chemical additives in order to soften the wax, separate the flask halves after about 5 minutes and remove the remaining wax with clean hot water. Thinly insulate the plaster at the contact area to the denture plastic twice (allowing each layer to dry) with Aislar® (prevent formation of pools).

After isolation of the plaster roughen the basal area of the teeth with a coarse diamond (remove sanding dust). To improve the bond between the teeth and the denture plastic, apply Palabond® adhesive twice with a brush and allow each layer to dry for 30 seconds. After the second layer the adhesive will remain active for 10 minutes.

Close both flask halves firmly and insert the flask into the Palajet.

NOTE

The sealing surface of the vent seal of the Palajet apparatus should not contact the opening of the flask's ventilation closure.

7.4 Preparing acrylic

Prior to the mixing procedure unscrew the cylinder cover from the fill cylinder and insert the fill cylinder cartridge. In order to receive the acrylic dough, the seal disk with the sealing lip showing upwards must be pressed into the fill cylinder cartridge (as far as it will go). The cover cartridge will be inserted into the cylinder cover.

Product	Mixing ratio
PalaXpress® ultra	2 g : 1 ml 30 g : 15 ml
PalaXpress®	2 g : 1 ml 30 g : 15 ml
Paladon® 65	5 g : 2 ml 30 g : 12 ml

Pour the powder quickly into the reserved liquid, mix well for 30 seconds. Transfer acrylic dough to the filling cylinder immediately and cover.

Remove air bubbles by swinging the mixing cup. When working with PalaXpress® ultra or PalaXpress® pour the acrylic dough in a thin jet without bubbles into the prepared fill cylinder cartridge.

NOTE

To avoid inhomogeneities do not wipe off the mixing cup on the surface of the fill cylinder cartridge.

7.5 Injecting

PalaXpress® / PalaXpress® ultra:

After filling the fill cylinder cartridge wait for the moment of injection.

This moment is reached when the dough maintains a dull surface.

Close the fill cylinder with the cylinder cover and turn it into the Palajet.
Lock the flask in position with the flask holder knob. After pulling the pressure module as far as it will go, the switch is reversed while the fill cylinder is inserted and the denture acrylic is injected with help of the pressure plunger.

As soon as material is pressed out of the ventilation channel, close the vent seal.
This ensures that the pressure is built up in the flask.
Leave the flask in the Palajet under pressure for 5 minutes (switch up).
Set the timer to 5 minutes.

Paladon® 65:

From the start of the mixing process the waiting period before injection at 23°C (73°F) is 20 minutes.
Waiting times are dependent on the room temperature and the mixed volume.
Close the fill cylinder with the cylinder cover and turn it into the Palajet.
Lock the flask in position with the flask holder knob.

After pulling the pressure module as far as it will go, the switch is reversed while the fill cylinder is inserted and the denture acrylic is injected with help of the pressure plunger.
Leave the flask in the Palajet under pressure for 25 minutes (switch up).
Set the timer to 25 minutes.

NOTE

For operation of the timer see the separate instructions.

7.6 Polymerising

For the polymerisation of the material the flask must be removed from the injection apparatus. Place the switch in the starting position (down); the pressure plunger retracts from the injection cylinder. Once it is completely down, push the pressure module to the rear as far as it will go.

After unscrewing the fill cylinder the flask can be removed, placed on the flask filling basket and then inserted into the Palamat® elite polymerisation unit. Disassemble the filling cylinder immediately and clean it. Acrylic residues can be removed much more easily from the filling cylinder and cover inserts after they have completely cured.

(See paragraph 7.8 Preparation for following processes.)

PalaXpress® / PalaXpress® ultra:

The time required for polymerisation is 40 min at 55°C (131°F) with a pressure of 2 bar in Palamat® elite / Palamat® practic.

Paladon® 65:

Rapid polymerisation takes 20 min at 100°C (212°F).
Kulzer water-bath polymerisation takes approx. 2.5 h at 70°C – 100°C (158°F – 212°F).
Slow polymerisation takes 10 h at 90°C (194°F).
For more detailed information concerning the polymerisation please see instruction for use of acrylic curing material.

7.7 Deflasking

After polymerisation, the flask is removed together with the flask filling basket from the Palamat® elite.
Prior to deflasking the denture allow the flask to cool to room temperature, then deflask the denture.
Remove the plaster from the flask halves with a plastic hammer.

NOTE

Do not use a metal hammer! Metal tools will deform the flask; this can result in inaccurate fitting.

NOTE

The final occlusion must be checked prior to removing the denture from the model!

7.8 Preparation for following processes

Prepare flask before reusing it and if necessary lubricate with Vaseline. Clean the sealing surfaces as well as the injection and ventilation channel from plaster and plastic remains in order to avoid raising the bite. Turn back the vent seal and remove any acrylic remains.

The fill cylinder and cover cartridge are constructed for multiple use (approximately 5 times).

NOTE

For better sealing and a safe pressure application a new sealing disk must be used for each injection procedure!

As noted above, acrylic residues can be removed from the filling cylinder and cover inserts much more easily after they have completely cured.

NOTE

Take care that the cartridges are not damaged!

Operating instruction – method and instructions for use

Injection process for manufacturing dentures

8 Cleaning – Maintenance – Repairs

Clean the surface of the apparatus with a moist cloth. Make sure that no water can enter. Do not use aggressive cleaning agents. Check the water separator weekly and remove accumulated water as required. Check the compressed-air supply and the compressor if leaks occur.

WARNING

The apparatus may be damaged by moisture!

Inspection of the pneumatic components is recommended every two years (see paragraph 12 Service).

WARNING

Use only approved original spare parts and accessories.

The use of parts from other sources involves unknown risks and is not permitted under any circumstances.

WARNING

The function and safety of the apparatus can only be guaranteed if all required tests, services and repairs are conducted by Kulzer service contractors or by technicians trained by the manufacturer.

Kulzer GmbH will not accept liability for damages resulting from damage to the apparatus resulting from incorrect repair work that was not conducted by Kulzer service contractors or by technicians trained by the manufacturer, or if original spare parts or accessories were not used to replace parts or accessories.

8.1 Disposal

The apparatus is designed for a service life of 10 years.

For the disposal of spare parts, or the unit, please contact the Kulzer subsidiary directly in your country.

9 Technical data

Compressed air connection	5 – 10 bar
Nominal pressure	4 bar
Plunger speed	0 – 50 mm / sec
Compression force on the fill cylinder	approx. 2 KN
Height	approx. 400 mm
Width	approx. 290 mm
Depth	approx. 265 mm
Weight	approx. 13.5 kg

10 Sources of error

Error sources	Remedy
Incorrect operating Pressure	<ul style="list-style-type: none"> – Pull setting knob on the regulator- / filter unit upwards until raster snaps out – Adjust regulation of operating pressure to 4 bar and press setting knob downwards until it snaps in
Too fast or too slow injection speed	<ul style="list-style-type: none"> – Remove protective stirrup above the throttling check valve – Remove locknut on the throttling check valve – Adjust knurled head screw according to the desired speed and fasten the locknut – Replace protective stirrup – Hoisting speed = upper throttling check valve – Sink speed = lower throttling check valve

11 Name plate

Kulzer GmbH		2017-07-1000
Leipziger Str. 2, 63450 Hanau		
P = 5-10 bar	PN = 4 bar	
66020450	Palajet®	
		

12 Service

We appreciate your comments, feedback, and suggestions.

Contact in the countries and more information are available by the displayed QR code or at our website www.kulzer-dental.com

13 Document history

2011-09	First edition.
2013-05	Editorial revision.
2017-02	Editorial revision paragraph 1.1, 1.2, 5.1, 6.1, 7.6, 11 and 12.
2017-07	Document version 11. Changed name manufacturer, new address and logo of the manufacturer.

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