D

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GB

**Manufacturer Information on Reprocessing of Medical Devices according to
DIN EN ISO 17664:2004**

Manufacturer: HAGER & WERKEN GmbH & Co. KG
Ackerstraße 1
47269 Duisburg, Germany
T +49 (203) 99269-0
F +49 (203) 299283
www.hagerwerken.de

Valid as of: 11/2016 (Rev: 2)

Reflect-Rhod mirror **REF 355 480 – 355 483**
Silicon handle separately **REF 355 487**

Warning notices: No special notices

Limitation of reprocessing: Frequent reprocessing has only a small impact on this instrument. The instrument's shelf-life is mainly determined by wear and damage by use.

Place of handling: Clean surface with a disposable cloth or paper towel.

Storage and transport: No special requirements.

Préparation for cleaning: Remove silicone handle.

Cleaning: Only use suitable cleaning agents, carefully following the manufacturer's instructions. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93°C in mechanical cleaning.

Disinfection: Only use suitable disinfectant solution, carefully following the manufacturer's instructions. After disinfectant, remove disinfectant solution thoroughly under running water. Do not exceed 93°C in mechanical disinfection.

Maintenance: No special requirements.

Packaging: Put the product into the enclosed fabric bag before sterilisation, in order to protect the mirror face from damages (scratches). Standardised packaging material may and should be used for sterilisation process.

Sterilisation: Vapour sterilisation at 134°C for 5 minutes or
vapour sterilisation at 121°C for 15 minutes.

Control / Functional check: Sight check on damages, wear and deformation.

Storage: No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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Fabricant: HAGER & WERKEN GmbH & Co. KG