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POLYPROPYLENE SUTURE MATERIAL**• SERAPREN®****WARNINGS AND PRECAUTIONS**

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).
- The product is intended for single use only.
- Do not resterilize or reuse the product as critical changes in performance and safety are possible.
- Do not use after expiration date.
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or accidentally opened prior to use.
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

STERILE EO

**DESCRIPTION AND INTENDED PURPOSE**

SERAPREN® is a synthetic, monofilament, nonabsorbable surgical suture. The raw material for the suture is composed of isotactic polypropylene. SERAPREN® is dyed with phthalocyanine blue, Colour Index Nr. 74160. It is available in various gauge sizes and lengths.

Polypropylene sutures are supplied in combination with various types of needles made of surgi-

cal steel. Pledgets to reinforce the body's own structures may be enclosed with SERAPREN®. It may also be provided in combination with specific components for intracutaneous sutures.

The suture material complies with the current edition of the European Pharmacopoeia monograph (EP) "Sterile nonabsorbable sutures (fila non resorbilia sterilis)".

SERAPREN® is indicated for use in soft-tissue approximation and/or for ligation in cases where nonabsorbable suture material is indicated, including use in cardiovascular and neurological procedures as well as in microsurgery and ophthalmic procedures. Polypropylene sutures can also be used as holding sutures and for tagging.

INSTRUCTIONS FOR USE

Due account should be taken of the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound when selecting and using the suture material. Standard knot-tying techniques should be used to ensure knot security. In order to prevent damage to the needle while suturing it is recommended that the needle be grasped towards the swaged end at a point between the end of the first third and the mid-point.

ACTION

The function of the surgical suture in wound care is to create tissue-to-tissue approximation or to attach tissue to allogeneic material. The ligature, or tie, as a special form of suture serves to occlude lumens and hollow organs.

After an initial minimal tissue reaction, polypropylene suture material is slowly and progressively encapsulated by connective tissue. It is nonabsorbable and therefore retains its tensile strength indefinitely. Due to its monofilament structure SERAPREN® may also be used in tissue at risk of infection.

CONTRAINDICATIONS

The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The demands placed on the suture material and the required techniques will vary with the indication. The user should therefore be familiar with the various surgical techniques and procedures and should consider the characteristics when selecting the suture material. As with any foreign body SERAPREN® may give rise to calculus formation with prolonged contact with the biliary or urinary tracts. The use of suitable ancillary measures is recommended when suturing within infected tissue. The surgical instruments should be handled with care to avoid damage to the suture. Bending the needles or grasping them outside the recommended area may result in needle breakage. Needles should be handled and disposed of with particular care due to the risk of infection from needle stick injuries.

ADVERSE REACTIONS

Even with correct use of SERAPREN® the following adverse effects may occur:
Transient localized irritation, transient inflammatory foreign-body reactions, the formation of fistula or granuloma.
Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

INSTRUCTIONS FOR REMOVAL FROM THE PACK

Minibox: For trouble-free removal the suture should be withdrawn from the minibox vertically and not at an angle.

Minibox – microsutures in protective tube: Detach the tube from the minibox fixation; grasp needle with needle-holder and withdraw the suture from the protective tube; if a **shorter suture** is needed, it should be cut to the required length together with the protective tube; with **double-armed sutures** 2 separate sutures can be removed, if required, by cutting the suture loop at the end of the protective tube.

MEANING OF THE SYMBOL

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG