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CE 0482

## Angular Clasps SPEIKO®

one-armed left, size 1; one-armed right, size 2

### Composition:

Stainless steel (type 1.4301), contains nickel, thickness: 1 mm Ø

### Fields of application:

- For manufacturing one-armed retaining elements in partial dentures
- For repairing/expanding one-armed retaining elements in partial dentures

### Contra-indications:

- When processing, pay attention to the pairing of materials, as adverse metal pairings can lead to corrosion or metal discolorations.
- Metal allergy

### Side effects/interactions:

None known if used according to intended purpose.

### Dosage, type and duration of treatment:

The acute angle of the prefabricated angle brackets is placed proximally between the natural tooth and the denture saddle/denture tooth. One end is bent in on the natural tooth, while the other free end is bent so it can be secured/cast in a torsion-resistant way.

Sharp edges can occur if the bracket is manipulated with the corresponding pliers. These sharp edges must be deburred and rounded. Any adhesion of the retaining element is ruled out for the thermal processing and machining.

### Notes:

Before inserting the complete assembly in the patient, it must be visually inspected to verify its perfect condition, otherwise the risk of injury with the corresponding danger of infection can occur with the dentist, dental technician and/or patient.

### Trading formats:

Angular Clasps SPEIKO®, one-armed left, size 1	Content: 10 pieces	Art. no. 1804 (HEKO no. 5521.0001)
Angular Clasps SPEIKO®, one-armed right, size 2	Content: 10 pieces	Art. no. 1805 (HEKO no. 5521.0002)

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