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1. Produktbeschreibung

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Instruction for use cmf etch

1. Product description

Dentine-preserving, buffered, phosphoric acid-containing etching gel for the total-etch technique.

2. Composition

water, phosphoric acid (35 % wt), phosphoric salt, gel former, colorant

3. Indication

Acid etching of enamel and dentin surfaces prior to the application of primers, adhesives, or sealers, etc.

4. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

5. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

6. Interactions

Underfilling materials based on glass ionomer and calcium hydroxide are etched by cmf etch.

7. Processing stages

7.1. Drying

Due to the following application of the filling material, a drainage of the working field is required. Avoid any contamination of the cavity with blood or saliva. A rubber dam is recommended.

7.2. Cavity preparation

Prepare the cavity as usual. Clean and dry the cavity.

7.3. Etching of tooth substance (total-etch technique)

7.3.1. Use a new syringe tip for each patient. Before using in the mouth, put a small amount of the gel on a cotton wool pad to familiarise yourself with the viscosity.

7.3.2. Apply cmf etch to enamel and dentine. Allow to have an effect on enamel and dentine for 30 seconds. The buffered etching gel reduces the risk of excessive demineralisation. Rinse and dry thoroughly. Avoid overdrying. Sufficiently etched enamel appears matt and white.

7.4. Adhesion

The adhesive has to be conditioned and applied in correspondence to the instruction for use of the product used. SAREMCO recommends cmf adhesive system, els duobond, els unibond and apt unibond.

7.5. Refilling the cmf etch syringe 2.5 ml

If necessary, cmf etch 2.5 ml can be refilled with cmf etch 25 ml. To do that, remove the caps from both syringes. Place the adapter supplied firmly on the 25 ml syringe and on the 2.5 ml syringe and screw tight. Once the syringes are firmly in place, slowly withdraw the plunger of the 2.5 ml syringe to collect the desired amount of material.

Avoid unintentional leakage of the material. Carefully separate the syringes from the adapter. When not in use, always keep the syringes closed with the caps.

8. Storage

cmf etch has been developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Bring to room temperature before use.

9. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

10. Precautionary measures

Close containers after each use with the right lid. cmf etch is corrosive! Avoid contact with skin, oral mucosa, eyes, clothes, metal etc. Keep out of reach of children. The principles of good personal and safe working practice apply.

11. Emergency measures

In case of direct contact with skin, oral mucosa or eyes, rinse immediately and thoroughly under running water. In case of contact with the eyes, consult an eye specialist.

12. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination.

13. Warranty

Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value is replaced. For further damages, namely that caused by non-compliance with the instructions for use of other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

14. Production / distribution

SAREMCO Dental AG
Gewerbstrasse 4
CH-9445 Rebstein / Switzerland
Tel: +41 (0) 71 775 80 90
Fax: +41 (0) 71 775 80 99
info@saremco.ch
www.saremco.ch

Edited 09/2020 | D600165

Class IIa medical devices
Class III medical devices (Canada)

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