


TOPDENT PIRANHA

CARBIDFRÄSER/ CARBIDE BURS



Die Carbid Fräser eignen sich hervorragend zum Anrauen von Kunststoffoberflächen aller Art.

Die einzigartige Ausrichtung und Formgebung der Carbid-Oberflächenpunkte ermöglichen durch geringen Abtrag ein größtmögliches Volumen (Retentionen) an der Oberfläche sowie ein besseres, laufruhigeres Abtragen als mit herkömmlichen Fräsern. Auch sehr gut zum Beschleifen der Gipsränder von Modellen geeignet. Blaue Schutzlackschicht kann sich bei Verwendung lösen! Dies hat keinen Einfluss auf Qualität, Funktion bzw. Schleifleistung.




Best suitable for roughen acrylic surfaces. Due to slight abrasion the unique alignment and design of the carbide surface points create a maximum volume (retentions) on the surface as well as a better abrasion compared to common burs. Also suited for grinding plaster margins of models. Attention: the blue protection layer may detach during use! This does not affect quality, function or grinding performance.

www.kkd-topdent.de



- für die Ausarbeitung und Vorbereitung von Unterfütterungen und Erweiterung von Prothesen
- erzeugen raue Oberfläche mit hervorragender Retention
- auch geeignet für Gips
- Reduktion zusätzlicher Retentionsmittel
- extrem hohe Standzeit
- laufruhig auch bei starker Beanspruchung
- problemlose Reinigung -> einfach abdampfen

- for preparation of relinings and enlargements of dentures
- produce a rough surface with excellent retention properties
- also suited for plasters
- reduction of additional retention agents
- extremely long lasting
- smooth running also during heavy use
- easy to clean -> just evaporate



Form Shape	Kugel Ball	Torpedo Torpedo	Kegel spitz Taper	Zylinder flach Cylinder	Zylinder rund Ball Nose	Flamme Flame	Flamme lang Bud
REF	20544	20542*	20546	20545	20543*	20541*	20547
Ø Kopf/ Head Dia.	8 mm	5 mm	10 mm	8 mm	8 mm	8 mm	8 mm
Kopflänge/ Head Length	8 mm	20 mm	16 mm	14 mm	16 mm	15 mm	19 mm
Gesamtlänge/ Overall Length	37 mm	51 mm	46 mm	44 mm	45 mm	44 mm	49 mm

Set*:
je/ ea. 1 St./ pc.
20541, 20542,
20543
REF 20540S